


Autoclavable Temperature Probes

400-AC

A complete range of Autoclavable Temperature Probes for accurate measurement of Body Surface and Core Temperatures



	YSI 400	Philips
		
Application/Size	P/N	P/N
Paediatric Oesophageal / Rectal  3mm Tip	0212010	0212019
Adult Oesophageal / Rectal  4mm Tip	0212011	0212020
Paediatric / Adult Skin / Surface  10mm Disc	0212022	0212031

- YSI 400 Compatible
- CE0086 Certified

- Latex Free
- Brand: Viamed

- Accuracy: +/- 0.1°
- Length: 2.8m