

Digital Thermometers

Hi-Care Digital Thermometer - Flex Tip



CODE : THRDFT

BAR CODE: 6009705230938

Specifications:

- Temperature range: 32°C - 42°C
- Display L°C: temperature below 32°C
- Display H°C: temperature above 42°C
- Accuracy: $\pm 0.1^{\circ}\text{C}$
- Display mode: LCD (3.5)
- Power: 0.15mW (working)
- Battery: AG3 1.5V, ordinary button battery
- Alarm: thermometer beeps after temperature taken (approx 60 seconds)
- Soft Flexible tip
- Waterproof
- Environmentally friendly - containing no mercury

Packing Details:

Single Unit – 1 complete set in Cardboard Pack
(180mm x 47mm x 15mm), 40 grams

Carton – 25 Pieces (207mm x 100mm x 190mm), 1kg


Digital Thermometers



prevention • monitoring • therapy • professional

Flexible Tip Thermometer

Model : TG380 Quite



- ≥ 10 seconds
- Flexible tip
- Fever alarm
- Waterproof
- Auto shut-off
- $^{\circ}\text{C}$ / $^{\circ}\text{F}$ switchable
- Low battery indicator
- Self-diagnosis for malfunction

10

≥ 10 seconds
measurement



Storage Case



Model	Display box	Qty per carton	Carton volume
TG380 Quite	15 pcs	270 pcs	6.045 cm ³

Digital Thermometers

Hi-Care Digital Thermometer - Flex Tip



CODE : THACD

BAR CODE: 6009705230921

Specifications:

- Temperature range: 32°C - 42°C
- Display $L^{\circ}\text{C}$: temperature below 32°C
- Display $H^{\circ}\text{C}$: temperature above 42°C
- Accuracy: $\pm 0.1^{\circ}\text{C}$
- Display mode: LCD (3.5)
- Power: 0.15mW (working)
- Battery: AG3 1.5V, ordinary button battery
- Alarm: thermometer beeps after temperature taken
- Environmentally friendly - containing no mercury
- Memory
- Automatic shut-off

Packing Details:

Single Unit – 1 complete set in Cardboard Pack
(175mm x 68mm x 17mm), 40 grams

Carton – 360 Pieces (710mm x 320mm x 385mm), 14.7kg