



- Full compatibility of the lead wires enable to be used with original ECG trunk Cables or devices, no extra adapter required.
- Built-in ribbon cable for patient comfort and ease of use
- Easy to read labeling facilities easy connection and lead placement
- Unique light blue color indicated the single use

3.1 Disposable ECG Lead Wires

Descriptions	Probe Connector	Lead #	REF# Snap(AHA/IEC)	REF# Grabber(AHA/IEC)	Packaging	Unit code
Covidien						
Disposable 3-lead shielded Ribbon cable, 310cm		3	CDM3-310DS CDM3-310DS-I	CDM3-310DP CDM3-310DP-I	10	B2
Disposable 5-lead shielded Ribbon cable, 90cm		5	CDM5-90DS CDM5-90DS-I	CDM5-90DP CDM5-90DP-I	10	B4
Disposable 5-lead shielded Ribbon V or C-lead cable, 310cm,		5	CDMV5-310DS CDMV5-310DS-I	CDMV5-310DP CDMV5-310DP-I	10	B5
Disposable 3-lead MRI cable, 90cm		3	NA	CDMA3-90DP CDMA3-90DP-I	10	B1
Disposable 5-lead MRI cable, 90cm		5	NA	CDMA5-90DP CDMA5-90DP-I	10	B4
Datascope						
Disposable 3-lead shielded Ribbon cable, 90cm		3	DT3-90DS DT3-90DS-I OEM REF:040-000748-01	DT3-90DP DT3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	DT5-90DS DT5-90DS-I OEM REF:040-000746-01	DT5-90DP DT5-90DP-I	10	B5
Drager /Siemens						
Disposable 3-lead shielded Ribbon cable, 90cm		3	DG3-90DS DG3-90DS-I	DG3-90DP DG3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	DG5-90DS DG5-90DS-I	DG5-90DP DG5-90DP-I	10	B5
Drager /Siemens						

③ Disposable ECG Monitoring Lead Wires



Patient Safety



Patient Comfort














Electrical Performance



Ease of Use



Single patient use (disposable) from patient to monitor, Reduces the risk of cross contamination from reusable cables and leads, improving patient care and cost control.

Descriptions	Probe Connector	Lead #	REF# Snap(AHA/IEC)	REF# Grabber(AHA/IEC)	Packaging	Unit code
Drager /Siemens						
Disposable 6-lead shielded Ribbon cable, 90cm		6	DG6-90DS DG6-90DS-I	DG6-90DP DG6-90DP-I	10	B6
Disposable 3-lead shielded Ribbon cable, 90cm		3	SM3-90DS SM3-90DS-I	SM3-90DP SM3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	SM5-90DS SM5-90DS-I	SM5-90DP SM5-90DP-I	10	B5
GE Multi-link						
Disposable 3-lead shielded Ribbon cable, 90cm		3	MQ3-90DS MQ3-90DS-I	MQ3-90DP MQ3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	MQ5-90DS MQ5-90DS-I	MQ5-90DP MQ5-90DP-I	10	B5
Disposable 5-lead shielded Ribbon cable, V or C lead ,90cm		5	MV5-90DS MV5-90DS-I	MV5-90DP MV5-90DP-I	10	B5
Disposable 6-lead shielded Ribbon cable for Apex, 90cm		6	MQB6-90DS MQB6-90DS-I	MQB6-90DP MQB6-90DP-I	10	B6
Mindray						
Disposable 3-lead shielded Ribbon cable, 90cm		3	MR3-90DS MR3-90DS-I	MR3-90DP MR3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	MR5-90DS MR5-90DS-I	MR5-90DP MR5-90DP-I	10	B5
Nihon Kohden						
Disposable 3-lead shielded Ribbon cable, 90cm		3	NKB3-90DS NKB3-90DS-I	NKB3-90DP NKB3-90DP-I	10	B3
Disposable 6-lead shielded Ribbon cable, 90cm		6	NKB6-90DS NKB6-90DS-I	NKB6-90DP NKB6-90DP-I	10	B5
Philips Veridia Telemetry						
Disposable 3-lead shielded Ribbon cable, 90cm with combiner		3	AT3-90DS AT3-90DS-I	AT3-90DP AT3-90DP-I	10	B5
Disposable 5-lead shielded Ribbon cable, 90cm with combiner		5	AT5-90DS AT5-90DS-I	AT5-90DP AT5-90DP-I	10	B6

Descriptions	Probe Connector	Lead #	REF# Snap(AHA/IEC)	REF# Grabber(AHA/IEC)	Packaging	Unit code
Philips IntelliVue						
Disposable 3-lead shielded Ribbon cable, 90cm		3	HP3-90DS HP3-90DS-I	HP3-90DP HP3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	HPA5-90DS HPA5-90DS-I	HPA5-90DP HPA5-90DP-I	10	B5
Disposable 5-lead shielded Ribbon cable, V or C lead ,90cm		5	HPB5-90DS HPB5-90DS-I	HPB5-90DP HPB5-90DP-I	10	B5
Disposable 3-lead shielded Ribbon cable, 90cm		3	HT3-90DS HT3-90DS-I	HT3-90DP HT3-90DP-I	10	B4
Disposable 5-lead shielded Ribbon V or C-lead cable, 90cm		5	HT5-90DS HT5-90DS-I	HT5-90DP HT5-90DP-I	10	B6
Philips Twin Pin(AA style)						
Disposable 3-lead shielded Ribbon cable, 90cm		3	AAB3-90DS AAB3-90DS-I	AAB3-90DP AAB3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	AAB5-90DS AAB5-90DS-I	AAB5-90DP AAB5-90DP-I	10	B6
Philips DIN style						
Disposable 3-lead shielded Ribbon cable, 90cm		3	HPC3-90DS HPC3-90DS-I	HPC3-90DP HPC3-90DP-I	10	Snap: B2 Grabber: B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	HPC5-90DS HPC5-90DS-I	HPC5-90DP HPC5-90DP-I	10	Snap: B5.3 Grabber: B6.5
Disposable 6-lead shielded Ribbon cable, 90cm		6	HPC6-90DS HPC6-90DS-I	HPC6-90DP HPC6-90DP-I	10	Snap: B6.5 Grabber: B7
Spacelabs						
Disposable 3-lead shielded Ribbon cable, 90cm	NA	3	SLB3-90DS SLB3-90DS-I	SLB3-90DP SLB3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm		5	SL5-90DS SL5-90DS-I	SL5-90DP SL5-90DP-I	10	B5

3.2 Disposable One Piece ECG Lead Wires




OEM & Models	ECG machine Plug	Lead #	REF# Snap(AHA/IEC)	REF# Grabber(AHA/IEC)	Packaging	Unit code
AAMI						
Merlin 78352C/78354C/ 78834CM1001A		3	2340DS 2340DS-I	2340DP 2340DP-I	8	B6
		5	2540DS 2540DS-I	2540DP 2540DP-I	8	B7
Covidien						
适配Covidien心电图转接头		3	23103DS-280 23103DS-280-I	23103DP-280 23103DP-280-I		B5
		5	25103DS-280 25103DS-280-I	25103DP-280 25103DP-280-I		B6.5
Philips						
Merlin 78352C/78354C/ 78834CM1001A/		3	2385DS 2385DS-I	2385DP 2385DP-I	8	B6
		5	2585DS 2585DS-I	2585DP 2585DP-I	8	B7
GE Multi-link						
Eagle,Solar,Dash monitor MAC-Lab Cath Lab Datex-Ohmeda M-Prestn & “E” modules		3	2386DS 2386DS-I	2386DP 2386DP-I	8	B5
		5	2586DS 2586DS-I	2586DP 2586DP-I	8	B6
Nihon Kohden						
BSM-2300 lifescopel I BSM-4100 iflescope P BSM-5100		3	2309DS 2309DS-I	2309DP 2309DP-I	8	B6
		6	2609DS 2609DS-I	2609DP 2609DP-I	8	B7
Mindray / Datascope						
PM5000 PM6000 BeneView T5/T6/T8		3	2312DS 2312DS-I	2312DP 2312DP-I	8	B6
		5	2512DS 2512DS-I	2512DP 2512DP-I	8	B7
Datascope-Mindray V12, V21		3	2326DS 2326DS-I	2326DP 2326DP-I	8	B6
		5	2526DS 2526DS-I	2526DP 2526DP-I	8	B7
Spacelabs						
Ultraview 90492, 91220, 90496, 91496		3	2396DS 2396DS-I	2396DP 2396DP-I	8	B6
		5	2596DS 2596DS-I	2596DP 2596DP-I	8	B7

3.3 Disposable ECG Radiolucent Lead Wires



Descriptions	Probe Connector	Lead #	REF# AHA (Snap/Grabber)	REF# IEC (Snap/Grabber)	Package	Unit price code
AAMI Disposable Radiolucent Lead Wires,Din Type/0.9m		3	D3-90XDS	D3-90XDS-I	1 PC	B1
		3	D3-90XDP	D3-90XDP-I	1 PC	A8
		5	D5-90XDS	D5-90XDS-I	1 PC	B4.1
		5	D5-90XDP	D5-90XDP-I	1 PC	B3
		10	D10-90XDS	D10-90XDS-I	1 PC	C1
		10	D10-90XDP	D10-90XDP-I	1 PC	B6.6
Siemens Disposable Radiolucent Lead Wires,Euro Type/0.9m		3	SME3-90XDS	SME3-90XDS-I	1 PC	B3
		3	SME3-90XDP	SME3-90XDP-I	1 PC	B1
		5	SME5-90XDS	SME5-90XDS-I	1 PC	B6
		5	SME5-90XDP	SME5-90XDP-I	1 PC	B4.1
		10	SME10-90XDS	SME10-90XDS-I	1 PC	C3
		10	SME10-90XDP	SME10-90XDP-I	1 PC	B8.1
Siemens Disposable Radiolucent Lead Wires,SMCTA/SMCTB Type/0.9m		5(Limb Leads)	SMCTA-90XDS	SMCTA-90XDS-I	1 PC	B4.1
		5(Limb Leads)	SMCTA-90XDP	SMCTA-90XDP-I	1 PC	B2
		5(chest leads)	SMCTB-90XDS	SMCTB-90XDS-I	1 PC	B4.1
		5(chest leads)	SMCTB-90XDP	SMCTB-90XDP-I	1 PC	B2

Descriptions	Probe Connector	Lead #	REF# AHA (Snap/Grabber)	REF# IEC (Snap/Grabber)	Package	Unit price code
GE Disposable Radiolucent Lead Wires,Vs Type/0.9m		3	MQD3-90XDS	MQD3-90XDS-I	1 PC	B4.1
		3	MQD3-90XDP	MQD3-90XDP-I	1 PC	B2
		5	MQD5-90XDS	MQD5-90XDS-I	1 PC	B5.3
		5	MQD5-90XDP	MQD5-90XDP-I	1 PC	B4
		10	MQD10-90XDS	MQD10-90XDS-I	1 PC	C2
		10	MQD10-90XDP	MQD10-90XDP-I	1 PC	B8
Spacelabs Disposable Radiolucent Lead Wires,LL Type/0.9m		3	LL3-90XDS	LL3-90XDS-I	1 PC	B2
		3	LL3-90XDP	LL3-90XDP-I	1 PC	B1
		5	LL5-90XDS	LL5-90XDS-I	1 PC	B5
		5	LL5-90XDP	LL5-90XDP-I	1 PC	B4
		10	LL10-90XDS	LL10-90XDS-I	1 PC	C1
		10	LL10-90XDP	LL10-90XDP-I	1 PC	B6.7
Spacelabs Disposable Radiolucent Lead Wires,SL Type/0.9m		3	SLB3-90XDS	SLB3-90XDS-I	1 PC	B5
		3	SLB3-90XDP	SLB3-90XDP-I	1 PC	B4
		5	SLB5-90XDS	SLB5-90XDS-I	1 PC	B7.9
		5	SLB5-90XDP	SLB5-90XDP-I	1 PC	B6.6
		10	SLB10-90XDS	SLB10-90XDS-I	1 PC	C3.8
		10	SLB10-90XDP	SLB10-90XDP-I	1 PC	C1
Philips Disposable Radiolucent Lead Wires,AA Type/0.9m		3	AAC3-90XDS	AAC3-90XDS-I	1 PC	B4
		3	AAC3-90XDP	AAC3-90XDP-I	1 PC	B2
		5	AAC5-90XDS	AAC5-90XDS-I	1 PC	B6
		5	AAC5-90XDP	AAC5-90XDP-I	1 PC	B5
		10	AAC10-90XDS	AAC10-90XDS-I	1 PC	C3.8
		10	AAC10-90XDP	AAC10-90XDP-I	1 PC	C1
COVIDIEN Disposable Radiolucent Lead Wires,CDM Style/0.9M		3	CDMA3-90XDS	CDMA3-90XDS-I	1 PC	B2
		3	CDMA3-90XDP	CDMA3-90XDP-I	1 PC	B1
		5	CDMA5-90XDS	CDMA5-90XDS-I	1 PC	B4.1
		5	CDMA5-90XDP	DMA5-90XDP-I	1 PC	B3
COVIDIEN Disposable Radiolucent Lead Wires,CDM Style/1.8m		5	CDMA5-180XDS	CDMA5-180XDS-I	1 PC	B6
		5	CDMA5-180XDP	CDMA5-180XDP-I	1 PC	B4.1
COVIDIEN Disposable Radiolucent Lead Wires,CDM Style/3.1m(V2-V6,C2-C6)		5	CDMV5-310XDS	CDMV5-310XDS-I	1 PC	B6.6
		5	CDMV5-310XDP	CDMV5-310XDP-I	1 PC	B6



4 ECG Lead Adapter

The unique design of the ECG lead wires adapters enable the hospitals to use the one type lead wires to fit major brands patient monitoring system running in one hospital by attaching the proper ECG adapter to the trunk cable. No confusion on lead wires any more !

Specification

Cable material	ABS Jacket, Light gray housing
Connector &Yoke	Nylon plug & Golden plated pins
Link Tie	12cm length with trunk cable clip

Compatibility

Major brands ECG monitoring leads system