


We also supply Invasive Blood Pressure (IBP) Interface Cable

- Adaptor Cable
- Convert Cable
- Transducer adaptor

Compatibility

- GE Marquette/Datex-Ohmeda
- Philips / HP
- Siemens/ Drager
- Nihon Kohden
- Spacelabs
- Mindray/ Datascope
- ...



Application Diagram

For detail information , please contact our sales or send emails to sales@unimed.cn.



Unimed Medical Supplies Inc.
Floor 5&6, Bld#8, Nangang 3rd Industrial Park,
Tangtou, Shiyan, Shenzhen 518108, China
T 86 0755 2669 5165 info@unimed.com
F 86 0755 2669 7984 www.unimed.cn

深圳市尤迈医疗用品有限公司
广东省深圳市宝安区石岩街道南岗
第三工业园8栋5-6楼 邮编:518108
T 0755 2669 5165 info@unimed.com
F 0755 2669 7984 www.unimed.cn



High Quality Transducers

High performance, cost effective.



Disposable Invasive Blood Pressure Transducer

How and Where To Use


- Intended for arterial and venous blood pressure monitoring and blood sampling testing, providing assistance for clinical diagnosis and treatment.
- A sterile and disposable medical device exposing in the vascular lumen for a short time. It can continuously monitor instantaneous hemodynamic changes, and shows the pressure waveforms, record and store values, providing accurate information for immediate treatment.

Features

- Continuously auto-monitoring the instant hemodynamic changes.
- Great sensitivity, wide- monitor range, and accurate measurment.
- Five types of popular connectors, suitable for most monitors in the market.



- Advanced CM OS chip guarantees superior sensitivity.



Technical Specification

Connector	Standard lock connector
Flush speed	≥3ml±1ml/1h/300mmHg
Pressure measurement range	-50~300mmHg
Sensitivity	5.0μv/v/mmHg±5%
Leakage current	<10μA
Marker	Artery-red Venous-blue


Quality control

- Accuracy
- Airproof

REF


Transducer Types

Type A: Single-channel set




REF	Description	Qty/Case	Price
DP-UT	Disposable Blood Pressure Transducer, UTAH connector, single channel	50	B4.5
DP-MX	Disposable Blood Pressure Transducer, MEDEX connector, single channel	50	B4.5
DP-ED	Disposable Blood Pressure Transducer, EDWARDS connector, single channel	50	B4.5
DP-BD	Disposable Blood Pressure Transducer, BD connector, single channel	50	B4.5
DP-BB	Disposable Blood Pressure Transducer, B.BRAUN connector, single channel	50	B4.5

Type B: Double-channel set



REF	Description	Qty/Case	Price
DP-UT/2	Disposable Blood Pressure Transducer, UTAH connector, double channels	50	C1
DP-MX/2	Disposable Blood Pressure Transducer, MEDEX connector, double channels	50	C1
DP-ED/2	Disposable Blood Pressure Transducer, EDWARDS connector, double channels	50	C1
DP-BD/2	Disposable Blood Pressure Transducer, BD connector, double channels	50	C1
DP-BB/2	Disposable Blood Pressure Transducer, B.BRAUN connector, double channels"	50	C1

Type C: Triple-channel set



REF	Description	Qty/Case	Price
DP-UT/3	Disposable Blood Pressure Transducer, UTAH connector, triple channels	50	C5
DP-MX/3	Disposable Blood Pressure Transducer, MEDEX connector, triple channels	50	C5
DP-ED/3	Disposable Blood Pressure Transducer, EDWARDS connector, triple channels	50	C5
DP-BD/3	Disposable Blood Pressure Transducer, BD connector, triple channels	50	C5
DP-BB/3	Disposable Blood Pressure Transducer, B.BRAUN connector, driple channels	50	C5