

Disposable Neonatal NIBP Cuffs

Product Feature


产品介绍

- Non-woven fabrics and TPU material for maximum patient comfort;
 - Latex free & allergy free;
 - Enviromentally protective materials, stable and durable
 - Stable and adhensive hook & loop tapes, easy to peel & fix;
 - Color coded cuff sizes.
- Cuffs with particular connectors to fit major brands
- monitoring systems
- Meeting different customization needs of customers





1.a) Disposable Neonate TPU Cuff(single & double tube)

一次性新生儿TPU血压袖带(单&双管)

Size# 尺寸号码	Limb Circumference 臂周长	Printing Color 袖套颜色	REF# Single tube 型号(单管)	Cuff Connector 连接端	REF # Double tube 型号(双管)	Cuff Connector 连接端	Package 包装	Price (Box) 单价(盒)
1	3-6cm	Orange	U1681S-C05-10		U1681D-C0505-10		10pcs/box	C3
2	4-8cm	Red	U1682S-C05-10		U1682D-C0505-10		10pcs/box	C4
3	6-11cm	Green	U1683S-C05-10		U1683D-C0505-10		10pcs/box	C5
4	7-13cm	Blue	U1684S-C05-10		U1684D-C0505-10		10pcs/box	C6
5	8-15cm	Brown	U1685S-C05-10		U1685D-C0505-10		10pcs/box	C6

1.b) Disposable Neonate Non-wovon Fabrics Cuff(single & double tube)

一次性新生儿无纺布血压袖带(单&双管)

Size# 尺寸号码	Limb Circumference 臂周长	PrintingColor 袖套颜色	REF# Single tube 型号(单管)	Cuff Connector 连接端	REF # Double tube 型号(双管)	Cuff Connector 连接端	Package 包装	Price (Box) 单价(盒)
1	3-6cm	Orange	U1781S-C05-10		U1781D-C0505-10		10pcs/box	C2
2	4-8cm	Red	U1782S-C05-10		U1782D-C0505-10		10pcs/box	C2
3	6-11cm	Green	U1783S-C05-10		U1783D-C0505-10		10pcs/box	C3
4	7-13cm	Blue	U1784S-C05-10		U1784D-C0505-10		10pcs/box	C3.8
5	8-15cm	Brown	U1785S-C05-10		U1785D-C0505-10		10pcs/box	C4