

Disposable Phototherapy Eye Shields



Features

- Made from very soft velvet fabric foam
- Very comfortable & gentle on baby's skin.
- Adjustable band for optimum comfort & a secure fit
- Ergonomically designed eye protection pad to reduce pressure and increase light protection
- Latex Free



Size 尺寸		Head Circumference 头围	REF# 型号	Color code 色标		Package 包装
Small	小号	20-25cm	UBE-S01-10	Green	绿色	10pcs/box 盒
Medium	中号	26-32cm	UBE-M01-10	Orange	橙色	10pcs/box 盒
Large	大号	33-38cm	UBE-L01-10	Blue	蓝色	10pcs/box 盒