



Disposable Temperature Probes

OEM Technology 原厂技术	REF# 型号	OEM REF# OEM型号	Application 应用	Cable length 线长	Probe adapter 探头转接头端	Adapter cable 连接线	Package 包装	Price 价格代码
400 series	TMQ-D-20	M1024222	Skin 体表	75cm/90cm		See the next page of adapter cable 见下页连接线	20	D5
	4009-ER See Next Page		General 成人体腔					
Airshield	TAS30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TAS30-DAG-10		General 成人体腔					
Artema/S&W	TSW30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TSW30-DAG-10		General 成人体腔					
Biolight	TBLT30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TBLT30-DAG-10		General 成人体腔					
Beijing Julongsanyou Infant Radiant Warmer BNT-1000/BNT-2000 /BNT-2000B	T2015-D-10		Skin 体表	150cm		No need 不需要连接线	10	D4
	T2030-D-10		Skin 体表	300cm				
Beijing Julongsanyou Infant Incubator B-800/B-1000/ B-2000/B-3000	TBNT15-D-10		Skin 体表	150cm		No need 不需要连接线	10	D4
	TBNT30-D-10		Skin 体表	300cm				
Drager/Siemens	TDG-D-20		Skin 体表	75cm/90cm		NA	20	D5
	TDG-DAG-20		General 成人体腔	75cm/90cm				
Mindray T5/T8	TMR30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TMR30-DAG-10		General 成人体腔					
Mindray - Datascope	TDT30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TDT30-DAG-10		General 成人体腔					
Ohmeda,10K	TOH30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TOH30-DAG-10		General 成人体腔					
Philips Medical	THP-D-20	21091A	Skin 体表	75cm/90cm		See the next page of adapter cable 见下页连接线	20	D5
	THP-DAG-20	21090A,M1837A	General 成人体腔	75cm/90cm				
Philips Medical	TPM30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TPM30-DAG-10		General 成人体腔					
Siemens/Drager	TSM16-D-10	MP00995(1.6m) MP00997(1.6m)	Skin 体表	160cm		No need 不需要连接线	10	D4
	TSM30-D-10	MP00996(3.0m) MP00998(3.0m)	Skin 体表	300cm				
	TSM16-DAG-10	MP00991(1.6m)	Adult General	160cm				
	TSM30-DAG-10	MP00992(3.0m)	成人体腔	300cm				
	TSM16-DPG-10	MP00993(1.6m)	Pediatric General	160cm				
	TSM30-DPG-10	MP00994(3.0m)	儿童体腔	300cm				
Spacelabs	TSL30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D5.2.5
	TSL30-DAG-10		General 成人体腔					
YSI400	YTSI30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	YTSI30-DAG-10		General 成人体腔					

Disposable Temperature Probes

#4009-ER 400 Series General Purpose/Esophageal/Rectal Temperature Probes



Product Features

- Smooth rounded beveled tip for non-traumatic insertion and withdrawal
- Available in 9 and 12F to cover both pediatric and adult patients
- Compatible with all hypo-hyperthermia equipment providing an accurate, low cost alternative

Product Specifications

- Individually pouched sterile single use probe
- Insert molded moisture resistant "Molex" Style connector
- Temperature accuracy: $\pm 0.1^{\circ}\text{C}$ @ 37°C , $\pm 0.2^{\circ}\text{C}$ @ 25°C to 45°C
- Completely enclosed lead wires for improved stability and hygiene

Ordering Information

Catalog Number	Description	Box Qty
4009-ER	Esophageal/Rectal Probe with 400 series thermistor, 9F	35 each

Disposable Temperature Probes

#10-1700-090 700 Series Skin Temperature Probes



- Reliable temperature accuracy $\pm 0.2^{\circ}\text{C}$.
- Extended length (32") lead wire.
- "Perfect fit" molded connector ensures a trouble-free connection and continuous temperature monitoring reliability.
- Twisted pair wiring and double strain relief provide safety reinforcement.
- Utilizes non-invasive, one-inch adhesive foam disk.
- Not made with natural rubber latex, hypoallergenic foam and adhesive.
- Safe for even the most sensitive skin.