



1 ECG Monitoring Cables& Lead Wires

Specification

Cable material	TPU Jacket, Light gray
Connector &Yoke	Nylon plug & Golden plated pins
Main Cable length	2.5m
Lead wire length	0.9m

Compatibility

Major brands ECG or multi-parameter monitoring system



Shielded Cable & Lead Wire	Photo	REF#	REF # Rule Description
One Piece Fixed Cable (Shielded)		2586P	2: Shielded(1=unshielded) 5: 5-lead(3= 3 lead) 86: factory/model code N: with neonate clip end(Default is adult) See the page below for the figure S: snap end(P=grabber end) I: IEC code (Default is AHA)
Trunk Cable (Shielded)		MQ-2586	MQ: cable yoke type, up to lead wire 2: Shielded(1=unshielded) 5: 5-lead(3= 3 lead) 86: factory/model code Default : AHA code
Lead Wires (Shielded)		MQ5-90P	MQ: lead wire type(Acc to yoke) 6: 6-lead (3= 3 lead,5=5lead) 90: 90cm(wire length) S: snap end(P=grabber end) I: IEC code(Default : AHA code)

Unshielded Cable & Lead Wire	Photo	REF#	REF # Rule Description
Din trunk Cable (1.5mm, Unshielded) Used with Din lead Wires especially for neonate ECG monitoring		D-1540-I	D: Din style cable 1: Unshielded 5: 5-lead(3= 3 lead) 40: Factory/model code I: IEC code (Others see below)
Din lead wires (1.5mm, Unshielded) Used with Din trunk cable only		D5-90S-I	D: Din lead wire 5: 5-lead (3= 3 lead) 90: 90cm(wire length) S: Snap end(P=grabber end) I: IEC code (Others see below)
Neonate lead wires (1.5mm, Unshielded)		D3-90NP-I	D: Din lead wire 3: 3-lead (5= 5 lead) 90: 90cm(wire length) NP=mini grabber end only I: IEC code (Others see below)




ECG Color Code on Electrodes (Regular grabber/Mini grabber/Regular snap)						
Code	AHA(A)	Alternate AHA(H)	Standard IEC(I)	Alternate IEC(IH)	AHA(A)	Standard IEC(I)
Grabber						
Mini Grabber						
Snap						





1.1 ECG Cables & Lead Wires Cross Reference




OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead Yoke Conn.	Lead Wire/w Snap	Lead Wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable
AAMI									
Without resistance BCI Datascope MDE/Invivo Mindray Spacelabs WelchAllyn Protocol		3	D-1340		D3-90S	D3-90P	2340S	2340P	
			D-1340-I		D3-90S-I	D3-90P-I	2340S-I	2340P-I	
		Price	C1		A8	B3	C1	C1	
		5	D-1540		D5-90S	D5-90P	2540S	2540P	
With 1K resistance Same as above			D-1540-I		D5-90S-I	D5-90P-I	2540S-I	2540P-I	
		Price	C2		B4	B6	C2	C2	
		3	D-1340R		D3-90S	D3-90P	2340RS	2340RP	
			D-1340R-I		D3-90S-I	D3-90P-I	2340RS-I	2340RP-I	
		Price	C1		A8	B3	C1	C1	
		5	D-1540R		D5-90S	D5-90P	2540RS	2540RP	
			D-1540R-I		D5-90S-I	D5-90P-I	2540RS-I	2540RP-I	
		Price	C2		B4	B6	C2	C2	
LL style		3	LL-2340		LL3-90S	LL3-90P	2340S	2340P	D-1340
			LL-2340-I		LL3-90S-I	LL3-90P-I	2340S-I	2340P-I	D-1340-I
		Price	C6		B7	B7	C1	C1	C1
		5	LL-2540		LL5-90S	LL5-90P	2540S	2540P	D-1540
			LL-2540-I		LL5-90S-I	LL5-90P-I	2540S-I	2540P-I	D-1540-I
		Price	C6		C2	C2	C2	C2	C2

OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable
Artema/S&W									
8000,Diascope1,-2, Diascope NT, Athena		3	NA	NA	NA	NA	2397S	2397P	D-1397
		Price					2397S-IH	2397P-IH	D-1397-IH
BCI See AAMI									
Biolight(China)									
M7000,M9500 AS Modular		5	NA	NA	NA	NA	2557S	2557P	D-1557
		Price					2557S-I	2557P-I	D-2557-I
AnyView A2/A3/A5/A6/A8 Q3/Q5/Q7/V6 Vital Sign Monitor		3	NA	NA	NA	NA	2367S	2367P	D-1367
		Price					2367S-I	2367P-I	D-1367-I
Bionet(Korea)									
		5	NA	NA	NA	NA	2514S	2514P	D-1514
		Price					2514S-I	2514P-I	D-2514-I
Biosys(Korea)									
		3	NA	NA	NA	NA	2335S	2335P	D-1335
		Price					2335S-I	2335P-I	D-1335-I
Bruker									
		3	NA	NA	NA	NA	2324S	2324P	D-1324
		Price					2324S-I	2324P-I	D-1324-I
		5	NA	NA	NA	NA	2524S	2524P	D-1524
		Price					2524S-I	2524P-I	D-1524-I
Cardell									
		3	NA	NA	NA	NA	2373S	2373P	D-1373
		Price					2373S-I	2373P-I	D-1373-I
CAS MEDICAL									
		3	NA	NA	NA	NA	2339S	2339P	D-1339
		Price					2339S-IH	2339P-IH	D-1339-IH
		3	NA	NA	NA	NA	2364S	2364P	D-1364
		Price					2364S-I	2364P-I	D-1364-I
COMEN									
STAR-8000E STAR-8000H STAR—5000C		3	SM-2337		SM3-90S	SM3-90P	2337S	2337P	D-1337P
		Price	SM-2337-I		SM3-90S-I	SM3-90P-I	2337S-I	2337P-I	D-1337P-I
		5	SM-2537		SM5-90S	SM5-90P	2537S	2537P	D-1537P
		Price	SM-2537-I		SM5-90S-I	SM5-90P-I	2537S-I	2537P-I	D-1537P-I
C100,/C100A/C90 C80/C70/C60/50 C30/C21/C20 NC3/NC5/NC8 NC12/NC10/NC19		3	AA-2347		AA3-90S	AA3-90P	2347S	2347P	D-1347
		Price	AA-2347-I		AA3-90S-I	AA3-90P-I	2347S-I	2347P-I	D-1347-I
		5	AA-2547		AA5-90S	AA5-90P	2547S	2547P	D-1547
		Price	AA-2547-I		AA5-90S-I	AA5-90P-I	2547S-I	2547P-I	D-1547-I





OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable
Colin/Omron(Japan)									
BP 308/308E/408	See AAMI								
BP88S		3	CL-2520		CL3-90S	CL3-90P	2320S	2320P	D-1320
			CL-2520-I		CL3-90S-I	CL3-90P-I	2320S-I	2320P-I	D-1320-I
		Price	C6		B7	B7	C4	C4	C4
		5	CL-2520		CL5-90S	CL5-90P	2520S	2520P	D-1520
			CL-2520-I		CL5-90S-I	CL5-90P-I	2520S-I	2520P-I	D-1520-I
		Price	C6	C2	C2	C5	C5	C5	
BP88 BP306		3	CL-2521		CL3-90S	CL3-90P	2321S	2321P	D-1321
			CL-2521-I		CL3-90S-I	CL3-90P-I	2321S-I	2321P-I	D-1321-I
		Price	C6		B7	B7	C4	C4	C4
		5	CL-2521		CL5-90S	CL5-90P	2521S	2521P	D-1521
			CL-2521-I		CL5-90S-I	CL5-90P-I	2521S-I	2521P-I	D-1521-I
		Price	C6	C2	C2	C5	C5	C5	
Criticare See AAMI									
602A-A, Poet 2, 504US,Scholar,506, 507,508,1100,1508, 8100									
Creative(China)									
		3	NA	NA	NA	NA	2330S	2330P	D-1330
							2330S-I	2330P-I	D-1330-I
		Price					C1	C1	C1
		5	NA	NA	NA	NA	2530S	2530P	D-1530
							2530S-I	2530P-I	D-1530-I
		Price					C2	C2	C2
Datascope /Mindcray See AAMI									
Datascope -Mindray V12, V21		3	DT-2526		DT3-90S	DT3-90P	2326S	2326P	D-1326
			DT-2526-I		DT3-90S-I	DT3-90P-I	2326S-I	2326P-I	D-1326-I
		Price	C6		C2	C2	C4	C4	C4
		5	DT-2526		DT5-90S	DT5-90P	2526S	2526P	D-1526
			DT-2526-I		DT5-90S-I	DT5-90P-I	2526S-I	2526P-I	D-1526-I
		Price	C6	C2	C2	C6	C6	C6	
C2DPM6 DPM7 Passport 8 Passport 12 T1		3	DT-2512 009-003652-00		DT3-90S 0012-00-1503-04	DT3-90P 0012-00-1514-06	2312S 040-000965-00	2312P 040-000964-00	D-1340
			DT-2512-I		DT3-90S-I	DT3-90P-I	2312S-I	2312P-I	D-1340-I
		Price	C6		C2	C2	C4	C4	C1
		5	DT-2512		DT5-90S 0012-00-1503-03	DT5-90P 0012-00-1514-03	2512S 040-000961-00	2512P 040-000961-00	D-1540
			DT-2512-I		DT5-90S-I	DT5-90P-I	2512S-I	2512P-I	D-1540-I
		Price	C6	C2	C2	C6	C6	C2	
Passport 2 Spectrum Spectrun OR		3	DT-2540		DT3-90S 0012-00-1503-04	DT3-90P 0012-00-1514-06	2340S 0010-30-43117	2340P	D-1340
			DT-2540-I		DT3-90S-I	DT3-90P-I	2340S-I	2340P-I	D-1340-I
		Price	C6		C2	C2	C1	C1	C1
		5	DT-2540 0012-00-1502-01		DT5-90S 0012-00-1503-03	DT5-90P 0013-00-1514-03	2540S 0010-30-43116	2540P	D-1540
			DT-2540-I		DT5-90S-I	DT5-90P-I	2540S-I	2540P-I	D-1540-I
		Price	C6	C2	C2	C2	C2	C2	
Dixtal									
Ep12			NA	NA	NA	NA	21070S	21070P	NA
							21070S-I	21070P-I	
		Price					D3	D3	
Drager									
UM1/3/3.1,RM1, PM8010,PM8014, Babyguard8000, Cato,PB8800, Series 8000,Cicero, Parameterbox		3	NA	NA	NA	NA	2377S	2377P	D-1377
							2377S-I	2377P-I	D-1377-I
		Price					C4	C4	C4

OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable
DRE									
WavelineEZ, Touch,Pro		3			MQ3-90S	MQ3-90P	2380S	2380P	D-1380R
							2380S-I	2380P-I	D-1380R-I
		Price			C2	C2	C4	C4	C4
		5			MQ5-90S-I	MQ5-90P-I	2580S	2580P	D-1580R
							2580S-I	2580P-I	D-1580R-I
		Price			C2	C2	C6	C6	C6
		5			MQ5-90S	MQ5-90P	2580ES	2580EP	
							2580ES-I	2580EP-I	
		Price		C2	C2				
ECAFIX									
		5			NA	NA	2552S 2552S-I	2552P 2552P-I	NA
							C2	C2	
Fukuda Danshi									
		3			SM5-90S	SM5-90P	2599S	2599P	D-1599
					SM5-90S-I	SM5-90P-I	2599S-I	2599P-I	D-1599-I
		Price			C2	C2	C6	C6	C6
		12			MQ5-90S MV5-90S	MQ5-90P MV5-90P	21099S	21099P	NA
					MQ5-90S-I MV5-90S-I	MQ5-90P-I MV5-90P-I	21099S-I	21099P-I	
		Price			2XC2	2XC2	D3	D3	
DS5100E,DS5300W, DS5400,DS7001 Mindray/Datascope /Expert		5		NA	NA	NA	2559S 2559S-I	2559P 2559P-I	D-1559 D-1559-I
		Price					C6	C6	C6
GE/Critikon									
Dianmap Pro 1000 MPS Select Protable monitors		3			DX3-90S	DX3-90P	2340RS	2340RP	D-1340R
					DX3-90S-I	DX3-90P-I	2340RS-I	2340RP-I	D-1340R-I
		Price			B7	B7	C1	C1	C1
GE-Datexx-Ohmeda									
AS/3, CS/3, Cardiacap II, Satlite Plus S5 Without Resister for respiration		3			DX3-90S	DX3-90P	2392S	2392P	D-1392
					DX3-90S-I	DX3-90P-I	2392S-I	2392P-I	D-1392-I
		Price			B7	B7	C4	C4	C4
		5			DX5-90S	DX5-90P	2592S	2592P	D-1592
					DX5-90S-I	DX5-90P-I	2592S-I	2592P-I	D-1592-I
AS/3, CS/3, Cardiacap II, Satlite Plus S5 With resisters		Price			C2	C2	C6	C6	C6
		3			DX3-90S	DX3-90P	2395S	2395P	D-1395
					DX3-90S-I	DX3-90P-I	2395S-I	2395P-I	D-1395-I
		Price			B7	B7	C4	C4	C4
		5			DX5-90S	DX5-90P	2595S	2595P	D-1595
					DX5-90S-I	DX5-90P-I	2595S-I	2595P-I	D-1595-I
		Price			C2	C2	C6	C6	C6
		12			MQ5-90S MV5-90S	MQ5-90P MV5-90P	21095S	21095P	NA
					MQ5-90S-I MV5-90S-I	MQ5-90P-I MV5-90P-I	21095S-I	21095P-I	
				2XC2	2XC2	NA	NA		



OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable	
GE/Hellige										
Eagle 1000 Def.. SCP 900/ Cardioserv SMC,SMS,SMK, SMU,SMV		3	DX-2355		DX3-90S	DX3-90P	2355S	2355P	D-1355	
			DX-2355-I		DX3-90S-I	DX3-90P-I	2355S-I	2355P-I	D-1355-I	
		Price	C6		B7	B7	C4	C4	C4	
		5	DX-2555		DX5-90S	DX5-90P	2555S	2555P	D-1555	
			DX-2555-I		DX5-90S-I	DX5-90P-I	2555S-I	2555P-I	D-1555-I	
		Price	C6		C2	C2	C6	C6	C6	
Tram 7200		5	MQ-2584		MQ5-90S	MQ5-90P	2584S	2584P	D-1584	
			MQ-2584-I		MQ5-90S-I	MQ5-90P-I	2584S-I	2584P-I	D-1584-I	
		Price	C6		C2	C2	C6	C6	C6	
GE/Marquette										
Eagle,Solar, Dash monitor MAC-Lab Cath Lab Datex-Ohmeda M-Prestn & “E” modules	 	3	MQ-2586		MQ3-90S	MQ3-90P	2386S	2386P	D-1386	
			MQ-2586-I		MQ3-90S-I	MQ3-90P-I	2386S-I	2386P-I	D-1386-I	
		Price	C6		C2	C2	C4	C4	C4	
		5	MQ-2586		MQ5-90S	MQ5-90P	2586S	2586P	D-1586	
			MQ-2586-I		MQ5-90S-I	MQ5-90P-I	2586S-I	2586P-I	D-1586-I	
		Price	C6		C2	C2	C6	C6	C6	
			6	MQ-2686		MQ6-90S	MQ6-90P	2686S	2686P	NA
				MQ-2686-I		MQ6-90S-I	MQ6-90P-I	2686S-I	2686P-I	
			Price	C8,2		C4	C4	D3	D3	
			12	MQ-21086		MQ5-90S MV5-90S	MQ5-90P MV5-90P	21086S	21086P	NA
				MQ-21086-I		MQ5-90S-I MV5-90S-I	MQ5-90P-I MV5-90P-I	21086S-I	21086P-I	
			Price	D3		2XC2	2XC2	D3	D3	
		3	MQA-2382		MQA3-90S	MQA3-90P	2382S	2382P	NA	
			MQA-2382-I		MQA3-90S-I	MQA3-90P-I	2382S-I	2382P-I		
		Price	C6		C2	C2	D1	D1		
		3	MQB-2381		MQB3-90S	MQB3-90P	2381S	2381P	NA	
			MQB-2381-I		MQB3-90S-I	MQB3-90P-I	2381S-I	2381P-I		
		Price	C6		C2	C2	D1	D1		
(*ESU-Proof)		3	MQ-2586EP		MQ3-90S	MQ-90P	2386ES	2386EP	D-1386EP	
					2386ES-I	2386EP-I	D-1386EP-I			
		Price	C6		C2	C2	D1	D1	D1	
Corometrics120 Patient Maternal Fetal		5	MQ-2586EP		MQ5-90S-I	MQ5-90P-I	2586ES	2586EP	D-1586EP	
					2586ES-I	2586EP-I	D-1586EP-I			
		Price	C6		C2	C2	D3	D3	D3	
		MECG	ADE-001-MF		2319S	2319P	NA	NA	NA	
		FECG	1552AAO		2319S-I	2319P-I				
		Price	C8		C4	C4				
		3	NA		NA	NA	2319S 1442AAO	2319P	D-1319	
							2319S-I	2319P-I	D-1319-I	
		Price					C4	C4	C4	
		3	MQ-2319 1554AAO		MQ3-90S	MQ3-90P	2319S	2319P	D-1319	
			MQ-2319I		MQ3-90S-I	MQ3-90P-I	2319S-I	2319P-I	D-1319-I	
		Price	C4		C2	C2	C4	C4	C4	
GE/Dash, Dinamap		5	AP-2586		AP5-90S	AP5-90P	2586S	2586P	D-1586	
			AP-2586-I		AP5-90S-I	AP5-90P-I	2586S-I	2586S-I	D-1586-I	
		Price	C5		C7	C7	C6	C6	C6	


OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable
Innocare									
IM200:36281		5	NA		NA	NA	2550S	2550P	NA
							2550S-I	2550P-I	
		Price					C2	C2	
Kontron									
Micromon, Minimon Vitalmon, Supermon Kolormon		3	NA		NA	NA	2398S	2398P	D-1398
							2398S-IH	2398P-IH	D-1398-IH
		Price					C4	C4	C4
105, 107, 108, 125, 128, VSM510. Old generation		3	NA		NA	NA	2393S	2393P	D-1393
							2393S-I	2393P-IH	D-1393-IH
		Price					C4	C4	C4
Lohmerier									
M010,M011/M211 Defi 501		3	LM-2358		LM3-90S	LM3-90P	2358S	2358P	D-1358
			LM-2358-I		LM3-90S-I	LM3-90P-I	2358S-I	2358P-I	D-1358-I
		Price	C6		B7	B7	C4	C4	C4
		5	LM-2558		LM5-90S	LM5-90P	NA	NA	NA
			LM-2558-I		LM5-90S-I	LM5-90P-I			
		Price	NA		C2	C2			
Medtronic Physio Control									
Lifepak 11		3	NA		NA	NA	2327S 3006218-02	2327P 1111P-000029	D-1380R
Lifepak 12		Price					2327S-I	2327P-I	D-1380R-I
Lifepak 15		5	NA		NA	NA	C4	C4	NA
Lifepak 20		Price					2527S 1111P-000066	2527P-I	
		12	MD-21027		MDA4-120S MDB6-90S MDA4-120S-I MDB6-90S-I C2/C6	MDA4-120P MDB6-90P MDA4-120P-I MDB6-90P-I C2/C6	21027S-I	21027P-I	NA
		Price	MD-21027-I				21027S-I	21027P-I	NA
Lifepak 11		4	NA		NA	NA	24108S 1111P-000018	24108P	NA
Lifepak 12		Price					24108S-I	24108P-I	
Lifepak 15		6	NA		NA	NA	26107S 1111P-000022	26107P	NA
Lifepak 20		Price					26107S-I	26107P-I	
							C6	C6	
MEDIANA									
ME306LK		3	AA-23104		AA3-90S	AA3-90P	23104S	23104P	D-13104
			AA-23104-I		AA3-90S-I	AA3-90P-I	23104S-I	23104P-I	D-13104-I
		Price	C4		B7	B7	C4	C4	C4

OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable	
MEK										
MP500, MP600, MP1000		3	NA		NA	NA	2310S	2310P	D-1310	
		Price				2310S-I	2310P-I	D-1310-I		
		5	NA		NA	NA	2510S	2510P	D-1510	
		Price				2510S-I	2510P-I	D-1510-I		
		3	NA		NA	NA	2311S	2311P	D-1311	
		Price				2311S-I	2311P-I	D-1311-I		
		5	NA		NA	NA	2511S	2511P	D-1511	
		Price				2511S-I	2511P-I	D-1511-I		
Mennen										
Series 9000, 1100,2000		3	NA		NA	NA	2374S	2374P	D-1374	
		Price				2374S-I	2374P-I	D-1374-I		
		5	NA		NA	NA	2574S	2574P	D-1574	
		Price				2574S-I	2574P-I	D-1574-I		
Envoy		3	AA-2375		MN3-90S	MN3-90P	2375S	2375P	D-1375	
		Price				MN3-90S-I	MN3-90P-I	2375S-I	2375P-I	D-1375-I
		5	AA-2575			B7	B7	C4	C4	C4
		Price					MN5-90S	MN5-90P	2575S	2575P
			AA-2575-I		MN5-90S	MN5-90P-I	2575S-I	2575P-I	D-1575-I	
			C6		C2	C2	C6	C6	C6	
Mindray										
PM7000/8000/9000	See AAMI									
PM5000 PM6000 BeneView T5/T6/T8		3	MR-2512		MR3-90S	MR3-90P	2312S	2312P	D-1312	
		Price				MR3-90S-I	MR3-90P-I	2312S-I	2312P-I	D-1312-I
		5	MR-2512		B7	B7	C4	C4	C4	
		Price				MR5-90S	MR5-90P	2512S	2512P	D-1512
			MR-2512-I		MR5-90S-I	MR5-90P-I	2512S-I	2512P-I	D-1512-I	
			C6		C2	C2	C6	C6	C6	
DPM4, DPM5 MEC 200, MEC 1000, MEC 1200, MEC 2000, PM 7000, PM 8000, PM 9000		3	MR-2340 0010-30-43126		MRB3-90S 0010-30-43251	MRB3-90P 0010-30-43252	2340S 0010-30-43117	2340P	D-1340R 040-000753-01	
		Price				MRB3-90S-I	MRB3-90P-I	2340S-I	2340P-I	D-1340R-I
		5	MR-2540 0010-30-42782		C1	C1	C1	C1	C1	
		Price				MR5-90S 0010-30-43145	MR5-90P 0010-30-43145	2540S 0010-30-43116	2540P	D-1540R
			MR-2540-I		MR5-90S-I	MR5-90P-I	2540S-I	2540P-I	D-1540R-I	
			C4		C2	C2	C2	C2	C2	
DPM6 DPM7 Passport 8 Passport 12 T1		3	MR-2512EP 0010-30-43140 MR-2512EP-I 0010-30-42724		MR3-90S 0010-30-43240	MR3-90P 0010-30-43240		NA	NA	
		Price				MR3-90S-I	MR3-90P-I		NA	NA
		5	MR-2512EP 0010-30-43140 MR-2512EP-I 0010-30-42724		B8	B8		NA	NA	
		Price				MR5-90S 0010-30-43145	MR5-90P 0010-30-43144		NA	NA
					MR5-90S-I	MR5-90P-I				
			D3		C2	C2				
BeneVision, BeneView, IPM, IMEC Series Monitor, Beneheart fibrillator		10	MR-21012		MRL5-90S MRP5-90S	MRL5-90P MRP5-90P		NA	NA	
		Price				MRL5-90S-I MRP5-90S-I	MRL5-90P-I MRP5-90P-I		NA	NA
Mindray										


OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable	
Mindray										
Passport 2 Spectrum Spectrun OR (*ESU-Proof)		3	MR-2540EP MR-2540EP-I		MR3-90S 0010-30-43251	MR3-90P 0010-30-43251	NA	NA	NA	
		Price	C7		MR3-90S-I	MR3-90P-I				
		5	MR-2540EP MR-2540EP-I		MR5-90S 0010-30-43251	MR5-90P 0010-30-43251	NA	NA	NA	
		Price	D1		MR5-90S-I	MR5-90P-I				
NEC										
YCE205		3	NA		NECA3-90S	NECA3-90P	NA	NA	NA	
		Price			NECA3-90S-I	NECA3-90P-I				
47502		3	NA		NECB3-90S	NECB3-90P				
		Price			NECB3-90S-I	NECB3-90P-I				
47504		3	NA		NECC3-90S	NECC3-90P	NA	NA	NA	
		Price	NECC3-90S-I		NECC3-90P-I					
			B7		B7					
Nihon Kohden										
Def. Cardiolife Lifescope serieos 6-8-9-11 6102-6205-6120		3	NK-2508 NK-2508-I		NKA3-90S	NKA3-90P	2308S	2308P	D-1308	
		Price	C6		NKA3-90S-I	NKA3-90P-I	2308S-I	2308P-I	D-1308-I	
		5	NK-2508 NK-2508-I		NKA5-90S	NKA5-90P	2508S	2508P	D-1508	
		Price	C6	C2	C2	C6	C6	C6		
			3	NK-2507 NK-2507-I		NKA3-90S	NKA3-90P	2307S	2307P	D-1307
			Price	C6		NKA3-90S-I	NKA3-90P-I	2307S-I	2307P-I	D-1307-I
5	NK-2507 NK-2507-I	NKA5-90S	NKA5-90P	2507S		2507P	D-1507			
TEC-7621C/TEC-7631C TEC-7721C/TEC-7731C 7700		3	NA		NA	NA	2313S	2313P		
		Price					2313S-I	2313P-I		
			C6		C6					
BSM-2300 lifescope I BSM-4100 lifescope P BSM-5100 lifescope A BSM-9510 lifescope M BSM-9800 lifescope S		3	NK-2509 NK-2509-I		NKB3-90S	NKB3-90P	2309S	2309P	D-1309	
		Price	C8		NKB3-90S-I	NKB3-90P-I	2309S-1	2309P-I	D-1309-I	
		B7	B7		C6	C6	C6			
		Price	6	NK-2509 NK-2509-I	NKB6-90S	NKB6-90P	2609S	2609P		
			C8	C3	C3	C8	C8			
			NKB6-90S-I	NKB6-90P-I	2609S-1	2609P-I				
Life Scope: OPV-1500		3	(JC-103TA) (JC-103T)		NKB3-90S	NKB3-90P	2316S	2316P	D-1316	
		Price	NA		NKB3-90S-I	NKB3-90P-I	2316S-1	2316P-I	D-1316-I	
		3	NK-2540 NK-2540-I		NKA3-90S	NKA3-90P	2340S	2340P	D-1340	
		Price	C6		NKA3-90S-I	NKA3-90P-I	2340S-I	2340P-I	D-1340-I	
			B7		B7	C1	C1	C1		
		Price	5	NK-2540 NK-2540-I	NKA5-90S	NKA5-90P	2540S	2540P	D-1540	
			C6	NKA5-90S-I	NKA5-90P-I	2540S-I	2540P-I	D-1540-I		
			C2	C2	C2	C2	C2			

OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable
Philips/HP									
43100,43200A, 77020A/AC/AR Serie 78300, serie 78500 Serie 78600, M2400A,M2406A, M2424A		3	AA-23029		AA3-90S	AA3-90P	23029S	23029P	D-13029
		Price	AA-23029-I		AA3-90S-I	AA3-90P-I	23029S-I	23029P-I	D-13029-I
			C6		B7.9	B7.9	C4	C4	C4
		5	AA-25029		AA5-90S	AA5-90P	25029S	25029P	D-15029
		Price	AA-25029-I		AA5-90S-I	AA5-90P-I	25029S-I	25029P-I	D-15029-I
			C6		C2	C2	C6	C6	C6
Merlin 78352C/78354C/ 78834CM1001A/ M1002A, Virdia, Ominicare CMA, ML65A,66A Heartstream Codemaster XL ...		3	AA-2385		AA3-90S	AA3-90P	2385S	2385P	D-1385
		Price	AA-2385-I		AA3-90S-I	AA3-90P-I	2385S-I	2385P-I	D-1385-I
			C6		B7.9	B7.9	C4	C4	C4
		5	AA-2585		AA5-90S	AA5-90P	2585S	2585P	D-1585
		Price	AA-2585-I		AA5-90S-I	AA5-90P-I	2585S-I	2585P-I	D-1585-I
			C6		C2	C2	C6	C6	C6
Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863064, 863065, 863066, 863068, 863087, 863088		3	HP-2385		HP3-90S	HP3-90P	2385S	2385P	D-1385
		Price	HP-2385-I		HP3-90S-I	HP3-90P-I	2385S-I	2385P-I	D-1385-I
			C7		C2	C2	C4	C4	C4
		5	HP-2585		HPA5-90S	HPA5-90P	2585S	2585P	D-1585
		Price	HP-2585-I		HPA5-90S-I	HPA5-90P-I	2585S-I	2585P-I	D-1585-I
			C8		C4	C4	C6	C6	C6
Multi-Measurement Modules: M3000A, M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18; Patient Monitors: MP5: M8105A Opt C12, MP2: M8102A Opt C12, ECG Modules: M1001B, M1002B, 860335, 860338		10	HP-21085		HPA5-90S (blue)	HPA5-90P (Blue)	21085S	21085P	
					HPB5-90S (Red)	HPB5-90P (Red)			
					HPA5-90S-I (Blue)	HPA5-90P-I (Blue)	21085S-I	21085P-I	
		Price	HP-21085-I		HPB5-90S-I (Red)	HPB5-90P-I (Red)	21085S-I	21085P-I	
Physio Control									
Lifepack5,6,7,8,9,9P, VSM2,1440,Vingmed 5	See GE Critikon								
LIFEPAK12- LIFEPAK20E- LIFEPAK15		3					2327S	2327P	
		Price					2327S-I	2327P-I	
							C4	C4	
		5					2527S	2527P	
With resistors		Price					2527S-I	2527P-I	
							C6	C6	
		12	MD-21027		MDA4-120S MDB6-90S		21027S	21027P-I	
			MD-21027-I		MDA4-120S-I MDB6-90S-I		21027S-I	21027P-I	
Price	D3		D3		D3	D3			


OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable		
Schiller Siemens											
Sirecust series 400 Sirecust 341,610/ 620/630/720/730/ 732/900/960/961 1260/1281/1801, 400-500 ECG		3	SM-2376		SM3-90S	SM3-90P	2376S	2376P	D-1376		
		Price	SM-2376-I		SM3-90S-I	SM3-90P-I	2376S-I	2376P-I	D-1376-I		
			C6		B7	B7	C4	C4	C4		
		5	SM-2576		SM5-90S	SM5-90P	2576S	2576P	D-1576		
			SM-2576-I		SM5-90S-I	SM5-90P-I	2576S-I	2576P-I	D-1576-I		
		Price	C6		C2	C2	C6	C6	C6		
SC5000,SC6000, SC7000 SC8000,SC9000 Multi-Parameter Cable SM-STE5		3	SM-STE5	See the following page	SMB3-90S	SMB3-90P					
			SM-STE5		SMB3-90S-I	SMB3-90P-I					
		Price	D8		B7	B7					
			5		SM-STE5	SM5-90S	SM5-90P				
		SM-STE5			SM5-90S-I	SM5-90P-I					
		Price	D8		C2	C2					
		3	SM-2323			SM3-90S	SM3-90P	2323S	2323P	D-1323	
		Price	SM-2323-I			SM3-90S-I	SM3-90P-I	2323S-I	2323P-I	D-1323-I	
	Defigard4000 defibrillator		3					2354S	2354P	D-1354R	
			Price					2354S-I	2354P-I	D-1354R-I	
						2354S-IH	2354P-IH				
Sirecust series 400 Sirecust 341,610 620/630/720/730 732/900/960/961 1260/1281/1801,		3	SM-2376		SM3-90S	SM3-90P	2376S	2376P	D-1376		
		Price	SM-2376-I		SM3-90S-I	SM3-90P-I	2376S-I	2376P-I	D-1376-I		
			C6		B7	B7	C4	C4	C4		
		5	SM-2576		SM5-90S	SM5-90P	2576S	2576P	D-1576		
			SM-2576-I		SM5-90S-I	SM5-90P-I	2576S-I	2576P-I	D-1576-I		
		Price	C6		C2	C2	C6	C6	C6		
			Spacelabs								
		90701,90303, 90305,90308, 90309,90385	See AAMI 								
Ultraview 90492, 91220, 90496, 91496		3	SL-2396		SL3-90S	SL3-90P	2396S	2396P	D-1396		
			SL-2396-I		SL3-90S-I	SL3-90P-I	2396S-I	2396P-I	D-1396-I		
		Price	C6		B7	B7	C4	C4	C4		
		5	SL-2596		SL5-90S	SL5-90P	2596S	2596P	D-1596		
			SL-2596-I		SL5-90S-I	SL5-90P-I	2596S-I	2596P-I	D-1596-I		
		Price	C6		C2	C2	C6	C6	C6		
		12	MQ-21096		MQ5-90S MV5-90S	MQ5-90P MV5-90P	21096S	21096P			
			MQ-21096-I		MQ5-90S-I MV5-90S-I	MQ5-90P-I MV5-90P-I	21096S-I	21096P-I			
			Price		D3	2XC2	2XC2	D3	D3		
Spacelabs											

OEM & Models	Monitor Plug	Lead No.	Trunk Cable	Lead yoke Conn.	Lead wire/w Snap	Lead wire/w Grabber	One piece/w Snap	One piece/w Grabber	Din Trunk Cable
Spacelabs									
Ultraview 90492, 91220, 90496, 91496		3	HPC-2596 HPC-2596-I		HPC3-90S HPC3-90S-I	HPC3-90P HPC3-90P-I	2396S 2396S-I	2396P 2396P-I	D-1396 D-1396-I
		Price	C8		C2	C2	C4	C4	C4
		5	HPC-2596 HPC-2596-I		HPC5-90S HPC5-90S-I	HPC5-90P HPC5-90P-I	2596S 2596S-I	2596P 2596P-I	D-1596 D-1596-I
		Price	C8		C4	C4	C6	C6	C6
Welch allyn									
Protocol	See AAMI								
Micropaq Propaq LT		3	D-1353 D-1353-I		D3-90S D3-90S-I	D3-90P D3-90P-I	2353S 2353S-I	2353P 2353P-I	
		Price	C7		A8	B3	C1	C1	
		5	D-1553 D-1553-I		D5-90S D5-90S-I	D5-90P D5-90P-I	2553S 2553S-I	2553P 2553P-I	
		Price	C7		B4	B6	C2	C2	
Welch Allyn 1500		10	NA		NA	NA	21079S 21079S-I	21079P 21079P-I	
		Price					D3	D3	
Zoll									
R SERIES E SERIES M SERIES M CCT SERIES		5	NA		NA	NA	21069S 21069S-I	21069P 21069P-I	
		Price					D3	D3	

One piece cables with neonate electrode end for most brands are available upon request.
For neonete application, Din style neonate leads are also the alternatives.



1.2 Major types or brandsECG lead wires

Types of Major ECG lead wires	Lead No.	REF# Snap(AHA/IEC)	OEM REF#	Price code	REF# Grabber(AHA/IEC)	OEM REF#	Price code
Din unshielded Individual set leads							
	3	D3-90S D3-90S-I		A8	D3-90P D3-90P-I		B3
	5	D5-90S D5-90S-I		B4	D5-90P D5-90P-I		B6
	7	D7-90S/D7-90S-I		B6	D7-90P/D7-90P-I		B7
	10	D10-90S/D10-90S-I		B8	D10-90P/D10-90P-I		C2
Din unshielded individual set (neonate)							
	3				D3-90NP D3-90NP-I (Mini Grabber)		B4
	5				D5-90NP D5-90NP-I (Mini Grabber)		B6
LL style(shielded)							
	3	LL3-90S LL3-90S-I		B7	LL3-90P LL3-90P-I		B7
	5	LL5-90S LL5-90S-I		C2	LL5-90P LL5-90P-I		C2
Mindray/Datascope lead							
	3	MR3-90S MR3-90S-I		C1	MR3-90P MR3-90P-I		C1
	5	MR5-90S MR5-90S-I	012-00-1503-02	C2	MR5-90P MR5-90P-I		C2
	3	DT3-90S DT3-90S-I	0012-00-1503-04/5/6 45cm,61cm,100cm	C2	DT3-90P DT3-90P-I	0012-00-1514-04/5/6 45cm,61cm,100cm	C2
	5	DT5-90S DT5-90S-I	0012-00-1503-01/2/3 45cm,61cm,100cm	C2	DT5-90P DT5-90P-I	0012-00-1514-01/2/3 45cm,61cm,100cm	C2
Datex-GE Integrated leads (DX300 series)							
	3	DX3-90S DX3-90S-I	545327	B7	DX3-90P DX3-90P-I	545317 545315	B7
	5	DX5-90S DX5-90S-I	545328	C2	DX5-90P DX5-90P-I	545318 545316	C2