



Descriptions	Probe Connector	Lead #	REF# AHA (Snap/Grabber)	REF# IEC (Snap/Grabber)	Package	Unit price code
GE Disposable Radiolucent Lead Wires,Vs Type/0.9m		3	MQD3-90XDS	MQD3-90XDS-I	1 PC	B4.1
		3	MQD3-90XDP	MQD3-90XDP-I	1 PC	B2
		5	MQD5-90XDS	MQD5-90XDS-I	1 PC	B5.3
		5	MQD5-90XDP	MQD5-90XDP-I	1 PC	B4
		10	MQD10-90XDS	MQD10-90XDS-I	1 PC	C2
		10	MQD10-90XDP	MQD10-90XDP-I	1 PC	B8
Spacelabs Disposable Radiolucent Lead Wires,LL Type/0.9m		3	LL3-90XDS	LL3-90XDS-I	1 PC	B2
		3	LL3-90XDP	LL3-90XDP-I	1 PC	B1
		5	LL5-90XDS	LL5-90XDS-I	1 PC	B5
		5	LL5-90XDP	LL5-90XDP-I	1 PC	B4
		10	LL10-90XDS	LL10-90XDS-I	1 PC	C1
		10	LL10-90XDP	LL10-90XDP-I	1 PC	B6.7
Spacelabs Disposable Radiolucent Lead Wires,SL Type/0.9m		3	SLB3-90XDS	SLB3-90XDS-I	1 PC	B5
		3	SLB3-90XDP	SLB3-90XDP-I	1 PC	B4
		5	SLB5-90XDS	SLB5-90XDS-I	1 PC	B7.9
		5	SLB5-90XDP	SLB5-90XDP-I	1 PC	B6.6
		10	SLB10-90XDS	SLB10-90XDS-I	1 PC	C3.8
		10	SLB10-90XDP	SLB10-90XDP-I	1 PC	C1
Philips Disposable Radiolucent Lead Wires,AA Type/0.9m		3	AAC3-90XDS	AAC3-90XDS-I	1 PC	B4
		3	AAC3-90XDP	AAC3-90XDP-I	1 PC	B2
		5	AAC5-90XDS	AAC5-90XDS-I	1 PC	B6
		5	AAC5-90XDP	AAC5-90XDP-I	1 PC	B5
		10	AAC10-90XDS	AAC10-90XDS-I	1 PC	C3.8
		10	AAC10-90XDP	AAC10-90XDP-I	1 PC	C1
COVIDIEN Disposable Radiolucent Lead Wires,CDM Style/0.9M		3	CDMA3-90XDS	CDMA3-90XDS-I	1 PC	B2
		3	CDMA3-90XDP	CDMA3-90XDP-I	1 PC	B1
		5	CDMA5-90XDS	CDMA5-90XDS-I	1 PC	B4.1
		5	CDMA5-90XDP	DMA5-90XDP-I	1 PC	B3
COVIDIEN Disposable Radiolucent Lead Wires,CDM Style/1.8m		5	CDMA5-180XDS	CDMA5-180XDS-I	1 PC	B6
		5	CDMA5-180XDP	CDMA5-180XDP-I	1 PC	B4.1
COVIDIEN Disposable Radiolucent Lead Wires,CDM Style/3.1m(V2-V6,C2-C6)		5	CDMV5-310XDS	CDMV5-310XDS-I	1 PC	B6.6
		5	CDMV5-310XDP	CDMV5-310XDP-I	1 PC	B6



4 ECG Lead Adapter

The unique design of the ECG lead wires adapters enable the hospitals to use the one type lead wires to fit major brands patient monitoring system running in one hospital by attaching the proper ECG adapter to the trunk cable. No confusion on lead wires any more !

Specification

Cable material	ABS Jacket, Light gray housing
Connector &Yoke	Nylon plug & Golden plated pins
Link Tie	12cm length with trunk cable clip

Compatibility

Major brands ECG monitoring leads system

OEM Manufacturer	3-lead REF#	Trunk Cable end	Lead wire end	Price	5-lead REF#	Trunk Cable end	Lead wire end	Price code
Convert To Covidien lead Covidien								
Philips MX40 Telemetry					AD-PTS5-CDF5			C4
					AD-PTS5-CDFB5			C7
Philips					AD-HP-CDF5			B7
GE Telemetry					AD-AP- CDF5			B5.3
GE Marquette					AD-MQ5-CDF5			B6
					AD-MQ6-CDF6			B7
					AD-MQ-CDF5			B6
Artema/S&W					AD-SW-CDF5			B8
Convert To Din style lead wires								
Colin	AD-CL3-D3			B5	AD-CL5-D5			B5
GE Marquette	AD-MQ5-D3			B5	AD-MQ5-D5			B5
GE Datex-Ohmeda	AD-DX3-D3			B5	AD-DX5-D5			B5
Nihon Kohden	AD-NKA3-D3			B5	AD-NKA5-D5			B5
Nihon Kohden	AD-NKB3-D3			B5				
Philips medical	AD-AA3-D3			B5				
Philips medical	AD-HP3-D3			B5				
Siemens/Drager	AD-SM3-D3			B5	AD-SM5-D5			B5
Spacelabs	AD-SL3-D3			B5	AD-SL5-D5			B5


New Adapter



Application Figure



OEM Manufacturer	3-lead REF#	Trunk Cable end	Lead wire end	Price	5-lead REF#	Trunk Cable end	Lead wire end	Price code
Convert To VS/Siemens lead VS/Siemens								
Colin	AD-CL3-SM3			B5	AD-CL5-SM5			B5
GE Marquette	AD-MQ5-SM3			B5	AD-MQ5-SM5			B5
GE Datex-Ohmeda	AD-DX3-SM3			B5	AD-DX5-SM5			B5
Nihon Kohden	AD-NKA3-SM3			B5	AD-NKA5-SM5			B5
Nihon Kohden	AD-NKB3-SM3			B5				
Philips medical	AD-AA3-SM3			B5				
Philips medical	AD-HP3-SM3			B8				
Spacelabs	AD-SL3-SM3			B5	AD-SL5-SM5			B5
Convert To Philips/HP style lead Philips/HP								
Colin	AD-CL3-HP3			B6	AD-CL5-HP5			B7
GE Marquette	AD-MQ5-HP3			B6	AD-MQ5-HP5			B7
GE Datex-Ohmeda	AD-DX3-HP3			B6	AD-DX5-HP5			B7
Nihon Kohden	AD-NKA3-HP3			B6	AD-NKA5-HP5			B7
Nihon Kohden	AD-NKB3-HP3			B6				
Philips medical	AD-AA3-HP3			B6				
Siemens/Drager	AD-SM3-HP3			B6	AD-SM5-HP5			B7
Spacelabs	AD-SL3-HP3			B6	AD-SL5-HP5			B7
Convert To GE Multi-Link lead GE Multi-Link								
Colin					AD-CL5-MQ5			B5
GE Datex-Ohmeda					AD-DX5-MQ5			B5
Nihon Kohden					AD-NKA5-MQ5			B5
Siemens/Drager					AD-SM5-MQ5			B5
Spacelabs					AD-SL5-MQ5			B5



Datex 300 style lead



Spacelabs style lead

OEM Manufacturer	3-lead REF#	Trunk Cable end	Lead wire end	Price	5-lead REF#	Trunk Cable end	Lead wire end	Price code
Convert To Datex 300 style lead								
Colin	AD-CL3-DX3			B5	AD-CL5-DX5			B5
GE Marquette	AD-CL3-DX3			B5	AD-MQ5-DX5			B5
Nihon Kohden	AD-NKA3-DX3			B5	AD-NKA5-DX5			B5
Nihon Kohden	AD-NKB3-DX3			B5				
Philips medical	AD-AA3-DX3			B5				
Philips medical	AD-HP3-DX3			B6				
Siemens/Drager	AD-SM3-DX3			B5	AD-SM5-DX5			B5
Spacelabs	AD-SL3-DX3			B5	AD-SL5-DX5			B5
Convert To Spacelabs style lead								
Colin	AD-CL3-SL3			B5	AD-CL5-SL5			B5
GE Marquette	AD-MQ5- SL3			B5	AD-MQ5-SL5			B5
GE Datex-Ohmeda	AD-DX3- SL3			B5	AD-DX5-SL5			B5
Nihon Kohden	AD-NKA3- SL3			B5	AD-NKA5-SL5			B5
Nihon Kohden	AD-NKB3- SL3			B5				
Philips medical	AD-AA3- SL3			B5				
Philips medical	AD-HP3- SL3			B6				
Spacelabs	AD-SM3- SL3			B5	AD-SM5-SL5			B5



5 Diagnosis EKG Cables & Accessories

Specification

Cable material	TPU Jacket, Light gray
Connector &Yoke	Nylon plug & Golden plated pins
Trunk Cable length	2m
Lead wire length	1.2m(limb lead);0.8m(Chest lead)

Compatibility

Major brands EKG diagnosis system