


Application Diagram




Product name	Photo	REF# Example	Reference to REF# Descriptions
One piece Fixed Cable		E10R-HP-B/I	E: EKG diagnostic cable or lead wire 10: 10-lead(12=12 lead) R: built in resisters for protection(default: no resisters) HP: Philips/HP,factory/model code 1:Connector configure; B: banana end(N= needle, P=grabber, S=snap) I: IEC code (Default is AHA code)
Trunk Cable		E10R-MQ-I	E: EKG diagnostic cable or lead wire 10: 10-lead(12=12 lead) R: built in resisters for protection(default: no resisters) MQ: GE Marquette, factory/model code 1:Connector configure (default: only one connector available) I: IEC code (Default is AHA code)
Lead Wire		MQ10-LB-I	MQ: Marquette, factory/model code 10: 10-lead(12=12 lead) L: lead wires B: banana end(N= needle, P=grabber, S=snap) I: IEC code (Default is AHA code)










All one piece cable or lead wires are available with banana end, needle end, grabber end and snap end

Code	AHA(A)	Standard IEC(I)	Alternate IEC(IH)
4mm Banana end			
3mm Needle end			
Grabber End			
Snap End			

5.1 One Piece EKG Fixed Cable



OEM & Models	EKG machine Plug	Lead#	With 4mm banana	With 3mm needle	With grabber end	With Snap end	Price code
Burdick/Cardiac Science							
EK-10 and Elite, E350, E350i, E550, E560, Eclipse™ LE/LE.		10	E10-BK1-B E10R-BK1-B	E10-BK1-N E10R-BK1-N	E10-BK1-P E10R-BK1-P	E10-BK1-S E10R-BK1-S	D3
		10	E10-BK1-B/I E10R-BK1-B/I	E10-BK1-N/I E10R-BK1-N/I	E10-BK1-P/I E10R-BK1-P/I	E10-BK1-S/I E10R-BK1-S/I	D3
Burdick Quinton Eclipse 4,8,400,800, 850,850i,plus, Eclipse LE II and ATRIA PLUS Quinton		10	E10-BK2-B E10R-BK2-B	E10-BK2-N E10R-BK2-N	E10-BK2-P E10R-BK2-P	E10-BK2-S E10R-BK2-S	D3
		10	E10-BK2-B/I E10R-BK2-B/I	E10-BK2-N/I E10R-BK2-N/I	E10-BK2-P/I E10R-BK2-P/I	E10-BK2-S/I E10R-BK2-S/I	D3
Edan							
SE-18,SE-12Express SE-1200Express , SE-1201SE-601,SE-3, SE-301,SE-100,DP12		10	E10-ED-B E10R-ED-B	E10-ED-N E10R-ED-N	E10-ED-P E10R-ED-P	E10-ED-S E10R-ED-S	D3
		10	E10-ED-B/I E10R-ED-B/I	E10-ED-N/I E10R-ED-N/I	E10-ED-P/I E10R-ED-P/I	E10-ED-S/I E10R-ED-S/I	D3
Fukuda Denshi							
FX101, FX102F, FX328U, FCP2201		10	E10-FD1-B E10R-FD1-B	E10-FD1-N E10R-FD1-N	E10-FD1-P E10R-FD1-P	E10-FD1-S E10R-FD1-S	D3
		10	E10-FD1-B/I E10R-FD1-B/I	E10-FD1-N/I E10R-FD1-N/I	E10-FD1-P/I E10R-FD1-P/I	E10-FD1-S/I E10R-FD1-S/I	D3
Fukuda ME							
KP-500, KP-500D		10	E10-FD2-B E10R-FD2-B	E10-FD2-N E10R-FD2-N	E10-FD2-P E10R-FD2-P	E10-FD2-S E10R-FD2-S	D3
		10	E10-FD2-B/I E10R-FD2-B/I	E10-FD2-N/I E10R-FD2-N/I	E10-FD2-P/I E10R-FD2-P/I	E10-FD2-S/I E10R-FD2-S/I	D3

OEM & Models	EKG machine Plug	Lead#	With 4mm banana	With 3mm needle	With grabber end	With Snap end	Price code
GE-Marquette							
MAC 500, MAC 1100, MAC 1200		10	E10-MQ-B E10R-MQ-B	E10-MQ-N E10R-MQ-N	E10-MQ-P E10R-MQ-P	E10-MQ-S E10R-MQ-S	D3
		10	E10-MQ-B/I E10R-MQ-B/I	E10-MQ-N/I E10R-MQ-N/I	E10-MQ-P/I E10R-MQ-P/I	E10-MQ-S/I E10R-MQ-S/I	D3
Hellige / Siemens Hormann / Bosch							
MEGO COMP 60PLUS DEGO Cardiomid 3 Van Begg Biomat 3000 / 3600 / 3700/Biomat 6000 / 8000 Hellige: EK36/ EK41/ EK43 / EK 53/EK56/ EK403, FK413/EK432/EK433/EK433P, Cardiomart, CardioSys, Cora Siemens: Cardiotest 1 / 3 T / 703 / 710, Cardina C, Cardires 370A/3702, Cathcor, MegacardElema/B, R/R/E, Mingodisc Mingocard 3 / 4 / 7.		10	E10-HL-B E10R-HL-B	E10-HL-N E10R-HL-N	E10-HL-P E10R-HL-P	E10-HL-S E10R-HL-S	D3
		10	E10-HL-B/I E10R-HL-B/I	E10-HL-N/I E10R-HL-N/I	E10-HL-P/I E10R-HL-P/I	E10-HL-S/I E10R-HL-S/I	D3
Kenz							
103,106 CARDIOLINE (Italy)		10	E10-KZ1-B E10R-KZ1-B	E10-KZ1-N E10R-KZ1-N	E10-KZ1-P E10R-KZ1-P	E10-KZ1-S E10R-KZ1-S	D3
		10	E10-KZ1-B/I E10R-KZ1-B/I	E10-KZ1-N/I E10R-KZ1-N/I	E10-KZ1-P/I E10R-KZ1-P/I	E10-KZ1-S/I E10R-KZ1-S/I	D3
ECG 108/110 Cardico 302/601/ 1210/1211		10	E10R-KZ2-B	E10R-KZ-N	E10R-KZ2-P	E10R-KZ2-S	D3
		10	E10R-KZ2-B/I	E10R-KZ2-N/I	E10R-KZ2-P/I	E10R-KZ2-S/I	D3
Mortara							
ELI-200, ELI-50		10	E10-MT-B E10R-MT-B	E10-MT-N E10R-MT-N	E10-MT-P E10R-MT-P	E10-MT-S E10R-MT-S	D3
		10	E10-MT-B/I E10R-MT-B/I	E10-MT-N/I E10R-MT-N/I	E10-MT-P/I E10R-MT-P/I	E10-MT-S/I E10R-MT-S/I	D3
Nihon Kohden							
CaduoFax 6151, 6851,6353 CadioFax 8110		10	E10-NK1-B E10R-NK1-B	E10-NK1-N E10R-NK1-N	E10-NK1-P E10R-NK1-P	E10-NK1-S E10R-NK1-S	D3
		10	E10-NK1-B/I E10R-NK1-B/I	E10-NK1-N/I E10R-NK1-N/I	E10-NK1-P/I E10R-NK1-P/I	E10-NK1-S/I E10R-NK1-S/I	D3
CaduoFax Q ECG-9110K ECG-9130K Doing Juang Cardioline (New)		10	E10-NK2-B E10R-NK2-B	E10-NK2-N E10R-NK2-N	E10-NK2-P E10R-NK2-P	E10-NK2-S E10R-NK2-S	D3
		10	E10-NK2-B/I E10R-NK2-B/I	E10-NK2-N/I E10R-NK2-N/I	E10-NK2-P/I E10R-NK2-P/I	E10-NK2-S/I E10R-NK2-S/I	D3
Philips/HP							
M1770A,M1771A, M1772A Pagewrite100/200/ 200i		10	E10-HP-B E10R-HP-B	E10-HP-N E10R-HP-N	E10-HP-P E10R-HP-P	E10-HP-S E10R-HP-S	D3
		10	E10-HP-B/I E10R-HP-B/I	E10-HP-N/I E10R-HP-N/I	E10-HP-P/I E10R-HP-P/I	E10-HP-S/I E10R-HP-S/I	D3
Schiller							
AT1,AT2,AT2 Plus, AT4 Resistor 10K Long Screw		10	E10-SH1-B E10R-SH1-B	E10-SH1-N E10R-SH1-N	E10-SH1-P E10R-SH1-P	E10-SH1-S E10R-SH1-S	D3
		10	E10-SH1-B/I E10R-SH1-B/I	E10-SH1-N/I E10R-SH1-N/I	E10-SH1-P/I E10R-SH1-P/I	E10-SH1-S/I E10R-SH1-S/I	D3
AT3,AT6,CS6,AT5, AT10,AT60 Short Screw		10	E10-SH2-B E10R-SH2-B	E10-SH2-N E10R-SH2-N	E10-SH2-P E10R-SH2-P	E10-SH2-S E10R-SH2-S	D3
		10	E10-SH2-B/I E10R-SH2-B/I	E10-SH2-N/I E10R-SH2-N/I	E10-SH2-P/I E10R-SH2-P/I	E10-SH2-S/I E10R-SH2-S/I	D3