

Finger Tip Pulse Oximeters

UC-60B



Features

- It is lightweight, small in size and easy to take;
- Accurately measures SpO2 value and Pulse Rate value;
- Large bright LCD display of SpO2, PR and Pulse Bar;
- Innovative 2 directions display;
- Automatic Power Off;
- More than 50 hours continuous working time with 2 AAA alkaline batteries;
- Low battery voltage indication

Finger Tip Pulse Oximeters

UC-60B1



Features

- It is lightweight, small in size and easy to take;
- Accurately measures SpO2 value and Pulse Rate value;
- Dual color display of SpO2, PR, PI, Plethysmogram and Pluse Bar;
- Innovative 2 directions display;
- Automatic Power On/Off;
- More than 30 hours continuous working
- Audible & visible alarm function;