

Handheld Pulse Oximeters

UC-66B Handheld Pulse Oximeter



Features

- Light, small in size and easy to carry for long time monitoring;
- High resolution LCD to display SpO₂, PR, PI and plethysmogram;
- Adult/pediatric/neonate probes replaceable & optional;
- Stored data can be uploaded to PC for further review and analysis;
- Up to 70 hours of oximetry data storage (Record data once every 1/2/4/8 seconds);
- Audible & visual alarm function;
- 3AA alkaline or rechargeable batteries or AC/DC adapter (optional)

Handheld Pulse Oximeters

UC-68B Wrist Pulse Oximeter



Features

- Small, light and easy to wear for long time monitoring
- Bright color OLED display with plethysmogram
- Accurate measuring SpO₂, Pulse Rate and Perfusion Index
- Audible and visual alarm with alarm limits adjustable
- Real time clock display
- Up to 72 hours data storage with selectable record intervals (1s/4s/8s)
- Data record can be uploaded to PC for future review and analyse