

UC-66B Handheld VET Pulse Oximeter



Features

- Light, small in size and easy to carry for long time monitoring;
- High resolution LCD to display SpO₂, PR, PI and plethysmogram;
- Adult/pediatric/neonate probes replaceable & optional;
- Stored data can be uploaded to PC for further review and analysis;
- Up to 70 hours of oximetry data storage (Record data once every 1/2/4/8 seconds);
- Audible & visual alarm function;
- 3AA alkaline or rechargeable batteries or AC/DC adapter (optional)

JPD-100B Handheld Digital Fetal Doppler



Features

- Safe to use
- Compact & Portable
- Large blue backlit LCD display
- Built in speaker
- 2.5MHz waterproof probe
- Rechargeable batteries.

JPD-200C Portable Digital Fetal Doppler



Features

- Safe to use
- Portable
- LCD display
- Built in speaker
- 2.2MHz or 3.3MHz waterproof probes
- Rechargeable 14.4V Nickel Hydrogen battery