




Multi-Parameter Patient Cable



REF#	Used with the OEM and models	Descriptions	Cable plug	Package	Price code
SM-STE5					
(3368391), 2.5m	Siemens/Drager Infinity 6000 / 7000 / 9000	Multimed 5 Pod patient cable, 2.5m or 1.5m 5-lead ECG, VS/Siemens Spo2/Temp/Resp. AHA & IEC		1pcs	D8
SM15-STE5	Drager Infinity Gamma/Gamma XL/Gamma X XL				
(5950192), 1.5m	Delta/Delta XL Kappa/Kappa XLT Vista/Vista XL				
SM-STE6					
(5191221)	Siemens/Drager Infinity 6000 7000 9000 Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/Delta XL	Multimed 6 Pod patient cable, 2.5m 6-lead ECG, VS/Siemens Spo2/Temp/Resp. AHA & IEC		1pcs	D8.5.8
SM-STE12					
(5589663)	Drager Infinity Delta/Delta XL Kappa/Kappa XLT	Multimed 12 Pod patient cable, 3.0m 12-lead ECG, VS/Siemens Spo2/AUX. AHA & IEC		1pcs	NA
SM-MX5					
(5590539)	Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/Delta XL Kappa/Kappa XLT Vista/Vista XL	NeoMed Pod Patient cable, 3.2m ECG/Resp./Spo2 Dual Temp/FiO2 AHA & IEC		1pcs	E1
DG-STE6					
(MS20093)	Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/Delta XL Kappa/Kappa XLT Vista/Vista XL	Multimed Plus patient cable, 2.5m Spo2/Temp 6-lead ECG/Resp. AHA & IEC		1pcs	E1
DG-STE6-OR					
(MS20094)	Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/Delta XL Kappa/Kappa XLT Vista/Vista XL	Multimed Plus OR patient cable, 2.5m Spo2/Temp 6-lead ECG/Resp. AHA & IEC		1pcs	NA
UNCC20-008					
(CareFusion 8120)	CareFusion Monitor	CareFusion/LM4-XX/2.0M		1pcs	D3
SM-STE6					
Y05-T091	Unimed compatible multi-parameter patient cable like SM-STE5.	The holder of the multi-parameter cable yoke to fix onto the pole or wall on bedside.		1pcs	C2
Y05-T092	Siemens/Drager original multi-parameter patient cable. Siemens/Drager	The purpose is the Same as above		1pcs	C6