

Oxygen Flow Meters

Ball & Tube Type

#FM4546-S
Single




#FM4546-D
Double



Features

- Polycarbonate Inner & Outer Tubes
- Calibrated at 50 psi inlet pressure
- Chromed brass body.
- Flow rate: 0-15LPM
- Adapter: Heyer



Extended versions also available

#FM4546-SE
Single

#FM4546-DE
Double

Oxygen Flow Meters

Click Type

#FM4000-S
Single




#FM4000-D
Double



Features

- Calibrated at 50 psi inlet pressure
- Anodized Aluminium body.
- Flow rate: 0-25LPM
- Adapter: Heyer



Extended versions also available

#FM4000-SE
Single

#FM4000-DE
Double