



1 Reusable SpO2 Sensor


Specification

Sensor Connector	OEM compatible, Nylon adapter with gold plated pins
Cable material & color	TPU jacket; light gray or blue
Cable length & diameter	1.1m/1.5m/3.0m; ϕ3.20mm/4.00mm or 5.00x2.50mmtwo line cable
Sensor house	ABS shell(finger clip); Silicon rubber boots(soft and multi-site)

Compatible With Major Patient Monitor Brands

Nellcor
Ohmeda
Philips
Marquette
Spacelabs
Novamatrix

Criticare
BCI
Datascopes
Datex
Nihon Kohden
Palco
.....



Series系列	U400	U100	U400S	U100S	U200S	U300
REF to Fig	A	B	C	D	E	F
Name	Adult Finger	Pediatric Finger	Adult soft tip	Pediatric Soft tip	Infant Soft tip	Neonate Wrap
Patient size	+40kg	15-40kg	+40kg	15-40kg	3-15kg	1-3kg
Cable color	Gray Cable	Gray Cable	Blue Cable	Gray two line	Gray two line	Gray two line
Series系列	U800	U900	U600	U600B	U800E	U800V
REF to Fig	G	H	I	J	K	L
Name	Multi-site	Adult Ear clip	Adult/Neo. Wrap	Pedia/Inf. Wrap	Multi-site	Forehead/Rectal
Patient size	All patients	+40kg	+40kg or -3kg	15-40kg	All patients	Adult
Cable color	Gray two line		Gray two line	Gray two line	Sensormate Used with DPT tapes	Gray Cable

How to Identify the REF# of Unimed SpO₂ Sensors

Reference Number Rule Description

U403S-01XE	
U	Unimed brand
4	adult sensor (1: pediatric; 2: infant;3: neonate; 5: disposable; 6 wrap; 7: ext-cable; 8: multi-site;9: ear clip)
03	1.1m cable length (03: 1.1m;05:1.5m;08: 2.2m;10: 3.0m)
S	Soft sensor (B:U style shell; Default: shell sensor or neonate wrap, alternative wrap, multi-site)
01	Nellcor code (02:Ohmeda;03:Novametrix; 04: Spacelabs; 05:Criticare; 06: BCI...)
X	Oximax technology(M: M technology, D: digital technology; Default: standard
E	Sensormate sensor (Default: non-Sensormate sensor)

Cross Reference

Ending "X" can be replaced by "XP" or "P" for fully compatible with all oximax oximetry

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
3F Medical (China)								
IRIS 300/IRIS 400 Phoebe PMS8300/8310 PhoebeProPMS8300/ 8310PMS8110A			U410-189	Adult finger /3.0m	C7	No need		
			U410S-189	Adult soft /3.0m	C8			
			U810-189	Multi-site Y/3.0m	C7			
Adecon (China)								
MCF-21L, DK-8000D,DK-8000L			U410-120	Adult finger /3.0m	C7	No need		
			U410S-120	Adult soft /3.0m	C8			
			U810-120	Multi-site Y/3.0m	C7			
Anke (China)								
ASC-553			U410-37	Adult finger /3.0m	C7	No need		
			U410S-37	Adult soft成/3.0m	C8			
			U810-37	Multi-site Y/3.0m	C7			
ASC-553 (*Oximax tech)			U410-78P	Adult finger /3.0m	C8	No need		
			U410S-78P	Adult soft /3.0m	D1			
			U810-78P	Multi-site Y/3.0m	C8			
ASC-553A			U410-82	Adult finger /3.0m	C7	No need		
			U410S-82	Adult soft /3.0m	C8			
			U810-82	Multi-site Y/3.0m	C7			
ASC-553A (*Ohmeda Tech)			U410-107	Adult finger /3.0m	C7	No need		
			U410S-107	Adult soft /3.0m	C8			
			U810-107	Multi-site Y/3.0m	C7			


OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Artema/S&W								
Diascope,8000, Athena (*Non-Oximax)		DS100A	U403-01	Adult finger/1.1m	C3	U708-80		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-80	Adult finger /3.0m	C8	No need		
			U410S-80	Adult soft /3.0m	D1			
			U810-80	Multi-site Y/3.0m	C8			
S&W			U410-19	Adult finger /3.0m	C8	No need		
			U410S-19	Adult soft成 /3.0m	D1			
			U810-19	Multi-site Y/3.0m	C8			
Baxter/Simed								
ASAT, S-50/100/100a			U410-12	Adult finger /3.0m	C8	No need		
			U410S-12	Adult soft /3.0m	D1			
			U810-12	Multi-site Y/3.0m	C8			
Monarch/Explorer (*Oximax)		DS100A	U403-01	Adult finger/1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft/1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
BCI								
3300,3301,3302, 3303,3304,3401 (*BCI)		3044	U403-06	Adult finge /1.1m	C3	Depends U708-06		
			U103-06	Pediatric finger /1.1m	C5			
			U403S-06	Adult soft /1.1m	C3			
			U103S-06	Pediatric soft /1.1m	C5			
			U203S-06	Infant soft /1.1m	C5			
			U303-06	Neonate wrap /1.1m	C5			
		3043	U803-06	Multi-site Y/1.1m	C3			
			U603-06	Adult/Neo Wrap/1.1m	C3			
			U603B-06	Pedia/infant wrap/1.1m	C3			
		3078	U903-06	Adult ear clip /1.1m	C5			
3300,3301 6100,9100 Advisor,Autocorr, MiniToor, Mini Torr plus (BCI*)		See above *	See the above sensor*	1.1m	See above sensors*	U708-76		
			U410-76	Adult finger /3.0m	C7	No need		
			U410S-76	Adult soft /3.0m	C8			
			U810-76	Multi-site Y/3.0m	C7			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Criticare (CSI)								
510,503,504, 507,508,581, 1100,8100,Poet, Scholar (CSI*)		934-10LN	U410-05	Adult finger /3.0m	C8	No need		
			U110-05	Pediatric finger /3.0m	D1			
		975AD-10	U410S-05	Adult soft /3.0m	D1			
			U110S-05	Pediatric soft /3.0m	D1			
			U210S-05	Infant soft /3.0m	D1			
			U310-05	Neonate wrap /3.0m	D1			
			U610-05	Adult/Neo wrap /3.0m	C8			
			U610B-05	Pedia/infant wrap/3.0m	C8			
			U810-05	Multi-site Y/3.0m	C8			
			U910-05	Adult ear clip /3.0m	D1			
			U403-05	Adult finger /1.1m	C3	U708-05		
			U103-05	Pediatric finger /1.1m	C4			
			U403S-05	Adult soft /1.1m	C4			
			U103S-05	Pediatric soft /1.1m	C4			
			U203S-05	Infant soft /1.1m	C4			
			U303-75	Neonate wrap /1.1m	C4			
			U603-05	Multi-site Y/1.1m	C3			
			U603B-05	Adult/Neo/1.1m	C3			
			U803-05	Pedia/infant wrap/1.1m	C4			
			U903-05	Adult ear clip /1.1m	C4			
503SPOT,504Plus, 504DX,506DX (*CSI)		934-10DN	U410-75	Adult finger /3.0m	C7	No need		
			U110-75	Pediatric finger /3.0m	C8			
			U410S-75	Adult soft /3.0m	C8			
			U110S-75	Pediatric soft /3.0m	C8			
			U210S-75	Infant soft /3.0m	C8			
			U310-75	Neonate wrap /3.0m	C8			
			U610-75	Adult/Neo wrap/3.0m	C7			
			U610B-75	Pedia/infant wrap/3.0m	C7			
			U810-75	Multi-site Y/3.0m	C7			
			U910-75	Adult ear clip /3.0m	C8			
503SPOT,503DX,504Plus, 504DX,506DX,506DXNP (*CSI)			U410-75N	Adult finger/3.0m	C7	No need		
			U410S-75N	Adult soft/3.0m	C8			
			U810-75N	Multi-site Y/3.0m	C7			
Digicare								
Digimax5500 (*BCI)		3044	U403-06	Adult finger /1.1m	C3	U708-76		
			U403S-06	Adult soft /1.1m	C5			
		3043	U803-06	Multi-site Y/1.1m	C3			
			U410-76	Adult finger /3.0m	C7	No need		
			U410S-76	Adult soft /3.0m	C8			
			U810-76	Multi-site Y/3.0m	C7			
			U410-52	Adult finger /3.0m	C7	No need		
			U410S-52	Adult soft /3.0m	C8			
			U810-52	Multi-site Y/3.0m	C7			


OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.	
Datascope									
Accusat, Accutor 3/4 sat Passport, Accutor plus (*Datascope)			U403-07	Adult finger /1.1m	C4	U708-07			
			U103-07	Pediatric finger /1.1m	C6				
			U403S-07	Adult soft /1.1m	C6				
			U103S-07	Pediatric soft /1.1m	C6				
			U203S-07	Infant soft /1.1m	C6				
			U303-07	Neonate wrap /1.1m	C6				
			U603-07	Adult/Neo wrap/1.1m	C5				
			U603B-07	Pedia/infant wrap/1.1m	C5				
			U803-07	Multi-site Y/1.1m	C5				
		U903-07	Adult ear clip /1.1m	C4					
			U410-07	Adult finger /3.0m	C8	No need			
			U410S-07	Adult soft /3.0m	D1				
		U810-07	Multi-site Y/3.0m	C8					
Accutor plus (*M Tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R			
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6				
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4				
		LNOP DC-12	U410-15	Adult finger /3.0m	C8	No need			
		LNOP DC-12	U410S-15	Adult soft /3.0m	D1				
			U810-15	Multi-site Y/3.0m	C8				
	PassportII, Expert Pssport XG, Accutor plus (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-69		
			DS100A	U403S-01	Adult soft /1.1m	C5			
D-YS			U803-01	Multi-site Y/1.1m	C3				
Accutor-Trio (*Non-Oximax)		See the above	See the above	See the above	See the above	U708-01			
Trio Nouvelle Generation (* Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710-70			
		DS100A	U403S-01XP	Adult soft /1.1m	C7				
		D-YS	U803-01XP	Multi-site Y/1.1m	C5				
Dima (Italia)									
(*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01			
		DS100A	U403S-01	Adult soft /1.1m	C5				
		D-YS	U803-01	Multi-site Y/1.1m	C3				
Dixon									
D6 D8			U410-205	Adult finger /3.0m	C8	No need			
			U410S-205	Adult soft/3.0m	D1				
			U810-205	Multi-site Y/3.0m	C8				
			U403-62D	Adult finger /1.1m	C8	U708-208			
			U403S-62D	Adult soft/1.1m	D1				
			U803-62D	Multi-site Y/1.1m	C8				
Dixon									


OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Nonin								
8500 2500 9840 Avant 4000,9600, 9700 Avant 2120 PulseSense Model 7500,7500FO (Nonin*)		8000AA	U403-08	Adult finger /1.1m	C3	No need		
		8000AP	U103-08	Pediatric finger /1.1m	C5			
		8000SL	U403S-08	Adult soft /1.1m	C5			
		8000SM	U103S-08	Pediatric soft /1.1m	C5			
			U203S-08	Infant soft /1.1m	C5			
			U303-08	Neonate wrap /1.1m	C5			
			U603-08	Adult/Neowrap /1.1m	C5			
			U603B-08	Pedia/infant wrap /1.1m	C3			
			U803-08	Multi-site Y /1.1m	C3			
			8000Q2	U903-08	Adult ear clip /1.1m			
Avant 4000 Wearable Digital Pulse Oximetry			U410-08L	Adult finger /3.0m	C5	No need		
			U401-146	Adult finger /0.3m	D1			
			U101-146	Pediatric finger /0.3m	D1			
			U401S-146	Adult soft /0.3m	D1			
8600		The same as above	The same as above	1.1m	The same as above	NA		
8604	 	The same as above	The same as above	1.1m	The same as above	U708-08		
			U410-08	Adult finger /3.0m	C8	No need		
Novametrix								
512/513/710/715, MARS			U403-73	Adult finger /1.1m	C5	No need		
			U103-73	Pediatric finger /1.1m	C7			
			U403S-73	Adult soft /1.1m	C7			
			U103S-73	Pediatric soft /1.1m	C7			
			U203S-73	Infant soft /1.1m	C7			
			U303-73	Neonate wrap /1.1m	C7			
			U603-73	Adult/Neowrap /1.1m	C5			
			U603B-73	Pedia/infant wrap /1.1m	C5			
			U803-73	Multi-site Y /1.1m	C5			
			U903-73	Adult ear clip /1.1m	C7			
			U403-103	Adult finger /1.1m	C8	NA		
			U103-103	Pediatric finger /1.1m	D1			
			U403S-103	Adult soft /1.1m	D1			
			U103S-103	Pediatric soft /1.1m	D1			
			U203S-103	Infant soft /1.1m	D1			
			U303-103	Neonate wrap /1.1m	D1			
			U603-103	Adult/Neowrap /1.1m	C8			
			U603B-103	Pedia/infant wrap /1.1m	C8			
			U803-103	Multi-site Y /1.1m	C8			
			U903-103	Adult ear clip /1.1m	D1			
Novametrix								

Novametrix




OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Novametrix								
505/510/511 /515/B/C/520 7000/7100,Oxypleth (Novametrix*)			U410-03	Adult finger /3.0m	C8	No need		
			U110-03	Pediatric finger /3.0m	D1			
			U410S-03	Adult soft /3.0m	D1			
			U110S-03	Pediatric soft /3.0m	D1			
			U210S-03	Infant soft /3.0m	D1			
			U310-03	Neonate wrap /3.0m	D1			
			U610-03	Adult/Neowrap /3.0m	C8			
			U610B-03	Pedia/infant wrap /3.0m	C8			
			U810-03	Multi-site Y/3.0m	C8			
			U910-03	Adult ear clip /3.0m	D1			
PaceTech/Generra								
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-38		
			U410-38	Adult finger /3.0m	C7	No need		
			U410S-38	Adult soft /3.0m	C8			
			U810-38	Multi-site Y/3.0m	C7			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-39		
			U410-39	Adult finger /3.0m	C8	No need		
			U410S-39	Adult soft /3.0m	D1			
			U810-39	Multi-site Y/3.0m	C8			
Mediaid/Palco								
M300/M305/340, M5305/M5340			U410-18	Adult finger /3.0m	C8	No need		
			U110-18	Pediatric finger /3.0m	D1			
			U410S-18	Adult soft /3.0m	D1			
			U110S-18	Pediatric soft /3.0m	D1			
			U210S-18	Infant soft /3.0m	D1			
			U310-18	Neonate wrap /3.0m	D1			
			U610-18	Adult/Neowrap /3.0m	C8			
			U610B-18	Pedia/infant wrap /3.0m	C8			
			U810-18	Multi-site Y/3.0m	D1			
M30/40,M100,M130, M900,M960P, M6 Aster,M4 Nova, M1000			U410-148	Adult finger /3.0m	D1.3	No need		
			U110-148	Pediatric finger /3.0m	D1.6			
			U410S-148	Adult soft /3.0m	D1.6			
			U110S-148	Pediatric soft /3.0m	D1.6			
			U210S-148	Infant soft /3.0m	D1.6			
			U310-148	Neonate wrap /3.0m	D1.6			
			U610-148	Adult/Neowrap /3.0m	D1.3			
			U610B-148	Pedia/infant wrap /3.0m	D1.3			
			U810-148	Multi-site Y/3.0m	D1.6			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
☆UTAS								
☆UM-300			U410-129	Adult finger /3.0m	C8			
			U410S-129	Adult soft /3.0m	D1			
			U810-129	Adult ear clip /3.0m	C8			
Unicare(China)								
☆F8/S			U410-120	Adult finger /3.0m	C7			
			U410S-120	Adult soft /3.0m	C8			
			U810-120	Adult ear clip /3.0m	C7			
Risingmed								
RPM-9000A			U410-217	Adult finger /3.0m	C7	No need		
			U410S-217	Adult soft/3.0m	C8			
			U810-217	Multi-site Y/3.0m	C7			
		DS100A	U403-01	Adult finger /1.1m	C3	U708-217		
		DS100A	U403S-01	Adult soft/1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
RGB medical								
(*BCI)		3044	U403-06	Adult finger /1.1m	C3	U708-76		
			U403S-06	Adult soft /1.1m	C3			
		3078	U803-06	Adult ear clip /1.1m	C3			



OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Schiller ODM/Bruker								
Propaq CS, Encore Series, Atlas (* non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
TM 910-Defigard 500(M) (* non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m/3.0m	See Nellcor	U708-63		
Argus VCM (*Oximax)		DS100A	U403-01XP	Adult finger /1.1m	C5	U710-70		
		DS100A	U403S-01XP	Adult soft /1.1m	C7			
		D-YS	U803-01XP	Multi-site Y/1.1m	C5			
CM8 SM783,784,785, TM910		See RGB medical	See RGB medical	See RGB medical	See RGB medical	U708-17		
			U410-17	Adult finger /3.0m	C8	No need		
			U410S-17	Adult soft /3.0m	D1			
			U810-17	Multi-site Y/3.0m	C8			
TM7,Argus LCM, (*M tech)		1863/ LNCS DCI	U403-49	Adult finger /1.1m	C4	U708M-17		
		1863/ LNCS DCI	U403S-49	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49	Multi-site Y/1.1m	C4			
			U410-17M	Adult finger /3.0m	C8	No need		
			U410S-17M	Adult soft /3.0m	D1			
			U810-17M	Multi-site Y/3.0m	C8			
Propaq CS, Encore, Micropaq (New) (M) (*M tech)		See M	See M	1.1m	See M	U708M-15R		
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
			U810-15R	Multi-site Y/3.0m	C8			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Spacelabs								
90419/01,90651A/07, 90501, Ultraview 90496 (*Novametrix)			U410-03	Adult finger /3.0m	C8	No need		
			U110-03	Pediatric finger /3.0m	D1			
			U410S-03	Adult soft /3.0m	D1			
			U110S-03	Pediatric soft /3.0m	D1			
			U210S-03	Infant soft /3.0m	D1			
			U310-03	Neonate wrap /3.0m	D1			
			U610-03	Adult/Neo wrap/3.0m	C8			
			U610B-03	Pedia/infant wrap/3.0m	C8			
			U810-03	Multi-site Y/3.0m	C8			
			U910-03	Adult ear clip /3.0m	D1			
90351-0/6,90465, 90466, 90456,PC exoress (*Non-Oximax)		DS100A	U403-01	Adult finger /1.1m	C3	U708-04		
		DS100A	U403S-01	Adult soft /1.1m	C5			
		D-YS	U803-01	Multi-site Y/1.1m	C3			
			U410-04	Adult finger /3.0m	C8	No need		
			U110-04	Pediatric finger /3.0m	D1			
			U410S-04	Adult soft /3.0m	D1			
			U110S-04	Pediatric soft /3.0m	D1			
			U210S-04	Infant soft /3.0m	D1			
			U310-04	Neonate wrap /3.0m	D1			
			U610-04	Adult/Neo wrap/3.0m	C8			
			U610B-04	Pedia/infant wrap/3.0m	C8			
			U810-04	Multi-site Y/3.0m	C8			
			U910-04	Adult ear clip /3.0m	D1			
Ultraview (*Non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-74		
			U410-74	Adult finger /3.0m	D1			
			U410S-74	Adult soft /3.0m	D2	No need		
			U810-74	Multi-site Y/3.0m	D1			
Ultraview (M Tech)			U410-74MR	Adult finger /3.0m	D1	No need		
			U410S-74MR	Adult soft /3.0m	D2			
			U810-74MR	Multi-site Y/3.0m	D1			
Ultraview (Oximax Tech)			U410-74P	Adult finger /3.0m	D1	No need		
			U410S-74P	Adult soft /3.0m	D2			
			U810-74P	Multi-site Y/3.0m	D1			
PRO2 (*BCI)		3044	U403-06	Adult finger /1.1m	C3	No need		
			U403S-06	Adult soft /1.1m	C5			
		3043	U803-06	Multi-site Y/1.1m	C3			

OEM Manufacturer & Models	Sensor Connector	REF# OEM	REF# UNIMED	Sensor Descriptions	Price Code	Ext - Cable	Cable Plug	Sensor Conn.
Takaoka								
Takaoka 9504, 9504ECG, Oxipani , Oxifast			U410-99	Adult finger /3.0m	C7	No need		
			U410S-99	Adult soft /3.0m	C8			
			U810-99	Multi-site Y/3.0m	C7			
Uray								
UR-8000 UR-9000			U410-143	Adult finger /3.0m	C7	No need		
			U410S-143	Adult soft /3.0m	C8			
			U810-143	Multi-site Y/3.0m	C7			
Zoll								
M series AED (*M tech)		1863/ LNCS DCI	U403-49R	Adult finger /1.1m	C4	U708M-15R		
		1863/ LNCS DCI	U403S-49R	Adult soft /1.1m	C6			
		2258/ LNCS YI	U803-49R	Multi-site Y/1.1m	C4			
		LNOP DC-12	U410-15R	Adult finger /3.0m	C8	No need		
		LNOP DC-12	U410S-15R	Adult soft /3.0m	D1			
			U810-15R	Multi-site Y/3.0m	C8			
Zondan (China)								
Apollo N2, Apollo N4 Apollo N5 (M tech optional) DS120C			U410-45	Adult finger /3.0m	C7	No need		
			U410S-45	Adult soft /3.0m	C8			
			U810-45	Multi-site Y/3.0m	C7			
ZD120B ZD120D ZD120E (*Digital Tech)			U410-66D	Adult finger /3.0m	C8	No need		
			U410S-66D	Adult soft /3.0m	D1			
			U810-66D	Multi-site Y/3.0m	C8			
ZD120D/ZD120E			U410-45	Adult finger /3.0m	C7	No need		
			U410S-45	Adult soft /3.0m	C8			
			U810-945	Multi-site Y/3.0m	C7			
Apollo N1/N2/N3/N4/N5			U410-210M	Adult finger /3.0m	D1	No need		
			U410S-210M	Adult soft /3.0m	D2			
			U810-210M	Multi-site Y/3.0m	D1			
Unicare(China)								
(*Digital tech)			U410-66D	Adult finger /3.0m	C8	No need		
			U410S-66D	Adult soft /3.0m	D1			
			U810-66D	Multi-site Y/3.0m	C8			

Please be Reminded : The other types of sensors like U100, U100S, U200S, U300, U600, U600B and U900 of the listed respected brands are available upon demands.