




1 Reusable Temp. Probe

Compatible With Major Patient Monitor Brands

Air Shielded(AS)

Datex-Ohmeda(DX)

Marquette0(MQ)

Philips(HP)

Spacelabs(SL)

Siemens(SM)

YSI(YSI)

Features 产品特点

- Accurate measurement
- Quick response
- Low noise signal transfer
- Fully compatible with OEM technology
- 12 month warranty


Specification

Probe type	General purpose; Skin
Cable length	3.0m/1.5m;unless requirement applied
Cable diameter	2.60mm(Skin purpose) 4.00mm(General purpose)
Resister type	NTC/YSI400/2.252K@25°C NTC/10K@25°C
Cable Material	TPU jacket ,Light gray
Connector	Golden plated pins


1.1 Reusable Temp. Probe

Reference Number Rule Description

TSM2-15AS	
T	Temperature probe
SM	Factory/model (See others left)
2	Dual channel(Default is single channel)
15	1.5m cable length(Default:=3.0m)
A	Adult(P= Pediatric)
S	Skin(G=general purpose)



OEM & Models	Unimed REF#	OEM REF#	Application	Cable length	Monitor end plug	Price code
Air Shield						
	TAS-15AS		Adult skin	1.5m		C3
	TAS-15PS		Pediatric skin	1.5m		C4
	TAS-15AG		Adult general	1.5m		C3
	TAS-15PG		Pediatric general	1.5m		C4
Artema-S&W						
	TSW-AS		Adult skin	3.0m		C3
	TSW-PS		Pediatric skin	3.0m		C4
	TSW-AG		Adult general	3.0m		C3
	TSW-PG		Pediatric general	3.0m		C4
ATOM						
Infant Warmer: Model V-3200D/V-505 Infant Incubator: Model V-2100G/V-2200 Transport Incubator:V-707/V-808	TAT-15PS		Pediatric skin	1.5m		C4
Biolight						
V6 Vital Sign AnyView A2/A3/A5/A6/A8	TBLT-AS		Adult skin	3.0m		C3
	TBLT-PS		Pediatric skin	3.0m		C4
	TBLT-AG		Adult general	3.0m		C3
	TBLT-PG		Pediatric general	3.0m		C4
Beijing Julongsanyou						
VInfant Radiant Warmer	T20-15PS		Pediatric skin	1.5m		B7
BNT-1000/BNT-2000/BNT-2000B	T20-PS		Pediatric skin	3.0m		B7
Beijing Julongsanyou						

OEM & Models	Unimed REF#	OEM REF#	Application	Cable length	Monitor end plug	Price code
Beijing Julongsanyou						
Infant Incubator	TBNT-15PS		Pediatric skin	1.5m		D1
B-800/B-1000/B-2000/B-3000	TBNT-PS		Pediatric skin	3.0m		D1
David						
HKN-90/HKN-93HKN-93A	TDV-15PS		Pediatric skin	1.5m		D1
HKN-93B/HKN-9010/HKN-2000	TDV-PS		Pediatric skin	3.0m		D1
GE Marquette						
Eagle, Tram, Solar, Mac-lab Cath lab	TMQ-AS		Adult skin	3.0m		C3
	TMQ-PS		Pediatric skin	3.0m		C4
	TMQ-AG		Adult general	3.0m		C3
	TMQ-PG		Pediatric general	3.0m		C4
	TMQ2-AS	Dual temp	Adult skin(2)	3.0m		D1
	TMQ2-AS/AG	Dual temp	Adult skin & general	3.0m		D1
Philips						
All	THP-AS	21078A	Adult skin	3.0m		C3
	THP-PS		Pediatric skin	3.0m		C4
	THP-AG	21075A,21076A	Adult general	3.0m		C3
	THP-PG		Pediatric general	3.0m		C4
Mindray						
T5/T8 DPM 6,DPM 7	TMR-AS	0011-30-37393	Adult skin	3.0m		C3
	TMR-PS	0011-30-37395	Pediatric skin	3.0m		C4
	TMR-AG	0011-30-37392	Adult general	3.0m		C3
	TMR-PG	0011-30-37394	Pediatric general	3.0m		C4
Datascope Passport V, V12, V21	TDT-AS	040-000057-00	Adult skin	3.0m		C4
	TDT-PS	040-000058-00	Pediatric skin	3.0m		C5
	TDT-AG	040-000055-00	Adult general	3.0m		C4
	TDT-PG	040-000056-00	Pediatric general	3.0m		C5
Siemens/Drager						
	TSM-15AS	4329822(1.5m)	Adult skin	1.5m		C4
	TSM-AS	5204669(3.0M)		3.0m		
	TSM-15PS		Pediatric skin	1.5m		C5
	TSM-PS			3.0m		
	TSM-15AG	4329889(1.5m)	Adult general	1.5m		C4
	TSM-AG	5204644(3.0m)		3.0m		
	TSM-15PG	4329848(1.5m)	Pediatric general	1.5m		C5
TSM-PG	5204651(3.0m)		3.0m			
Spacelabs						
	TSL-AS		Adult skin	3.0m		C7
	TSL-PS		Pediatric skin	3.0m		C8
	TSL-AG		Adult general	3.0m		C7
	TSL-PG		Pediatric general	3.0m		C8
	TSL2-AS		Adult skin(2)	3.0m		D3
	TSL2-AS/AG		Adult skin & general	3.0m		D3
YSI						
YSI 400 Series	T2252-AS	409B	Adult skin	3.0m		B6
	T2252-PS	429	Pediatric skin	3.0m		B7
	T2252-AG	401	Adult general	3.0m		B6
	T2252-PG	402	Pediatric general	3.0m		B7
YSI 700 Series	T700-AS	709B	Adult skin	3.0m		C3
	T700-PS	729	Pediatric skin	3.0m		C4
	T700-AG	701	Adult general	3.0m		C3
	T700-PG	702	Pediatric general	3.0m		C4

OEM & Models	Unimed REF#	OEM REF#	Application	Cable length	Monitor end plug	Price code
10K@25°C						
	T10-AS/T10-09AS		Adult skin	3.0m/0.9m		B6/B5
	T10-PS		Pediatric skin	3.0m		B7
	T10-AG		Adult general			B6
	T10-PG		Pediatric general	3.0m		B7
Ohmeda						
Ohmeda Incubator: Giraffe/Omnibed	TOH-PS		Neonate skin	3.0m		C4
Air Temperature						
	T100-GP		For air Temp. measurement	3.0m		C4

Unimed REF#	Descriptions & Application	PCS /Per box	Price code
T004-16	Disposable temperature probe cover Used with reusable skin temp. probe to fixed the probe onto the skin and reflect the ambient light to ensure the accurate measurement.	16	A5



The probe cover package not includes the reusable probe. Separated order required.

