







③ SpO2 Interface Cable




Specification




Cable material	TPU Jacket, Light gray, Blue
Connector	Nylon plug &Golden plated pins
Cable length	1.1m/2.2m/3.0m
Cable diameter	4.0mm/5.0mm
Package	1 pcs per bag








REF#	OEM REF#	Sensor Tech.	Cable length	Monitor end Conn.	Sensor end Conn.	Price code
Artema						
U708-80	59831-0	Nellcor,non-Oximax	2.2m			C3
BCI						
U708-76		BCI	2.2m			B8
U708-06		BCI	2.2m			B8
Biolight						
U710X-51		Nellcor Oximax	3.0m			C3
U708-68		Digital sensor	2.2m			C3
Bionet						
U708-58		Nellcor,non-Oximax	2.2m			B8
U708-108		Bionet	2.2m			C3
Criticare						
U708-05	518LD	CSI	2.2m			C3
U708-75	518DD	CSI	2.2m			B8
Datascope						
U708-07	0012-00-0516-02	DATASCOPE	2.2m			C3
Dolphin						
U708-50		Dolphin digital sensor	2.2m			C3
Drager						
U708-28	5720084	Nellcor,non-Oximax	2.2m			D4
EDAN						
U708X-191	MS1-30131 (M9,M9A,M8,M8A,M8B)	Nellcor-Oximax	2.2m			C3
U708-29	M13-30691 (M9,M9A,M8,M8A,M8B, M3,M3B)	EDAN Oximetry	2.2m			B8








REF#	OEM REF#	Sensor Tech.	Cable length	Monitor end Conn.	Sensor end Conn.	Price code
GE healthcare						
U708-21	2006644-001 E9004GE	Nellcor,non-Oximax	2.2M			C3
U710X-21	2021406-001	Nellcor,Oximax	3.0M			C7
U708-21M	M technology	Nellcor Non-Oximax	2.2M			C3
U708M-21R	LNC-10-GE 2027263-002 (REF 2016)	M tech.(Masimo LNCS)	2.2M			C3
U708M-21P	20169041-001	M tech.(Masimo LNOP)	2.2M			C7
U708-09	OXY-C3	Datex	2.2m			C3
U708-81	Datex>>Oxytip Used with Oxy-C3 cable	Oxytip	2.2m			C3
U708-55	30344358	Hellige(Nellcor)	2.2m			C3
U708-25	4033CAX	Corometrics(Nellcor)	2.2m			C3
U708-72*	All Ohmeda Oximeters	Ohmeda Used with Unimed N500 or U500 Disposable sensors	2.2m			C3
U708-160	OXY-MC3	GE OxyTip+ sensor	2.2m			C7
U708-102 U704-102	OXY-OL3	GE OxyTip+ sensor	2.2m 1.1m			C5
U708-209	OXY-SL3	GE OxyTip+ sensor	2.2m			C5
U708-121 U710-121	OXY-ES3	GE OxyTip+ sensor	2.2m 3.0m			C7
U708-81 U701-81	OXY-SLA OXY-SLC	GE OxyTip+ sensor	2.2m 0.3m			C5 C3
U708-111	TS-G3	GE OxyTip+ sensor	2.2m			C5
U708-112	TS-H3	GE OxyTip+ sensor	2.2m			C5
U708-116	TS-M3	GE OxyTip+ sensor	2.2m			C7
U708-119	TS-N3	GE OxyTip+ sensor	2.2m			C5
Goldway						
U708-11		Nellcor Non-Oximax	2.2m			C3
U708-37		Nellcor Non-Oximax	2.2m			B8

REF#	OEM REF#	Sensor Tech.	Cable length	Monitor end Conn.	Sensor end Conn.	Price code
Goldway						
U708-38		Nellcor Non-Oximax	2.2m			B8
U710X-38		Nellcor Oximax	3.0m			C3
U708-11		Nellcor Non-Oximax	2.2m			C3
U708-37		Nellcor Non-Oximax	2.2m			B8
Invivo						
U708-14		Nellcor,non-Oximax	2.2m			C3
U708-14M		(Masimo) Nellcor, non-Oximax	2.2m			C3
U708M-14		Masimo LNCS sensor	2.2m			C3
Kontron						
U708-13		Nellcor,non-Oximax	2.2m			C3
U708-53		Nellcor,non-Oximax	2.2m			B8
M Technology						
U708-15	M technology	(Masimo) Nellcor, non-Oximax	2.2m			C3
U708M-15R	LNCS Patient Cable LNC-4(REF 2017) LNC-10(REF 1814) LNC-14(REF 2013)	Masimo LNCS sensor	2.2m			C3
U708M-15PR	PC 04 / PC 08 / PC 12	Masimo LNOP sensor	2.2m			C3
U704M-49	LNC-4-Ext (REF 2021) For LNCS Cable	Masimo LNCS sensor	1.1m			B8
U701M-125	LNCS to LNOP LNCS to PC series (REF 1816)	Masimo LNCS sensor	0.3m			C3
MEK						
U708-56		MEK (Nellcor non-oximax)	2.2m			C3
Mennen						
U708-54		Nellcor,non-Oximax	2.2m			C3
U708-24		Nellcor,non-Oximax	2.2m			D4
U708-59M	M technology	M tech. (Nellcor, non-Oximax)	2.2m			C3
U708M-59	M technology	Masimo LNCS sensor	2.2m			C3
U710X-59		Nellcor Oximax	3.0m			C3

REF#	OEM REF#	Sensor Tech.	Cable length	Monitor end Conn.	Sensor end Conn.	Price code
Mindray						
U708-29	0010-20-42594 Accutorr V, DPM 3, DPM 4, DPM 5, Passport V	Nellcor,non-Oximax	2.2M			B8
U708-29M	M technology For old PM 6000 only	Nellcor,non-Oximax	2.2M			B8
U708M-29	M technology For old PM 6000 only	Masimo LNCS	2.2M			B8
U708-30M	M technology (0010-30-42625) DPM 3	Nellcor,non-Oximax	2.2M			B8
U708M-30	DPM 4 DPM 5	Masimo LNCS	2.2M			B8
U708-22	0515-30-11221	Nellcor,non-Oximax	2.2M			C3
U708-22M	M technology 0010-30-12452	non-Oximax Nellcor	2.2M			C3
U708M-22	M technology 0010-30-12452	Masimo LNCS	2.2M			C3
U708-31	CSI Technology	CSI	2.2M			C3
U708M-48R	0010-30-42738 T5/T8 DPM 6, DPM 7	Masimo LNCS	2.2m			C3
U708-48	0010-20-42710 (C) T5/T8 DPM 6, DPM 7 Nellcor Non-Oximax	Nellcor,non-Oximax	2.2m			C3
U710X-48	0010-20-42712 T5/T8, DPM 6, DPM 7 Nellcor Oximax	Nellcor Oximax	3.0m			C4
U708M-15R	(0012-00-1652),1.2m (0012-00-1599),3.1m (0012-00-1653),4.3m Accutorr V, Duo, Trio, Passport 2, PassportV, Spectrum/OR, V12, V21	Masimo LNCS				C3
U708-69	0012-00-1254 Trio, Passport 2, Spectrum	Nellcor,non-Oximax	2.2m			C3
U710-70	0012-00-1464 Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum/OR, V12, V21	Nellcor Oximax	3.0m			C6
U710X-180	0010-20-42595 DPM4	Nellcor Oximax	3.0m			C3
Nellcor/Coviden						
U708-69	SCP-10/MC-10	Nellcor, Non-Oximax	2.2m			C3
U710-85	M-200-13		3.0m			D7
U704-71	DEC-4	Nellcor-Oximax	1.1m			B8
U708-71	DEC-8		2.2m			B8
U710-70	DOC-10	Nellcor-Oximax	3.0m			C6
U710-70X	DOC-10 + Oximax adapter	Nellcor- Non Oximax	3.0m			C8

REF#	OEM REF#	Sensor Tech.	Cable length	Monitor end Conn.	Sensor end Conn.	Price code
Nellcor/Coviden						
U701-01X U708-01X	Oximax adapter	Nellcor, Non-Oximax	0.3m 2.2m			C3
Nihon Kohden						
U710-36	JL-900P BSM/PVM	NIHON KOHDEN	3.0m			C7
U710X-36	JL-650P	Oximax	3.0m			C7
U710M-36R	JL-631P	M Tech	3.0m			C7
U710-106	JL-302T OPV-1500,OLV	NIHON KOHDEN	3.0m			C7
Nonin						
U708-08		Nonin	2.2m			C3
Novametrix						
U708-103	8853	Novametrix	2.2m			C3
Pace tech/Generra						
U708-37		Nellcor non-Oximax	2.2m			B8
U708-38		Nellcor non-Oximax	2.2m			B8
U708-39		Nellcor non-Oximax	2.2m			C3
Philips/HP						
U708-40	M1940A/B	Philips	2.2m			C3
U708-41	M1941A/B		2.2m			C3
U708-20	M1900A/B	Nellcor non-Oximax	2.2m			C3
U708M-20P	M1020-61102 CMS M1020A module.	Masimo LNOP	2.2m			C3
U708-43 U710-43	M1943A M1943AL	Nellcor non-Oximax	2.2m 3.0m			C3
U710X-43	M1943NL	Nellcor,Oximax	3.0m			C7
U708-43MR	M Technology	Nellcor non-Oximax	2.2m			C3
U708M-43R	989803148221 LNC MP10 (REF 2281)	Masimo LNCS sensor	2.2m			C3
U708M-43PR	M1020-61100	Masimo LNOP sensor	2.2m			C7

4 Accessories

REF#	OEM REF#	Sensor Tech.	Cable length	Monitor end Conn.	Sensor end Conn.	Price code
Schiller						
U708-17		BCI	2.2m			C3
U708-17M	M tech.	Nellcor non-Oximax	2.2m			C3
U708M-17	M tech.	Masimo LNCS sensor	2.2m			C3
U708-63 ODAM/Bruker		BCI	2.2m			C3
Siemens/Dräger						
U704-23	3368433	Nellcor non-Oximax	1.1m			C6
U708-23	3375834	Nellcor non-Oximax	2.2m			C6
U710X-23		Nellcor Oximax	3.0m			C7
U704M-23P	MS18680	Masimo [®] LNOP	1.1m			C6
U708M-23PR	MS17041		2.2m			C6
U704M-23R	MS24303	Masimo [®] LNCS	1.1m			C6
U708M-23R	MS17522		2.2m			C6
U708-33	4533840	Nellcor non-Oximax	2.2m			C3
Spacelabs						
U708-04	700-0002-00 90351-D/G, 90465, 90466, 90467, 90489,M77	Nellcor non-Oximax	2.2m			C3
U708-74	700-0030-00 90496,91496-U	Nellcor non-Oximax	2.2m			C6
U710X-74	700-0792-00 91496-N,91220	Nellcor-Oximax	3.0m			C7
*U708M-74	700-0906-00 91496-M	Masimo LNCS	2.2m			C6
*U708M-74P	700-0789-00 91496-M	Masimo LNOP	2.2m			C6
U701-04NV	175-0646-00 90351-D/G,90465, 90466,90467, 90489, 90651-08)	Novametrix- Hypertronic Con	0.3m			C3
U701-74NV	700-0029-00 (90496,91496-U)	Novametrix- Hypertronic Con	0.3m			C6
Welchallyn						
U704X-104	001966-E (PIC 40TM Defibrillator PIC 50TM Defibrillator)	Nellcor Oximax	2.2m			C7
U704-104	001955-U (PIC40 Defibrillator PIC50 Defibrillator)	Nellcor	1.1m			C7

*Refurbish only

					
REF#	Product Name	Purpose of Use	Package	Price	Figures
Y04-E041	Ear sensor hanger	Used with Unimed ear clip sensor U900 series or Other brands ear sensor(Gray)	1 pcs	A5	A
UEC-01	Adult ear clip Or as animal small size lingua clip	Used with D-YS or Unimed U800 sensor(Black)	1 pcs	A5	B
UVC-01	Animal lingua clip,Big size	Used with D-YS or Unimed U800 sensor(Gray)	1 pcs	A8	C
UVC-02	Animal lingua clip,Middle size	Used with D-YS or Unimed U800 sensor(Gray)	1 pcs	A7	D
Y04-E019	Reusable Silicon Y sensor wrap tape	Used with Nellcor D-YS or Unimed U800 series sensors	1 pcs	A6	E
E023-12	Reusable 3 holes foam tape For All patient sizes	Used with Nellcor D-YS or Unimed U800 /U600 series sensors	12 pcs Per bag	B5/bag	F
E024-25	Disposable 3M adhesive textile tape For Adult and neonate	Used with Nellcor D-YS or Unimed U800 series sensors	25 pcs Per bag	B7/bag	G
E025-25	Disposable 3M adhesive textile tape for Pediatric and infant	Used with Nellcor D-YS or Unimed U800 series sensors	25 pcs Per bag	B7/bag	H
T002-12	Velcro foam wrap tape	To fix the sensor cable onto the patient wrist	12 pcs per bag	C1.4/bag	I