



1 Reusable SpO2 Sensor


Specification

|                         |                                                                    |
|-------------------------|--------------------------------------------------------------------|
| Sensor Connector        | OEM compatible, Nylon adapter with gold plated pins                |
| Cable material & color  | TPU jacket; light gray or blue                                     |
| Cable length & diameter | 1.1m/1.5m/3.0m; ϕ3.20mm/4.00mm or 5.00x2.50mmtwo line cable        |
| Sensor house            | ABS shell(finger clip); Silicon rubber boots( soft and multi-site) |

Compatible With Major Patient Monitor Brands

Nellcor  
Ohmeda  
Philips  
Marquette  
Spacelabs  
Novamatrix

Criticare  
BCI  
Datascope  
Datex  
Nihon Kohden  
Palco  
.....



| Series系列     | U400          | U100             | U400S           | U100S              | U200S                          | U300            |
|--------------|---------------|------------------|-----------------|--------------------|--------------------------------|-----------------|
| REF to Fig   | A             | B                | C               | D                  | E                              | F               |
| Name         | Adult Finger  | Pediatric Finger | Adult soft tip  | Pediatric Soft tip | Infant Soft tip                | Neonate Wrap    |
| Patient size | +40kg         | 15-40kg          | +40kg           | 15-40kg            | 3-15kg                         | 1-3kg           |
| Cable color  | Gray Cable    | Gray Cable       | Blue Cable      | Gray two line      | Gray two line                  | Gray two line   |
| Series系列     | U800          | U900             | U600            | U600B              | U800E                          | U800V           |
| REF to Fig   | G             | H                | I               | J                  | K                              | L               |
| Name         | Multi-site    | Adult Ear clip   | Adult/Neo. Wrap | Pedia/Inf. Wrap    | Multi-site                     | Forehead/Rectal |
| Patient size | All patients  | +40kg            | +40kg or -3kg   | 15-40kg            | All patients                   | Adult           |
| Cable color  | Gray two line |                  | Gray two line   | Gray two line      | Sensormate Used with DPT tapes | Gray Cable      |


How to Identify the REF# of Unimed SpO<sub>2</sub> Sensors

Reference Number Rule Description

| U403S-01XE |                                                                                                                   |
|------------|-------------------------------------------------------------------------------------------------------------------|
| U          | Unimed brand                                                                                                      |
| 4          | adult sensor (1: pediatric; 2: infant;3: neonate; 5: disposable; 6 wrap; 7: ext-cable; 8: multi-site;9: ear clip) |
| 03         | 1.1m cable length (03: 1.1m;05:1.5m;08: 2.2m;10: 3.0m)                                                            |
| S          | Soft sensor (B:U style shell; Default: shell sensor or neonate wrap, alternative wrap, multi-site)                |
| 01         | Nellcor code (02:Ohmeda;03:Novametrix; 04: Spacelabs; 05:Criticare; 06: BCI...)                                   |
| X          | Oximax technology( M: M technology, D: digital technology; Default: standard                                      |
| E          | Sensormate sensor (Default: non-Sensormate sensor)                                                                |



Cross Reference


Ending "X" can be replaced by "XP" or "P" for fully compatible with all oximax oximetry

| OEM Manufacturer & Models                                                     | Sensor Connector                                                                    | REF# OEM | REF# UNIMED | Sensor Descriptions | Price Code | Ext - Cable | Cable Plug | Sensor Conn. |
|-------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------|-------------|---------------------|------------|-------------|------------|--------------|
| 3F Medical (China)                                                            |                                                                                     |          |             |                     |            |             |            |              |
| IRIS 300/IRIS 400<br>Phoebe PMS8300/8310<br>PhoebeProPMS8300/<br>8310PMS8110A |  |          | U410-189    | Adult finger /3.0m  | C7         | No need     |            |              |
|                                                                               |                                                                                     |          | U410S-189   | Adult soft /3.0m    | C8         |             |            |              |
|                                                                               |                                                                                     |          | U810-189    | Multi-site Y/3.0m   | C7         |             |            |              |
| Adecon (China)                                                                |                                                                                     |          |             |                     |            |             |            |              |
| MCF-21L,<br>DK-8000D,DK-8000L                                                 |  |          | U410-120    | Adult finger /3.0m  | C7         | No need     |            |              |
|                                                                               |                                                                                     |          | U410S-120   | Adult soft /3.0m    | C8         |             |            |              |
|                                                                               |                                                                                     |          | U810-120    | Multi-site Y/3.0m   | C7         |             |            |              |
| Anke (China)                                                                  |                                                                                     |          |             |                     |            |             |            |              |
| ASC-553                                                                       |  |          | U410-37     | Adult finger /3.0m  | C7         | No need     |            |              |
|                                                                               |                                                                                     |          | U410S-37    | Adult soft成/3.0m    | C8         |             |            |              |
|                                                                               |                                                                                     |          | U810-37     | Multi-site Y/3.0m   | C7         |             |            |              |
| ASC-553<br>(*Oximax tech)                                                     |  |          | U410-78P    | Adult finger /3.0m  | C8         | No need     |            |              |
|                                                                               |                                                                                     |          | U410S-78P   | Adult soft /3.0m    | D1         |             |            |              |
|                                                                               |                                                                                     |          | U810-78P    | Multi-site Y/3.0m   | C8         |             |            |              |
| ASC-553A                                                                      |  |          | U410-82     | Adult finger /3.0m  | C7         | No need     |            |              |
|                                                                               |                                                                                     |          | U410S-82    | Adult soft /3.0m    | C8         |             |            |              |
|                                                                               |                                                                                     |          | U810-82     | Multi-site Y/3.0m   | C7         |             |            |              |
| ASC-553A<br>(*Ohmeda Tech)                                                    |  |          | U410-107    | Adult finger /3.0m  | C7         | No need     |            |              |
|                                                                               |                                                                                     |          | U410S-107   | Adult soft /3.0m    | C8         |             |            |              |
|                                                                               |                                                                                     |          | U810-107    | Multi-site Y/3.0m   | C7         |             |            |              |

| OEM Manufacturer & Models                                                            | Sensor Connector                                                                      | REF# OEM    | REF# UNIMED              | Sensor Descriptions    | Price Code            | Ext - Cable        | Cable Plug                                                                            | Sensor Conn.                                                                          |
|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-------------|--------------------------|------------------------|-----------------------|--------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Artema/S&W                                                                           |                                                                                       |             |                          |                        |                       |                    |                                                                                       |                                                                                       |
| Diascope,8000, Athena<br>(*Non-Oximax)                                               |    | DS100A      | <b>U403-01</b>           | Adult finger/1.1m      | C3                    | U708-80            |    |    |
|                                                                                      |                                                                                       | DS100A      | <b>U403S-01</b>          | Adult soft /1.1m       | C5                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       | D-YS        | <b>U803-01</b>           | Multi-site Y/1.1m      | C3                    |                    |                                                                                       |                                                                                       |
|                                                                                      |    |             | <b>U410-80</b>           | Adult finger /3.0m     | C8                    | No need            |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U410S-80</b>          | Adult soft /3.0m       | D1                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U810-80</b>           | Multi-site Y/3.0m      | C8                    |                    |                                                                                       |                                                                                       |
| S&W                                                                                  |    |             | <b>U410-19</b>           | Adult finger /3.0m     | C8                    | No need            |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U410S-19</b>          | Adult soft成 /3.0m      | D1                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U810-19</b>           | Multi-site Y/3.0m      | C8                    |                    |                                                                                       |                                                                                       |
| Baxter/Simed                                                                         |                                                                                       |             |                          |                        |                       |                    |                                                                                       |                                                                                       |
| ASAT,<br>S-50/100/100a                                                               |   |             | <b>U410-12</b>           | Adult finger /3.0m     | C8                    | No need            |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U410S-12</b>          | Adult soft /3.0m       | D1                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U810-12</b>           | Multi-site Y/3.0m      | C8                    |                    |                                                                                       |                                                                                       |
| Monarch/Explorer<br>(*Oximax)                                                        |  | DS100A      | <b>U403-01</b>           | Adult finger/1.1m      | C3                    | U708-01            |  |  |
|                                                                                      |                                                                                       | DS100A      | <b>U403S-01</b>          | Adult soft/1.1m        | C5                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       | D-YS        | <b>U803-01</b>           | Multi-site Y/1.1m      | C3                    |                    |                                                                                       |                                                                                       |
| BCI                                                                                  |                                                                                       |             |                          |                        |                       |                    |                                                                                       |                                                                                       |
| 3300,3301,3302,<br>3303,3304,3401<br>(*BCI)                                          |  | 3044        | <b>U403-06</b>           | Adult finge /1.1m      | C3                    | Depends<br>U708-06 |  |  |
|                                                                                      |                                                                                       |             | <b>U103-06</b>           | Pediatric finger /1.1m | C5                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U403S-06</b>          | Adult soft /1.1m       | C3                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U103S-06</b>          | Pediatric soft /1.1m   | C5                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U203S-06</b>          | Infant soft /1.1m      | C5                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U303-06</b>           | Neonate wrap /1.1m     | C5                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       | 3043        | <b>U803-06</b>           | Multi-site Y/1.1m      | C3                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U603-06</b>           | Adult/Neo Wrap/1.1m    | C3                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U603B-06</b>          | Pedia/infant wrap/1.1m | C3                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       | 3078        | <b>U903-06</b>           | Adult ear clip /1.1m   | C5                    |                    |                                                                                       |                                                                                       |
| 3300,3301<br>6100,9100<br>Advisor,Autocorr,<br>MiniToor,<br>Mini Torr plus<br>(BCI*) |  | See above * | See the<br>above sensor* | 1.1m                   | See above<br>sensors* | U708-76            |  |  |
|                                                                                      |                                                                                       |             | <b>U410-76</b>           | Adult finger /3.0m     | C7                    | No need            |                                                                                       |                                                                                       |
|                                                                                      |  |             | <b>U410S-76</b>          | Adult soft /3.0m       | C8                    |                    |                                                                                       |                                                                                       |
|                                                                                      |                                                                                       |             | <b>U810-76</b>           | Multi-site Y/3.0m      | C7                    |                    |                                                                                       |                                                                                       |

| OEM Manufacturer & Models                                            | Sensor Connector                                                                    | REF# OEM | REF# UNIMED | Sensor Descriptions | Price Code | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|----------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------|-------------|---------------------|------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Biosys                                                               |                                                                                     |          |             |                     |            |             |                                                                                       |                                                                                       |
| BPM200<br>BPM700<br>Guardian                                         |    | DS100A   | U403-01     | Adult finger /1.1m  | C3         | U708-01     |    |    |
|                                                                      |                                                                                     | DS100A   | U403S-01    | Adult soft /1.1m    | C5         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     | D-YS     | U803-01     | Multi-site Y/1.1m   | C3         |             |                                                                                       |                                                                                       |
|                                                                      |    |          | U410-27     | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-27    | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-27     | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| Biolight(China)                                                      |                                                                                     |          |             |                     |            |             |                                                                                       |                                                                                       |
| M6,M12,M69                                                           |    |          | U410-47     | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-47    | Adult soft/3.0m     | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-47     | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| M66,M69                                                              |    |          | U410-57     | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-57    | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-57     | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| (*Oximax)                                                            |    |          | U410-51P    | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-51P   | Adult soft/3.0m     | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-51P    | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| M700<br>(*Digital Tech)                                              |    |          | U403-62D    | Adult finger /1.1m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U403S-62D   | Adult soft /1.1m    | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U803-62D    | Multi-site Y/1.1m   | C8         |             |                                                                                       |                                                                                       |
| M8000,M9000,M8000A<br>M9000A,BTD-352A,<br>new M69<br>(*Digital Tech) |  |          | U410-62D    | Adult finger /3.0m  | D1         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-62D   | Adult soft /3.0m    | D2         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-62D    | Multi-site Y/3.0m   | D1         |             |                                                                                       |                                                                                       |
| M900,M7000,M9500<br>A8 Modular<br>(*Digital Tech)                    |  |          | U410-68D    | Adult finger /3.0m  | D1         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-68D   | Adult soft /3.0m    | D2         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-68D    | Multi-site Y/3.0m   | D1         |             |                                                                                       |                                                                                       |
| Anyview A2E<br>Anyview A3<br>V6vet                                   |  |          | U410-205D   | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-205D  | Adult soft/3.0m     | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-205D   | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| M800 handled vet                                                     |  |          | U403-62D    | Adult finger /1.1m  | C8         | U708-205    |  |  |
|                                                                      |                                                                                     |          | U403S-62D   | Adult soft/1.1m     | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U803-62D    | Multi-site Y/1.1m   | C8         |             |                                                                                       |                                                                                       |
| Anyview A3<br>(*Masimo tech)                                         |  |          | U410-87     | Adult finger /3.0m  | D2         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-87    | Adult soft/3.0m     | D3         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-87     | Multi-site Y/3.0m   | D2         |             |                                                                                       |                                                                                       |
| Bionet                                                               |                                                                                     |          |             |                     |            |             |                                                                                       |                                                                                       |
|                                                                      |  | DS100A   | U403-01     | Adult finger /1.1m  | C3         | U708-58     |  |  |
|                                                                      |                                                                                     | DS100A   | U403S-01    | Adult soft /1.1m    | C5         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     | D-YS     | U803-01     | Multi-site Y/1.1m   | C3         |             |                                                                                       |                                                                                       |
| BM3/BM5                                                              |  |          | U410-58     | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-58    | Adult soft /3.0m    | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-58     | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| BM3/BM5(New)                                                         |  |          | U410-108    | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U410S-108   | Adult soft /3.0m    | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |          | U810-108    | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |

| OEM Manufacturer & Models                               | Sensor Connector                                                                      | REF# OEM          | REF# UNIMED | Sensor Descriptions | Price Code | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|---------------------------------------------------------|---------------------------------------------------------------------------------------|-------------------|-------------|---------------------|------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Biocare (China)                                         |                                                                                       |                   |             |                     |            |             |                                                                                       |                                                                                       |
| BP-12B,BP-12C<br>(*Digital Tech)                        |    |                   | U410-61D    | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U410S-61D   | Adult soft/3.0m     | D1         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-61D    | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| BM-9000,BM-9000S<br>FM-801,FM-801V<br>(*Digital Tech)   |    |                   | U410-62D    | Adult finger /3.0m  | D1         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U410S-62D   | Adult soft/3.0m     | D2         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-62D    | Multi-site Y/3.0m   | D1         |             |                                                                                       |                                                                                       |
| BM-9000                                                 |    |                   | U410-29     | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U410S-29    | Adult soft/3.0m     | C8         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-29     | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| IM15                                                    |    |                   | U410-223    | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U410S-223   | Adult soft/3.0m     | D1         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-223    | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| PM-900                                                  |    |                   | U410-143    | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U410S-143   | Adult soft/3.0m     | C8         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-143    | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| CAS Medical                                             |                                                                                       |                   |             |                     |            |             |                                                                                       |                                                                                       |
| Oscillomate 9002<br>9303,9302S<br>(*Nonin)              |    | 8000AA            | U403-08     | Adult finger /1.1m  | C3         | U708-01     |    |    |
|                                                         |                                                                                       |                   | U403S-08    | Adult soft/1.1m     | C5         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U803-08     | Multi-site Y/1.1m   | C3         |             |                                                                                       |                                                                                       |
| 740-750(M)<br>(*M tech)                                 |  | 1863/<br>LNCS DCI | U403-49R    | Adult finger /1.1m  | C4         | U708M-15R   |  |  |
|                                                         |                                                                                       | 1863/<br>LNCS DCI | U403S-49R   | Adult soft/1.1m     | C6         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       | 2258/<br>LNCS YI  | U803-49R    | Multi-site Y/1.1m   | C4         |             |                                                                                       |                                                                                       |
|                                                         |  | LNOP<br>DC-12     | U410-15R    | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       | LNOP<br>DC-12     | U410S-15R   | Adult soft/3.0m     | D1         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-15R    | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| Charmcare                                               |                                                                                       |                   |             |                     |            |             |                                                                                       |                                                                                       |
| CX100<br>(*Non-Oximax)                                  |  | DS100A            | U403-01     | Adult finger /1.1m  | C3         | U708-01     |  |  |
|                                                         |                                                                                       | DS100A            | U403S-01    | Adult soft/1.1m     | C5         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       | D-YS              | U803-01     | Multi-site Y/1.1m   | C3         |             |                                                                                       |                                                                                       |
| ChoiceMed超思(China)                                      |                                                                                       |                   |             |                     |            |             |                                                                                       |                                                                                       |
| MMED3001                                                |  | DS100A            | U403-01     | Adult finger /1.1m  | C3         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       | DS100A            | U403S-01    | Adult soft/1.1m     | C5         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       | D-YS              | U803-01     | Multi-site Y/1.1m   | C3         |             |                                                                                       |                                                                                       |
| MMED100H,<br>MMED300A,<br>MMED6000DP<br>(*Digital tech) |  |                   | U410-84D    | Adult finger /3.0m  | D1         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U410S-84D   | Adult soft/3.0m     | D2         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-84D    | Multi-site Y/3.0m   | D1         |             |                                                                                       |                                                                                       |
| MMED6000DP<br>(*BCI Tech)                               |  |                   | U410-96     | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U410S-96    | Adult soft/3.0m     | C8         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-96     | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| MMED6000DP                                              |  |                   | U410-37     | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U410S-37    | Adult soft/3.0m     | C8         |             |                                                                                       |                                                                                       |
|                                                         |                                                                                       |                   | U810-37     | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| ChoiceMed(China)                                        |                                                                                       |                   |             |                     |            |             |                                                                                       |                                                                                       |



| OEM Manufacturer & Models                                                                                         | Sensor Connector                                                                    | REF# OEM          | REF# UNIMED   | Sensor Descriptions | Price Code    | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|-------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------|---------------|---------------------|---------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| ChoiceMed(China)                                                                                                  |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| MMED3001,MD300K<br>MA300K2,MD300M,<br>MMED6000DP-M3<br>(*BCI Tech)                                                |    |                   | U403-06       | Adult finger /1.1m  | C3            | No need     |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U403S-06      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U803-06       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| Contec康泰(China)                                                                                                   |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
|                                                                                                                   |    |                   | U410-29       | Adult finger /3.0m  | C7            | No need     |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U410S-29      | Adult soft /3.0m    | C8            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U810-29       | Multi-site Y/3.0m   | C7            |             |                                                                                       |                                                                                       |
| CMS60C/CW/D/DW,CMS60,<br>PM50,PM60D,CSM5000,<br>CMS5000C,CMS5100,<br>CMS6000,CMS7000/<br>8000/9000(*Digital Tech) |    |                   | U410-46D      | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U410S-46D     | Adult soft /3.0m    | D2            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U810-46D      | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| CMS50F/CMS50FW                                                                                                    |    |                   | U401-149      | Adult finger /0.3m  | C7            | No need     |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U401S-149     | Adult soft /0.3m    | C8            |             |                                                                                       |                                                                                       |
| Colin/Omron                                                                                                       |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| BP306,BP88,BP508<br>(Nellcor Non-Oximax)<br>M570                                                                  |    | DS100A            | U403-01       | Adult finger /1.1m  | C3            | U708-01     |    |    |
|                                                                                                                   |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| BP88S<br>(*Non-Oximax)                                                                                            |    | See the above     | See the above | See the above       | See the above | U708-69     |    |    |
|                                                                                                                   |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
|                                                                                                                   |  |                   | U410-69       | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U410S-69      | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U810-69       | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| BX 10NE<br>(* Oximax)                                                                                             |  | DS100A            | U403-01XP     | Adult finger /1.1m  | C5            | U710-70     |  |  |
|                                                                                                                   |                                                                                     | DS100A            | U403S-01XP    | Adult soft /1.1m    | C7            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     | D-YS              | U803-01XP     | Multi-site Y/1.1m   | C5            |             |                                                                                       |                                                                                       |
| BX 10MA<br>(* M tech)                                                                                             |  | 1863/<br>LNCS DCI | U403-49R      | Adult finger /1.1m  | C4            | U708M-15R   |  |  |
|                                                                                                                   |                                                                                     | 1863/<br>LNCS DCI | U403S-49R     | Adult soft /1.1m    | C6            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     | 2258/<br>LNCS YI  | U803-49R      | Multi-site Y/1.1m   | C4            |             |                                                                                       |                                                                                       |
|                                                                                                                   |  | LNOP<br>DC-12     | U410-15R      | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     | LNOP<br>DC-12     | U410S-15R     | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                                                                                   |  |                   | U810-15R      | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| Corpuls                                                                                                           |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| 08/16<br>(*BCI)                                                                                                   |  | 3044              | U403-06       | Adult finger /1.1m  | C3            | NA          |                                                                                       |  |
|                                                                                                                   |                                                                                     |                   | U403S-06      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     | 3043              | U803-06       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| 08/16(M)<br>(* Non-Oximax)                                                                                        |  | DS100A            | U403-01       | Adult finger /1.1m  | C3            | NA          |                                                                                       |  |
|                                                                                                                   |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| 08/16(M)<br>M Tech                                                                                                |  |                   | U410-152M     | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U410S-152M    | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                                                                                   |                                                                                     |                   | U810-152M     | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |

| OEM Manufacturer & Models                                                                                                          | Sensor Connector                                                                      | REF# OEM | REF# UNIMED | Sensor Descriptions | Price Code | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------|-------------|---------------------|------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Comen (China)                                                                                                                      |                                                                                       |          |             |                     |            |             |                                                                                       |                                                                                       |
| Start-8000A/8000B/8000C<br>Star-8000D,8000H<br>Star-5000B/5000C,C50<br>(*Digital Tech)                                             |    |          | U410-66D    | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-66D   | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-66D    | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| Start-8000E                                                                                                                        |    |          | U410-183    | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-183   | Adult soft /3.0m    | C8         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-183    | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| C60,C80<br>(*Oximax Tech)                                                                                                          |    |          | U410-187P   | Adult finger /3.0m  | D1         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-187P  | Adult soft /3.0m    | D2         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-187P   | Multi-site Y/3.0m   | D1         |             |                                                                                       |                                                                                       |
| C100/C100A/C21/C30<br>C50/C70/C80/NC3<br>NC5/NC8/NC12/NC19                                                                         |    |          | U410-208D   | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-208D  | Adult soft/3.0m     | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-208D   | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| Old C60<br>Masimo Rainbow                                                                                                          |    |          | U410-202M   | Adult finger /3.0m  | D2         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-202M  | Adult soft /3.0m    | D3         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-202M   | Multi-site Y/3.0m   | D2         |             |                                                                                       |                                                                                       |
| C20,C60,C80,C90<br>(*Nellcor OxiMax Tech)                                                                                          |    |          | U410-208P   | Adult finger /3.0m  | D1         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-208P  | Adult soft /3.0m    | D2         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-208P   | Multi-site Y/3.0m   | D1         |             |                                                                                       |                                                                                       |
| C20,C60,C80,C90<br>(Masimo Rainbow)                                                                                                |   |          | U410-208MR  | Adult finger /3.0m  | D2         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-208MR | Adult soft /3.0m    | D3         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-208MR  | Multi-site Y/3.0m   | D2         |             |                                                                                       |                                                                                       |
| C60,C80<br>(*Masimo Tech)                                                                                                          |  |          | U410-187M   | Adult finger /3.0m  | D2         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-187M  | Adult soft/3.0m     | D3         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-187M   | Multi-site Y/3.0m   | D2         |             |                                                                                       |                                                                                       |
| Creative(China)                                                                                                                    |                                                                                       |          |             |                     |            |             |                                                                                       |                                                                                       |
| UP-6000,PC-9000<br>(*Non-Oximax)                                                                                                   |  | DS100A   | U403-01     | Adult finger /1.1m  | C3         | U708-29     |  |  |
|                                                                                                                                    |                                                                                       | DS100A   | U403S-01    | Adult soft /1.1m    | C5         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       | D-YS     | U803-01     | Multi-site Y/1.1m   | C3         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |  |          | U410-29     | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-29    | Adult soft /3.0m    | C8         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-29     | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| UP-7000,UP-9000,PC-900,<br>PC-900A,PC-900S<br>PC-900SNPC-3000,<br>PC-5000A/B,PC-9000B,<br>PC-66A, PC-66B,<br>PC-66V(*Digital tech) |  |          | U410-61D    | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U410S-61D   | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                                    |                                                                                       |          | U810-61D    | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |

| OEM Manufacturer & Models                                            | Sensor Connector                                                                    | REF# OEM             | REF# UNIMED | Sensor Descriptions    | Price Code | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|----------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------------------|-------------|------------------------|------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Criticare (CSI)                                                      |                                                                                     |                      |             |                        |            |             |                                                                                       |                                                                                       |
| 510,503,504,<br>507,508,581,<br>1100,8100,Poet,<br>Scholar<br>(CSI*) |    | 934-10LN             | U410-05     | Adult finger /3.0m     | C8         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U110-05     | Pediatric finger /3.0m | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     | 975AD-10             | U410S-05    | Adult soft /3.0m       | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U110S-05    | Pediatric soft /3.0m   | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U210S-05    | Infant soft /3.0m      | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U310-05     | Neonate wrap /3.0m     | D1         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U610-05     | Adult/Neo wrap /3.0m   | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U610B-05    | Pedia/infant wrap/3.0m | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U810-05     | Multi-site Y/3.0m      | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U910-05     | Adult ear clip /3.0m   | D1         |             |                                                                                       |                                                                                       |
|                                                                      |    |                      | U403-05     | Adult finger /1.1m     | C3         | U708-05     |    |    |
|                                                                      |                                                                                     |                      | U103-05     | Pediatric finger /1.1m | C4         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U403S-05    | Adult soft /1.1m       | C4         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U103S-05    | Pediatric soft /1.1m   | C4         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U203S-05    | Infant soft /1.1m      | C4         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U303-75     | Neonate wrap /1.1m     | C4         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U603-05     | Multi-site Y/1.1m      | C3         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U603B-05    | Adult/Neo/1.1m         | C3         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U803-05     | Pedia/infant wrap/1.1m | C4         |             |                                                                                       |                                                                                       |
|                                                                      | U903-05                                                                             | Adult ear clip /1.1m | C4          |                        |            |             |                                                                                       |                                                                                       |
| 503SPOT,504Plus,<br>504DX,506DX<br>(*CSI)                            |  | 934-10DN             | U410-75     | Adult finger /3.0m     | C7         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U110-75     | Pediatric finger /3.0m | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U410S-75    | Adult soft /3.0m       | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U110S-75    | Pediatric soft /3.0m   | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U210S-75    | Infant soft /3.0m      | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U310-75     | Neonate wrap /3.0m     | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U610-75     | Adult/Neo wrap/3.0m    | C7         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U610B-75    | Pedia/infant wrap/3.0m | C7         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U810-75     | Multi-site Y/3.0m      | C7         |             |                                                                                       |                                                                                       |
|                                                                      | U910-75                                                                             | Adult ear clip /3.0m | C8          |                        |            |             |                                                                                       |                                                                                       |
| 503SPOT,503DX,504Plus,<br>504DX,506DX,506DXNP<br>(*CSI)              |  |                      | U410-75N    | Adult finger/3.0m      | C7         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U410S-75N   | Adult soft/3.0m        | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U810-75N    | Multi-site Y/3.0m      | C7         |             |                                                                                       |                                                                                       |
| Digicare                                                             |                                                                                     |                      |             |                        |            |             |                                                                                       |                                                                                       |
| Digimax5500<br>(*BCI)                                                |  | 3044                 | U403-06     | Adult finger /1.1m     | C3         | U708-76     |  |  |
|                                                                      |                                                                                     |                      | U403S-06    | Adult soft /1.1m       | C5         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     | 3043                 | U803-06     | Multi-site Y/1.1m      | C3         |             |                                                                                       |                                                                                       |
|                                                                      |  |                      | U410-76     | Adult finger /3.0m     | C7         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U410S-76    | Adult soft /3.0m       | C8         |             |                                                                                       |                                                                                       |
|                                                                      |  |                      | U810-76     | Multi-site Y/3.0m      | C7         | No need     |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U410-52     | Adult finger /3.0m     | C7         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U410S-52    | Adult soft /3.0m       | C8         |             |                                                                                       |                                                                                       |
|                                                                      |                                                                                     |                      | U810-52     | Multi-site Y/3.0m      | C7         |             |                                                                                       |                                                                                       |





| OEM Manufacturer & Models                                    | Sensor Connector                                                                      | REF# OEM                                                                              | REF# UNIMED          | Sensor Descriptions    | Price Code         | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |                                                                                       |
|--------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------------------|------------------------|--------------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Datascope                                                    |                                                                                       |                                                                                       |                      |                        |                    |             |                                                                                       |                                                                                       |                                                                                       |
| Accusat, Accutor 3/4 sat Passport, Accutor plus (*Datascope) |    |                                                                                       | U403-07              | Adult finger /1.1m     | C4                 | U708-07     |    |    |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U103-07              | Pediatric finger /1.1m | C6                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U403S-07             | Adult soft /1.1m       | C6                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U103S-07             | Pediatric soft /1.1m   | C6                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U203S-07             | Infant soft /1.1m      | C6                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U303-07              | Neonate wrap /1.1m     | C6                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U603-07              | Adult/Neo wrap/1.1m    | C5                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U603B-07             | Pedia/infant wrap/1.1m | C5                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U803-07              | Multi-site Y/1.1m      | C5                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       | U903-07                                                                               | Adult ear clip /1.1m | C4                     |                    |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |    |                                                                                       | U410-07              | Adult finger /3.0m     | C8                 | No need     |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U410S-07             | Adult soft /3.0m       | D1                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       | U810-07                                                                               | Multi-site Y/3.0m    | C8                     |                    |             |                                                                                       |                                                                                       |                                                                                       |
| Accutor plus (*M Tech)                                       |    | 1863/ LNCS DCI                                                                        | U403-49R             | Adult finger /1.1m     | C4                 | U708M-15R   |    |    |                                                                                       |
|                                                              |                                                                                       | 1863/ LNCS DCI                                                                        | U403S-49R            | Adult soft /1.1m       | C6                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       | 2258/ LNCS YI                                                                         | U803-49R             | Multi-site Y/1.1m      | C4                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |   | LNOP DC-12                                                                            | U410-15              | Adult finger /3.0m     | C8                 | No need     |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       | LNOP DC-12                                                                            | U410S-15             | Adult soft /3.0m       | D1                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U810-15              | Multi-site Y/3.0m      | C8                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              | PassportII, Expert Pssport XG, Accutor plus (*Non-Oximax)                             |  | DS100A               | U403-01                | Adult finger /1.1m | C3          | U708-69                                                                               |  |  |
|                                                              |                                                                                       |                                                                                       | DS100A               | U403S-01               | Adult soft /1.1m   | C5          |                                                                                       |                                                                                       |                                                                                       |
| D-YS                                                         |                                                                                       |                                                                                       | U803-01              | Multi-site Y/1.1m      | C3                 |             |                                                                                       |                                                                                       |                                                                                       |
| Accutor-Trio (*Non-Oximax)                                   |  | See the above                                                                         | See the above        | See the above          | See the above      | U708-01     |  |  |                                                                                       |
| Trio Nouvelle Generation (* Oximax)                          |  | DS100A                                                                                | U403-01XP            | Adult finger /1.1m     | C5                 | U710-70     |  |  |                                                                                       |
|                                                              |                                                                                       | DS100A                                                                                | U403S-01XP           | Adult soft /1.1m       | C7                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       | D-YS                                                                                  | U803-01XP            | Multi-site Y/1.1m      | C5                 |             |                                                                                       |                                                                                       |                                                                                       |
| Dima (Italia)                                                |                                                                                       |                                                                                       |                      |                        |                    |             |                                                                                       |                                                                                       |                                                                                       |
| (*Non-Oximax)                                                |  | DS100A                                                                                | U403-01              | Adult finger /1.1m     | C3                 | U708-01     |  |  |                                                                                       |
|                                                              |                                                                                       | DS100A                                                                                | U403S-01             | Adult soft /1.1m       | C5                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       | D-YS                                                                                  | U803-01              | Multi-site Y/1.1m      | C3                 |             |                                                                                       |                                                                                       |                                                                                       |
| Dixon                                                        |                                                                                       |                                                                                       |                      |                        |                    |             |                                                                                       |                                                                                       |                                                                                       |
| D6 D8                                                        |  |                                                                                       | U410-205             | Adult finger /3.0m     | C8                 | No need     |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U410S-205            | Adult soft/3.0m        | D1                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U810-205             | Multi-site Y/3.0m      | C8                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |  |                                                                                       | U403-62D             | Adult finger /1.1m     | C8                 | U708-208    |  |  |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U403S-62D            | Adult soft/1.1m        | D1                 |             |                                                                                       |                                                                                       |                                                                                       |
|                                                              |                                                                                       |                                                                                       | U803-62D             | Multi-site Y/1.1m      | C8                 |             |                                                                                       |                                                                                       |                                                                                       |
| Dixon                                                        |                                                                                       |                                                                                       |                      |                        |                    |             |                                                                                       |                                                                                       |                                                                                       |


| OEM Manufacturer & Models              | Sensor Connector                                                                    | REF# OEM      | REF# UNIMED   | Sensor Descriptions | Price Code        | Ext - Cable        | Cable Plug                                                                            | Sensor Conn.                                                                          |
|----------------------------------------|-------------------------------------------------------------------------------------|---------------|---------------|---------------------|-------------------|--------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Dixon                                  |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
| D2,D3,D5,D6,D8                         |    |               | U410-205D     | Adult finger /3.0m  | C8                | No need            |                                                                                       |                                                                                       |
| DS3,DS5,DS7                            |                                                                                     |               | U410S-205D    | Adult soft /3.0m    | D1                |                    |                                                                                       |                                                                                       |
| (Digital)                              |                                                                                     |               | U810-205D     | Multi-site Y/3.0m   | C8                |                    |                                                                                       |                                                                                       |
| M800 handled vet                       |    |               | U403-62D      | Adult finger /1.1m  | C8                | No need            |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U403S-62D     | Adult soft /1.1m    | D1                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U803-62D      | Multi-site Y/1.1m   | C8                |                    |                                                                                       |                                                                                       |
| DS3,DS5,DS7                            |    |               | U410-87       | Adult finger /3.0m  | D2                | No need            |                                                                                       |                                                                                       |
| (Masimo Rainbow)                       |                                                                                     |               | U410S-87      | Adult soft /3.0m    | D3                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U810-87       | Multi-site Y/3.0m   | D2                |                    |                                                                                       |                                                                                       |
| Storm 5000                             |    |               | U403-01       | Adult finger /1.1m  | C3                | No need            |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U403S-01      | Adult soft /1.1m    | C5                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U803-01       | Multi-site Y/1.1m   | C3                |                    |                                                                                       |                                                                                       |
| Storm 5300,Storm 5500                  |    |               | U410-29       | Adult finger /3.0m  | C7                | No need            |                                                                                       |                                                                                       |
| Storm 5800,Storm 5900                  |                                                                                     |               | U410S-29      | Adult soft /3.0m    | C8                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U810-29       | Multi-site Y/3.0m   | C7                |                    |                                                                                       |                                                                                       |
| Dixtal(Brazil)-Philips                 |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
| Full compatible with Novamatrix Sensas |    | DX-1NC01-0    | U410-113      | Adult finger /3.0m  | C7                | No need            |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U410S-113     | Adult soft /3.0m    | C8                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U810-113      | Multi-site Y/3.0m   | C7                |                    |                                                                                       |                                                                                       |
| Dolphin                                |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
| (*Digital Tech)                        |  |               | U403-50D      | Adult finger /1.1m  | C8                | U708-50            |  |  |
|                                        |                                                                                     |               | U403S-50D     | Adult soft /1.1m    | D1                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U803-50D      | Multi-site Y/1.1m   | C8                |                    |                                                                                       |                                                                                       |
| (*Digital Tech)                        |  |               | U410-50D      | Adult finger /3.0m  | D1                | No need            |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U410S-50D     | Adult soft /3.0m    | D2                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U810-50D      | Multi-site Y/3.0m   | D1                |                    |                                                                                       |                                                                                       |
| Drager/Siemens                         |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
| Oxipac 2000 (*BCI)                     |  | 3044          | U403-06       | Adult finger /1.1m  | C3                | No need            |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | U403S-06      | Adult soft /1.1m    | C5                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               | 3043          | U803-06             | Multi-site Y/1.1m |                    |                                                                                       |                                                                                       |
| Vitalert 1000,icero (*Non-Oximax)      |  | DS100A        | U403-01       | Adult finger /1.1m  | C3                | U708-28            |  |  |
|                                        |                                                                                     | DS100A        | U403S-01      | Adult soft /1.1m    | C5                |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     | D-YS          | U803-01       | Multi-site Y/1.1m   | C3                |                    |                                                                                       |                                                                                       |
| Narkomed-Oxysat 2 (*Non-Oximax)        |  | See the above | See the above | See the above       | See the above     | U710-85            |  |  |
|                                        |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
| Sirecust, Sirem Spo2 (*Non-Oximax)     |  | See the above | See the above | See the above       | See the above     | U708-33            |  |  |
|                                        |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
| SC6000, 7000,9000 (*Non-Oximax)        |  | See the above | See the above | See the above       | See the above     | U704-23<br>U708-23 |  |  |
|                                        |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
|                                        |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |
| Drager/Siemens                         |                                                                                     |               |               |                     |                   |                    |                                                                                       |                                                                                       |


| OEM Manufacturer & Models                                                                                            | Sensor Connector                                                                      | REF# OEM       | REF# UNIMED | Sensor Descriptions | Price Code | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|----------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------------|-------------|---------------------|------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Drager/Siemens                                                                                                       |                                                                                       |                |             |                     |            |             |                                                                                       |                                                                                       |
| Oxytrend (*M tech))                                                                                                  |    | 1863/ LNCS DCI | U403-49R    | Adult finger /1.1m  | C4         | U708M-15R   |    |    |
|                                                                                                                      |                                                                                       | 1863/ LNCS DCI | U403S-49R   | Adult soft /1.1m    | C6         |             |                                                                                       |                                                                                       |
|                                                                                                                      |    | 2258/ LNCS YI  | U803-49R    | Multi-site Y/1.1m   | C4         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       | LNOP DC-12     | U410-15R    | Adult finger /3.0m  | C8         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       | LNOP DC-12     | U410S-15R   | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-15R    | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| (*Oximax Tech)                                                                                                       |    |                | U410-23P    | Adult finger /3.0m  | D1         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-23P   | Adult soft /3.0m    | D2         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-23P    | Multi-site Y/3.0m   | D1         |             |                                                                                       |                                                                                       |
| (*M Tech)                                                                                                            |    |                | U410-23MR   | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-23MR  | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-23MRT  | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| Sirecust 1281 Monitor<br>SC 6002XL,SC7000,<br>SC9000XL,<br>( Masimo Rainbow)                                         |    |                | U410-23MR   | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-23MR  | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-23MR   | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| DRE                                                                                                                  |                                                                                       |                |             |                     |            |             |                                                                                       |                                                                                       |
| DRE Waveline EZ<br>Touch<br>Pro                                                                                      |    |                | U410-213    | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-213   | Adult soft /3.0m    | C8         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-213    | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| Waveline EZ MAX,<br>Waveline Pro<br>(*Oximax Tech)                                                                   |  |                | U410-212P   | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-212P  | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-212P   | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| Waveline,<br>Waveline Pro<br>Omni compatible, One                                                                    |  |                | U410-137    | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-137   | Adult soft /3.0m    | C8         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-137    | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
| Edan (China)                                                                                                         |                                                                                       |                |             |                     |            |             |                                                                                       |                                                                                       |
| M9,M99A,M8,M8A,<br>M8B, M3,M3B<br>(*Oximax)                                                                          |  | DS100A         | U403-01XP   | Adult finger /1.1m  | C5         | U708X-191   |  |  |
|                                                                                                                      |                                                                                       | DS100A         | U403S-01XP  | Adult soft /1.1m    | C7         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       | D-YS           | U803-01XP   | Multi-site Y/1.1m   | C5         |             |                                                                                       |                                                                                       |
| EDAN Oximetry<br>M3,M3A,M3B,<br>M8,M8A,M8B<br>M9,M9A,M50,M80<br>iM8,iM8A,iM8B<br>iM9,iM9A,iM50,iM80<br>(*Non-Oximax) |  | MS3-109069     | U410-29     | Adult finger /3.0m  | C7         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-29    | Adult soft /3.0m    | C8         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-29     | Multi-site Y/3.0m   | C7         |             |                                                                                       |                                                                                       |
|                                                                                                                      |  | MS-109079      | U403-01     | Adult finger /1.1m  | C3         | U708-29     |  |  |
|                                                                                                                      |                                                                                       | DS100A         | U403S-01    | Adult soft /1.1m    | C5         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       | D-YS           | U803-01     | Multi-site Y/1.1m   | C3         |             |                                                                                       |                                                                                       |
| M3,M3B,M8,M8A,M8B<br>M9,M9A,M50,M80<br>iM9,iM9A,iM50,iM80<br>(*Oximax)                                               |  |                | U410-191P   | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-191P  | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-191P   | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| (*Digital tech)                                                                                                      |  |                | U410-46D    | Adult finger /3.0m  | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U410S-46D   | Adult soft /3.0m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       |                | U810-46D    | Multi-site Y/3.0m   | C8         |             |                                                                                       |                                                                                       |
| H100B<br>(*Non Oximax)                                                                                               |  | MS-109079      | U403-01     | Adult finger /1.1m  | C3         | No need     |                                                                                       |                                                                                       |
|                                                                                                                      |                                                                                       | DS100A         | U403S-01    | Adult soft /1.1m    | C5         |             |                                                                                       |                                                                                       |
| Edan(China)                                                                                                          |                                                                                       |                |             |                     |            |             |                                                                                       |                                                                                       |


| OEM Manufacturer & Models                                                  | Sensor Connector                                                                    | REF# OEM          | REF# UNIMED   | Sensor Descriptions | Price Code         | Ext - Cable                | Cable Plug                                                                            | Sensor Conn.                                                                          |  |
|----------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------|---------------|---------------------|--------------------|----------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--|
| Edan (China)                                                               |                                                                                     |                   |               |                     |                    |                            |                                                                                       |                                                                                       |  |
| H100N<br>(*Oximax)                                                         |    | DS100A            | U403-01XP     | Adult finger /1.1m  | C5                 | No need                    |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     | DS100A            | U403S-01XP    | Adult soft /1.1m    | C7                 |                            |                                                                                       |                                                                                       |  |
| iT20,iM20,<br>elite V5/V6/V8                                               |    |                   | U410-48       | Adult finger /3.0m  | C7                 | No need                    |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   | U410S-48      | Adult soft /3.0m    | C8                 |                            |                                                                                       |                                                                                       |  |
| iT20,iM20,<br>elite V5/V6/V8<br>(*Nellcor OxiMax)                          |    |                   | U810-48       | Multi-site Y/3.0m   | C7                 | No need                    |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   | U410-48P      | Adult finger /3.0m  | C8                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   | U410S-48P     | Adult soft/3.0m     | D1                 |                            |                                                                                       |                                                                                       |  |
| iT20,iM20,<br>elite V5/V6/V8<br>(Masimo Rainbow)                           |    |                   | U810-48P      | Multi-site Y/3.0m   | C8                 | No need                    |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   | U410-48MR     | Adult finger /3.0m  | D1                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   | U410S-48MR    | Adult soft/3.0m     | D2                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   | U810-48MR     | Multi-site Y/3.0m   | D1                 |                            |                                                                                       |                                                                                       |  |
| Fukuda Denshi                                                              |                                                                                     |                   |               |                     |                    |                            |                                                                                       |                                                                                       |  |
| DS7100-DS7300<br>(*Oximax)                                                 |    | DS100A            | U403-01XP     | Adult finger /1.1m  | C5                 | U710-70                    |    |    |  |
|                                                                            |                                                                                     | DS100A            | U403S-01XP    | Adult soft /1.1m    | C7                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     | D-YS              | U803-01XP     | Multi-site Y/1.1m   | C5                 |                            |                                                                                       |                                                                                       |  |
| GE Critikon                                                                |                                                                                     |                   |               |                     |                    |                            |                                                                                       |                                                                                       |  |
| Dinamap 9700/9701<br>9702 PLUS, MPS<br>ORTABLE,/COMPACT,<br>(* Non-Oximax) |    | DS100A            | U403-01       | Adult finger /1.1m  | C3                 | U708-01<br>(EC-8)          |    |    |  |
|                                                                            |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3                 |                            |                                                                                       |                                                                                       |  |
| Dinamap,ProSeries V1<br>Pro 1000 V2<br>(Nellcor Optional)                  |  | See the above     | See the above | See the above       | See the above      | U708-69<br>(SCP-10)        |  |  |  |
| Dinamap,ProCare,<br>ProSeries V2Pro1000 V3<br>(*Oximax)                    |  | DS100A            | U403-01XP     | Adult finger /1.1m  | C5                 | U710-70<br>(DOC-10)        |  |  |  |
|                                                                            |                                                                                     | DS100A            | U403S-01XP    | Adult soft /1.1m    | C7                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     | D-YS              | U803-01XP     | Multi-site Y/1.1m   | C5                 |                            |                                                                                       |                                                                                       |  |
| Dinamap<br>8700/8710/8720Pulse,<br>Oxyshuttle11 et 2                       |  | DS100A            | U403-01       | Adult finger /1.1m  | C3                 | Not Available              |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |  |                   | U410-10       | Adult finger /3.0m  | C8                 | No need                    |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   | U410S-10      | Adult soft /3.0m    | D1                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   | U810-10       | Multi-site Y/3.0m   | C8                 |                            |                                                                                       |                                                                                       |  |
| Dinamap<br>ProCare<br>Series V2<br>Pro 1000 V3<br>(*M Tech)                |  | 1863/<br>LNCS DCI | U403-49R      | Adult finger /1.1m  | C4                 | U708M-15R<br>(2027263-001) |  |  |  |
|                                                                            |                                                                                     | 1863/<br>LNCS DCI | U403S-49R     | Adult soft /1.1m    | C6                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     | 2258/<br>LNCS YI  | U803-49R      | Multi-site Y/1.1m   | C4                 |                            |                                                                                       |                                                                                       |  |
|                                                                            |  |                   | LNOP<br>DC-12 | U410-15R            | Adult finger /3.0m | C8                         | No need                                                                               |                                                                                       |  |
|                                                                            |                                                                                     |                   | LNOP<br>DC-12 | U410S-15R           | Adult soft /3.0m   | D1                         |                                                                                       |                                                                                       |  |
|                                                                            |                                                                                     |                   |               | U810-15R            | Multi-site Y/3.0m  | C8                         |                                                                                       |                                                                                       |  |




| OEM Manufacturer & Models                                                                                                     | Sensor Connector                                                                      | REF# OEM             | REF# UNIMED          | Sensor Descriptions    | Price Code    | Ext - Cable       | Cable Plug                                                                            | Sensor Conn.                                                                          |
|-------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------------------|----------------------|------------------------|---------------|-------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| GE Corometrics                                                                                                                |                                                                                       |                      |                      |                        |               |                   |                                                                                       |                                                                                       |
| 188,129,129F, 128,128F (*Non Oximax)                                                                                          |    | DS100A               | U403-01              | Adult finger /1.1m     | C3            | U708-01           |    |    |
|                                                                                                                               |                                                                                       | DS100A               | U403S-01             | Adult soft /1.1m       | C5            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       | D-YS                 | U803-01              | Multi-site Y/1.1m      | C3            |                   |                                                                                       |                                                                                       |
| 120,510,511, 118,556 (*Non Oximax)                                                                                            |    | See the above        | See the above        | See the above          | See the above | U708-25 (4033CAX) |    |    |
|                                                                                                                               |                                                                                       |                      |                      |                        |               |                   |                                                                                       |                                                                                       |
| 128, 129 (*Non Oximax)                                                                                                        |    | See the above        | See the above        | See the above          | See the above | U708-21           |    |    |
|                                                                                                                               |                                                                                       |                      |                      |                        |               |                   |                                                                                       |                                                                                       |
|                                                                                                                               |    |                      | U410-21              | Adult finger /3.0m     | C8            | No need           |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U410S-21             | Adult soft /3.0m       | D1            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U810-21              | Multi-site Y/3.0m      | C8            |                   |                                                                                       |                                                                                       |
| SERIE 250 (*Oximax)                                                                                                           |    | DS100A               | U403-01XP            | Adult finger /1.1m     | C5            | U710X-21          |    |    |
|                                                                                                                               |                                                                                       | DS100A               | U403S-01XP           | Adult soft /1.1m       | C7            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       | D-YS                 | U803-01XP            | Multi-site Y/1.1m      | C5            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |    | See Oximax Marquette | See Oximax Marquette | See Oximax Marquette   | See GE        | No need           |                                                                                       |                                                                                       |
| GE Datex-Ohmeda                                                                                                               |                                                                                       |                      |                      |                        |               |                   |                                                                                       |                                                                                       |
| AS/3,CS/3, AS/5 (Built in M module) CH/S,CH/RSCG/2S, CARDIOCAPI-II Capnomac Ultima Light monitor Satlite, Oscar (Datex tech*) |  | OXY-F4-N             | U410-09              | Adult finger /3.0m     | C8            | No need           |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U110-09              | Pediatric finger /3.0m | D1            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U410S-09             | Adult soft /3.0m       | D1            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U110S-09             | Pediatric soft /3.0m   | D1            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U210S-09             | Infant soft /3.0m      | D1            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U310-09              | Neonate wrap /3.0m     | D1            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       | OXY-W4-N             | U610-09              | Adult/Neo wrap /3.0m   | C8            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U610B-09             | Pedia/infant wrap/3.0m | C8            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       | OXY-W4-N             | U810-09              | Multi-site Y/3.0m      | C8            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       | OXY-E4-N             | U910-09              | Adult ear clip /3.0m   | D1            |                   |                                                                                       |                                                                                       |
| M-NSAT Module (AS/3,CS/3,S/5) (* Non Oximax)                                                                                  |  | DS100A               | U403-01              | Adult finger /1.1m     | C3            | U708-69 (MC-10)   |  |  |
|                                                                                                                               |                                                                                       | DS100A               | U403S-01             | Adult soft /1.1m       | C5            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       | D-YS                 | U803-01              | Multi-site Y/1.1m      | C3            |                   |                                                                                       |                                                                                       |
| M-NSAT Module (i4) (*Oximax)                                                                                                  |  | DS100A               | U403-01XP            | Adult finger /1.1m     | C5            | U710-70 (DOC-10)  |  |  |
|                                                                                                                               |                                                                                       | DS100A               | U403S-01XP           | Adult soft /1.1m       | C7            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       | D-YS                 | U803-01XP            | Multi-site Y/1.1m      | C5            |                   |                                                                                       |                                                                                       |
| AS/3,CS/3, AS/5 (Built in M module) CH/S,CH/RSCG/2S, CARDIOCAPI-II Capnomac Ultima Light monitor Satlite, Oscar (Datex tech*) |  | OXY-F-DB             | U403-09              | Adult finger /1.1m     | C3            | U708-09           |  |  |
|                                                                                                                               |                                                                                       |                      | U103-09              | Pediatric finger /1.1m | C5            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U403S-09             | Adult soft /1.1m       | C5            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U103S-09             | Pediatric soft /1.1m   | C5            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U203S-09             | Infant soft /1.1m      | C5            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      | U303-09              | Neonate wrap /1.1m     | C5            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       | OXY-W-DB             | U603-09              | Adult/Neo wrap/1.1m    | C3            |                   |                                                                                       |                                                                                       |
|                                                                                                                               |                                                                                       |                      |                      |                        |               |                   |                                                                                       |                                                                                       |
| GE Datex-Ohmeda                                                                                                               |                                                                                       |                      |                      |                        |               |                   |                                                                                       |                                                                                       |


| OEM Manufacturer & Models                                                                                                                                   | Sensor Connector                                                                    | REF# OEM             | REF# UNIMED       | Sensor Descriptions    | Price Code        | Ext - Cable        | Cable Plug                                                                            | Sensor Conn.                                                                          |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------------------|-------------------|------------------------|-------------------|--------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| GE Datex-Ohmeda                                                                                                                                             |                                                                                     |                      |                   |                        |                   |                    |                                                                                       |                                                                                       |
| 3700,3710,3770, 3775,Biox3740,3760, 3800,3900,3900P, 4700,5250RGB, Cardiocap/5, S/5, AS/3 (built in Module N-XOSAT) CS/3 (Built in module M-OSAT) (*Ohmeda) |    |                      | U410-02           | Adult finger /3.0m     | C8                | No need            |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U110-02           | Pediatric finger /3.0m | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U410S-02          | Adult soft /3.0m       | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U110S-02          | Pediatric soft /3.0m   | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U210S-02          | Infant soft /3.0m      | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U310-02           | Neonate wrap /3.0m     | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U610-02           | Adult/Neo wrap /3.0m   | C8                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U610B-02          | Pedia/infant wrap/3.0m | C8                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U810-02           | Multi-site Y/3.0m      | C8                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             | U910-02                                                                             | Adult ear clip /3.0m | D1                |                        |                   |                    |                                                                                       |                                                                                       |
| Ohmeda Tuffsat (*Ohmeda)                                                                                                                                    |    |                      | U403-02           | Adult finger /1.1m     | C8                | No need            |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U403S-02          | Adult soft /1.1m       | D1                |                    |                                                                                       |                                                                                       |
| Ohmeda TruSat<br>Ohmeda TruSat TD (*+Oxytip+ Technology)                                                                                                    |    |                      | U410-60           | Adult finger /3.0m     | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U110-60           | Pediatric finger /3.0m | D2                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U410S-60          | Adult soft /3.0m       | D2                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U110S-60          | Pediatric soft /3.0m   | D2                |                    |                                                                                       |                                                                                       |
|                                                                            |    |                      | U210S-60          | Infant soft /3.0m      | D2                | No need            |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U310-60           | Neonate wrap /3.0m     | D2                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U610-60           | Adult/Neo wrap /3.0m   | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U610B-60          | Pedia/infant wrap/3.0m | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U810-60           | Multi-site Y/3.0m      | D1                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U910-60           | Adult ear clip /3.0m   | D2                |                    |                                                                                       |                                                                                       |
|                                                                          |  |                      | U603B-09          | Pedia/infant wrap/1.1m | C3                | U708-09            |  |  |
|                                                                                                                                                             |                                                                                     |                      | U803-09           | Multi-site Y/1.1m      | C3                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     | OXY-E-DB             | U903-09           | Adult ear clip /1.1m   | C5                |                    |                                                                                       |                                                                                       |
| GE Hellige                                                                                                                                                  |                                                                                     |                      |                   |                        |                   |                    |                                                                                       |                                                                                       |
| Cardioserv<br>Eagle1000-3000 (*Non Oximax)                                                                                                                  |  | DS100A               | U403-01           | Adult finger /1.1m     | C3                | U708-55 (30344358) |  |  |
|                                                                                                                                                             |                                                                                     | DS100A               | U403S-01          | Adult soft /1.1m       | C5                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     | D-YS                 | U803-01           | Multi-site Y/1.1m      | C3                |                    |                                                                                       |                                                                                       |
| Cardioserv (* Oximax)                                                                                                                                       |  | DS100A               | U403-01XP         | Adult finger /1.1m     | C5                | NA (2026751-001)   |  |  |
|                                                                                                                                                             |                                                                                     | DS100A               | U403S-01XP        | Adult soft /1.1m       | C7                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     | D-YS                 | U803-01XP         | Multi-site Y/1.1m      | C5                |                    |                                                                                       |                                                                                       |
| SMU,SMK,SMV (*Ohmeda)                                                                                                                                       |  | See above *Ohmeda    | See above *Ohmeda | See above *Ohmeda      | See above *Ohmeda | No need            |                                                                                       |                                                                                       |
| GE Marquette                                                                                                                                                |                                                                                     |                      |                   |                        |                   |                    |                                                                                       |                                                                                       |
| Dash3000/4000/5000<br>Solar8000/8000M<br>Solar 8000i/9500<br>B40i/ 60i/B650 (M Tech)                                                                        |  |                      | U410-21MR         | Adult finger /3.0m     | D2                | No need            |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U410S-21MR        | Adult soft /3.0m       | D3                |                    |                                                                                       |                                                                                       |
|                                                                                                                                                             |                                                                                     |                      | U810-21MR         | Multi-site Y/3.0m      | D2                |                    |                                                                                       |                                                                                       |
| 7000/7010,Tram 200/400/600/800 (*Ohmeda Tech.)                                                                                                              |  | See above *Ohmeda    | See above *Ohmeda | See above *Ohmeda      | See above *Ohmeda |                    |                                                                                       |                                                                                       |
| GE Marquette                                                                                                                                                |                                                                                     |                      |                   |                        |                   |                    |                                                                                       |                                                                                       |


| OEM Manufacturer & Models                                                                                                                                                                                        | Sensor Connector                                                                      | REF# OEM       | REF# UNIMED | Sensor Descriptions    | Price Code | Ext - Cable             | Cable Plug                                                                            | Sensor Conn.                                                                          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------------|-------------|------------------------|------------|-------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| GE Marquette                                                                                                                                                                                                     |                                                                                       |                |             |                        |            |                         |                                                                                       |                                                                                       |
| Eagle 1000/3000/4000<br>Dash 2000/3000/4000<br>Solar 8000/8000M<br>i/9500(Tram 451/451N /851/851N)<br>CardioLab/Mac-Lab, Corometrics 120/2120 (Tran 451/451N)<br>Responder3000 (nellcor optional) (* Non Oximax) |    | DS100A         | U403-01     | Adult finger /1.1m     | C3         | U708-21 (2006644-001)   |    |    |
|                                                                                                                                                                                                                  |                                                                                       | DS100A         | U403S-01    | Adult soft /1.1m       | C5         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       | D-YS           | U803-01     | Multi-site Y/1.1m      | C3         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |    |                | U410-21     | Adult finger /3.0m     | C8         | No need                 |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U110-21     | Pediatric finger /3.0m | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U410S-21    | Adult soft /3.0m       | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U110S-21    | Pediatric soft /3.0m   | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U210S-21    | Infant soft /3.0m      | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U310-21     | Neonate wrap /3.0m     | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U610-21     | Adult/Neo wrap/3.0m    | C8         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U610B-21    | Pedia/infant wrap/3.0m | C8         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U810-21     | Multi-site Y/3.0m      | C8         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U910-21     | Adult ear clip /3.0m   | D1         |                         |                                                                                       |                                                                                       |
| Dash2500/Dash 3000/ 4000/5000<br>Solar 8000/8000M<br>i/9500 (Tram 451N5/851N5)<br>Corometrics 250 (GE-Oximax*)<br>Ending “X” can be replaced by “P” for fully compatible with all Oximax oximetry                |    | DS100A         | U403-01XP   | Adult finger /1.1m     | C5         | U710X-21 (2021406-001)  |    |    |
|                                                                                                                                                                                                                  |                                                                                       | DS100A         | U403S-01XP  | Adult soft /1.1m       | C7         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       | D-YS           | U803-01XP   | Multi-site Y/1.1m      | C5         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |  |                | U410-21P    | Adult finger /3.0m     | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U110-21P    | Pediatric finger /3.0m | D2         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U410S-21P   | Adult soft /3.0m       | D2         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U110S-21P   | Pediatric soft /3.0m   | D2         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U210S-21P   | Infant soft /3.0m      | D2         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U310-21P    | Neonate wrap /3.0m     | D2         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U610-21P    | Adult/Neo wrap/3.0m    | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U610B-21P   | Pedia/infant wrap/3.0m | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U810-21P    | Multi-site Y/3.0m      | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U910-21P    | Adult ear clip /3.0m   | D2         |                         |                                                                                       |                                                                                       |
| Dash 3000, 4000,5000<br>Solar 8000/8000M<br>i/9500 (Tram 451M,851M)<br>Cardiolab/Mac-Lab (Tram 451M)<br>Corometrics 259 (*M tech)                                                                                |  | 1863/ LNCS DCI | U403-49     | Adult finger /1.1m     | C4         | U708M-21R (2027263-002) |  |  |
|                                                                                                                                                                                                                  |                                                                                       | 1863/ LNCS DCI | U403S-49    | Adult soft /1.1m       | C6         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       | 2258/ LNCS YI  | U803-49     | Multi-site Y/1.1m      | C4         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |  |                | U410-21MR   | Adult finger /3.0m     | C8         | No need                 |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U110-21MR   | Pediatric finger /3.0m | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U410S-21MR  | Adult soft /3.0m       | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U110S-21MR  | Pediatric soft /3.0m   | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U210S-21MR  | Infant soft /3.0m      | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U310-21MR   | Neonate wrap /3.0m     | D1         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U610-21MR   | Adult/Neo wrap/3.0m    | C8         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U610B-21MR  | Pedia/infant wrap/3.0m | C8         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U810-21MR   | Multi-site Y/3.0m      | C8         |                         |                                                                                       |                                                                                       |
|                                                                                                                                                                                                                  |                                                                                       |                | U910-21MR   | Adult ear clip /3.0m   | D1         |                         |                                                                                       |                                                                                       |


| OEM Manufacturer & Models                                                                                                                                                                                                              | Sensor Connector                                                                    | REF# OEM  | REF# UNIMED | Sensor Descriptions    | Price Code | Ext - Cable                                     | Cable Plug                                                          | Sensor Conn.                                                                          |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------|-------------|------------------------|------------|-------------------------------------------------|---------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| GE Ohmeda OxyTip+                                                                                                                                                                                                                      |                                                                                     |           |             |                        |            |                                                 |                                                                     |                                                                                       |
| GE S/5<br>(i4 E series module built in)                                                                                                                                                                                                |    | OXY-F4-GE | U410-121    | Adult finger /3.0m     | D1         | No need                                         |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U110-121    | Pediatric finger /3.0m | D2         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U410S-121   | Adult soft /3.0m       | D2         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U110S-121   | Pediatric soft /3.0m   | D2         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U210S-121   | Infant soft /3.0m      | D2         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U310-121    | Neonate wrap /3.0m     | D2         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U610-121    | Adult/Neo wrap/3.0m    | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U610B-121   | Pedia/infant wrap/3.0m | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U810-121    | Multi-site Y/3.0m      | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U910-121    | Adult ear clip /3.0m   | D2         |                                                 |                                                                     |                                                                                       |
| Ohmeda TruSat<br>3755,37XX series<br>3800,3900,3900P<br>4700, 5250<br>Cardiocap/5, S/5,AS/3<br>(N-XOSAT or M built in)<br>CS/3 (M-OSAT built in)<br>S/5 light, Satlite Puls II<br>Capnomac UltimaS/5<br>,Coro 259<br>(*TruSignal tech) |    | OXY-F-UN  | U403-127    | Adult finger /1.1m     | C8         | Cables required depending on the monitor models | Depends on the monitor models<br>See the page of Spo2 adapter cable |    |
|                                                                                                                                                                                                                                        |                                                                                     |           | U103-127    | Pediatric finger /1.1m | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | OXY-F-UN  | U403S-127   | Adult soft /1.1m       | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U103S-127   | Pediatric soft /1.1m   | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U203S-127   | Infant soft /1.1m      | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U303-127    | Neonate wrap /1.1m     | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | OXY-W-UN  | U603-127    | Adult/Neo wrap/1.1m    | C8         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | OXY-W-UN  | U603B-127   | Pedia/infant wrap/1.1m | C8         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | OXY-W-UN  | U803-127    | Multi-site Y/1.1m      | C8         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | OXY-E-UN  | U903-127    | Adult ear clip /1.1m   | D1         |                                                 |                                                                     |                                                                                       |
| GE Trusignal                                                                                                                                                                                                                           |                                                                                     |           |             |                        |            |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |  | TS-F-D    | U403-117    | Adult finger /1.1m     | C8         | Required                                        | Depends on the monitor models<br>See the page of Spo2 adapter cable |  |
|                                                                                                                                                                                                                                        |                                                                                     |           | U103-117    | Pediatric finger /1.1m | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | TS-SA-D   | U403S-117   | Adult soft /1.1m       | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U103S-117   | Pediatric soft /1.1m   | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U203S-117   | Infant soft /1.1m      | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U303-117    | Neonate wrap /1.1m     | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | TS-W-D    | U603-117    | Adult/Neo wrap/1.1m    | C8         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U603B-117   | Pedia/infant wrap/1.1m | C8         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | TS-W-D    | U803-117    | Multi-site Y/1.1m      | C8         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     | TS-E-D    | U903-117    | Adult ear clip /1.1m   | D1         |                                                 |                                                                     |                                                                                       |
| General (China)                                                                                                                                                                                                                        |                                                                                     |           |             |                        |            |                                                 |                                                                     |                                                                                       |
| G1B<br>(*Digital Tech)                                                                                                                                                                                                                 |  |           | U403-42D    | Adult finger /1.1m     | C5         | No need                                         |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U403S-42D   | Adult soft /1.1m       | C7         |                                                 |                                                                     |                                                                                       |
| G2A,G3A,G3C,G3D,<br>G3F,G3G,G3H,G3L,<br>G6A,G6A plus, G9L<br>(*Digital Tech)                                                                                                                                                           |  |           | U410-65D    | Adult finger /3.0m     | C8         | No need                                         |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U410S-65D   | Adult soft /3.0m       | D1         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U810-65D    | Multi-site Y/3.0m      | C8         |                                                 |                                                                     |                                                                                       |
| Guoteng (China)                                                                                                                                                                                                                        |                                                                                     |           |             |                        |            |                                                 |                                                                     |                                                                                       |
| PM300                                                                                                                                                                                                                                  |  | DS100A    | U403-01     | Adult finger /1.1m     | C3         | No need                                         |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U403S-01    | Adult soft /1.1m       | C5         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |  |           | U410-171    | Adult finger /3.0m     | C8         | No need                                         |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U410S-171   | Adult soft /3.0m       | D1         |                                                 |                                                                     |                                                                                       |
| GT3000,GT5000,<br>GT6000GT7000,<br>GT9000,GT9002                                                                                                                                                                                       |  |           | U410-179    | Adult finger /3.0m     | C7         | No need                                         |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U410S-179   | Adult soft /3.0m       | C8         |                                                 |                                                                     |                                                                                       |
|                                                                                                                                                                                                                                        |                                                                                     |           | U810-179    | Multi-site Y/3.0m      | C7         |                                                 |                                                                     |                                                                                       |


| OEM Manufacturer & Models                                            | Sensor Connector                                                                      | REF# OEM               | REF# UNIMED            | Sensor Descriptions | Price Code    | Ext - Cable | Cable Plug                                                                          | Sensor Conn.                                                                        |
|----------------------------------------------------------------------|---------------------------------------------------------------------------------------|------------------------|------------------------|---------------------|---------------|-------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Goldway(China)                                                       |                                                                                       |                        |                        |                     |               |             |                                                                                     |                                                                                     |
| UT4000A,UT4000B, UT4000F<br>UT4000ACS                                |    |                        | U410-37                | Adult finger /3.0m  | C7            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-37               | Adult soft /3.0m    | C8            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-37                | Multi-site Y/3.0m   | C7            |             |                                                                                     |                                                                                     |
|                                                                      |    | See Nellcor Non-Oximax | See Nellcor Non-Oximax | 1.1m                | See Nellcor   | U708-37     |  |  |
|                                                                      |    | See Pace tech          | See Pace tech          | See Pace tech       | See Pace tech | No need     |                                                                                     |                                                                                     |
| UT4000A,UT4000C, UT4000EUT4000F, UT4000APro<br>UT4000FPro, UT6000A   |    |                        | U410-11                | Adult finger /3.0m  | C8            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-11               | Adult soft /3.0m    | D1            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-11                | Multi-site Y/3.0m   | C8            |             |                                                                                     |                                                                                     |
|                                                                      |    | See Nellcor Non-Oximax | See Nellcor Non-Oximax | 1.1m                | See Nellcor   | U708-11     |  |  |
|                                                                      |    | See Pace tech          | See Pace tech          | See Pace tech       | See Pace tech | No need     |                                                                                     |                                                                                     |
| 3100                                                                 |    | Digital sensor<br>数字探头 | U410-61D               | Adult finger /3.0m  | C8            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-61D              | Adult soft /3.0m    | D1            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-67D               | Multi-site Y/3.0m   | C8            |             |                                                                                     |                                                                                     |
| UT4000B,UT4000F                                                      |  |                        | U410-67D               | Adult finger /3.0m  | C8            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-67D              | Adult soft /3.0m    | D1            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-67D               | Multi-site Y/3.0m   | C8            |             |                                                                                     |                                                                                     |
| (*Oximax)                                                            |  |                        | U410-38P               | Adult finger /3.0m  | D1            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-38P              | Adult soft /3.0m    | D2            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-38P               | Multi-site Y/3.0m   | D1            |             |                                                                                     |                                                                                     |
| UT6000A<br>G3<br>(*Philips tech)                                     |  |                        | U410-91                | Adult finger /3.0m  | C8            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-91               | Adult soft /3.0m    | D1            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-91                | Multi-site Y/3.0m   | C8            |             |                                                                                     |                                                                                     |
| FJBF(China)                                                          |                                                                                       |                        |                        |                     |               |             |                                                                                     |                                                                                     |
| NLB7000N/8000B/8000N<br>NLB_V3,NBL_K3<br>乐夫:LF-6000A/6000B/<br>6000C |  |                        | U410-134D              | Adult finger /3.0m  | C8            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-134D             | Adult soft /3.0m    | D1            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-134D              | Multi-site Y/3.0m   | C8            |             |                                                                                     |                                                                                     |
| Infinium Medical                                                     |                                                                                       |                        |                        |                     |               |             |                                                                                     |                                                                                     |
| OmniOne, Omni II<br>Omni III<br>Infitron                             |  |                        | U410-97D               | Adult finger /3.0m  | C8            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-97D              | Adult soft /3.0m    | D1            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-97D               | Multi-site Y/3.0m   | C8            |             |                                                                                     |                                                                                     |
| DRE Waveline Pro<br>(*Digital Tech)                                  |  |                        | U410-137D              | Adult finger /3.0m  | C8            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-137D             | Adult soft /3.0m    | D1            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-137D              | Multi-site Y/3.0m   | C8            |             |                                                                                     |                                                                                     |
| DRE Waveline Pro+                                                    |  |                        | U410-147               | Adult finger /3.0m  | C7            | No need     |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U410S-147              | Adult soft /3.0m    | C8            |             |                                                                                     |                                                                                     |
|                                                                      |                                                                                       |                        | U810-147               | Multi-site Y/3.0m   | C7            |             |                                                                                     |                                                                                     |



| OEM Manufacturer & Models                                                            | Sensor Connector                                                                    | REF# OEM          | REF# UNIMED   | Sensor Descriptions | Price Code    | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|--------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------|---------------|---------------------|---------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Innocare                                                                             |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| VET-Multi-Parameter                                                                  |    |                   | U410-204      | Adult finger /3.0m  | C7            | No need     |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U410S-204     | Adult soft /3.0m    | C8            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U810-204      | Multi-site Y/3.0m   | C7            |             |                                                                                       |                                                                                       |
| IVY                                                                                  |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| 405A,405D<br>(* Non Oximax)                                                          |    | DS100A            | U403-01       | Adult finger /1.1m  | C3            | U710-85     |    |    |
|                                                                                      |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| 405P,405T<br>(*Non Oximax)                                                           |    | See the above     | See the above | See the above       | See the above | U708-01     |    |    |
| Vital-Guard 450CN(M)<br>(*M tech)                                                    |    | 1863/<br>LNCS DCI | U403-49R      | Adult finger /1.1m  | C4            | U708M-15R   |    |    |
|                                                                                      |                                                                                     | 1863/<br>LNCS DCI | U403S-49R     | Adult soft /1.1m    | C6            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     | 2258/<br>LNCS YI  | U803-49R      | Multi-site Y/1.1m   | C4            |             |                                                                                       |                                                                                       |
|                                                                                      |    | LNOP<br>DC-12     | U410-15R      | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     | LNOP<br>DC-12     | U410S-15R     | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U810-15R      | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| 450C & CM5<br>(*Nellcr Oximax)                                                       |    | LNOP<br>DC-12     | U410-70P      | Adult finger /3.0m  | D1            | No need     |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     | LNOP<br>DC-12     | U410S-70P     | Adult soft /3.0m    | D2            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U810-70P      | Multi-site Y/3.0m   | D1            |             |                                                                                       |                                                                                       |
| Heal Force                                                                           |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| Classic 90,Classic 120,<br>Prince 100F                                               |   |                   | U410-61D      | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U410S-61D     | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U810-61D      | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| Jeteem                                                                               |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| BD-6000<br>BD-9000                                                                   |  |                   | U410-174      | Adult finger /3.0m  | C7            | No need     |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U410S-174     | Adult soft /3.0m    | C8            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U810-174      | Multi-site Y/3.0m   | C7            |             |                                                                                       |                                                                                       |
| JNR                                                                                  |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| DS-100                                                                               |  |                   | U410-188      | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U410S-188     | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U810-188      | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| Kontron                                                                              |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| 7840,7845,Minimon<br>7337Micromon,<br>Supermon,Kolormon,<br>Trakmon<br>(*Non Oximax) |  |                   | U410-13       | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U410S-13      | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U810-13       | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
|                                                                                      |  | DS100A            | U403-01       | Adult finger /1.1m  | C3            | U708-13     |  |  |
|                                                                                      |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| 7150,7850,<br>Ivy 2000(M)                                                            |  | See the above     | See the above | See the above       | See the above | U708-01     |  |  |
|                                                                                      |  | See the above     | See the above | See the above       | See the above | U708-53     |  |  |
| Lohmeier                                                                             |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
|                                                                                      |  |                   | U410-64       | Adult finger /3.0m  | D1            | No need     |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U410S-64      | Adult soft /3.0m    | D2            |             |                                                                                       |                                                                                       |
|                                                                                      |                                                                                     |                   | U810-64       | Multi-site Y/3.0m   | D1            |             |                                                                                       |                                                                                       |


| OEM Manufacturer & Models         | Sensor Connector                                                                      | REF# OEM                                                    | REF# UNIMED                                                                           | Sensor Descriptions    | Price Code | Ext - Cable | Cable Plug                                                                          | Sensor Conn.                                                                          |
|-----------------------------------|---------------------------------------------------------------------------------------|-------------------------------------------------------------|---------------------------------------------------------------------------------------|------------------------|------------|-------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| IBIO R&DT                         |                                                                                       |                                                             |                                                                                       |                        |            |             |                                                                                     |                                                                                       |
| RD9012                            |    |                                                             | U410-100                                                                              | Adult finger /3.0m     | C7         | No need     |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U410S-100                                                                             | Adult soft /3.0m       | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U810-100                                                                              | Multi-site Y/3.0m      | C7         |             |                                                                                     |                                                                                       |
| (*non Oximax)                     |    | DS100A                                                      | U410-01                                                                               | Adult finger /3.0m     | C5         | U708-56     |  |    |
|                                   |                                                                                       | FlexMAX                                                     | U410S-01                                                                              | Adult soft /3.0m       | C7         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | D-YS                                                        | U810-01                                                                               | Multi-site Y/3.0m      | C5         |             |                                                                                     |                                                                                       |
| *M Technology                     |                                                                                       |                                                             |                                                                                       |                        |            |             |                                                                                     |                                                                                       |
| Rad-8<br>Rad-5/Rad-5v<br>SatShare |    | 1269/<br>LNOP DCI                                           | U403-125R                                                                             | Adult finger /1.1m     | C7         | Required    |                                                                                     |    |
|                                   |                                                                                       | 1276/<br>LNOP DCIP                                          | U103-125R                                                                             | Pediatric finger /1.1m | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | 1269/<br>LNOP DCI                                           | U403S-125R                                                                            | Adult soft成 /1.1m      | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | 1276/<br>LNOP DCIP                                          | U103S-125R                                                                            | Pediatric soft /1.1m   | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U203S-125R                                                                            | Infant soft /1.1m      | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U303-125R                                                                             | Neonate wrap /1.1m     | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U603-125R                                                                             | Adult/Neo wrap/1.1m    | C7         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U603B-125R                                                                            | Pedia/infant wrap/1.1m | C7         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | 1504/<br>LNOP YI                                            | U803-125R                                                                             | Multi-site Y/1.1m      | C7         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | LNOP TC-I                                                   | U903-125R                                                                             | Adult ear clip /1.1m   | C8         |             |                                                                                     |                                                                                       |
|                                   |  | 1863/<br>LNCS DCI                                           | U403-49R                                                                              | Adult finger /1.1m     | C4         | Required    |                                                                                     |  |
|                                   |                                                                                       | 1864/<br>LNCS DCIP                                          | U103-49R                                                                              | Pediatric finger /1.1m | C6         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | 1863/<br>LNCS DCI                                           | U403S-49R                                                                             | Adult soft /1.1m       | C6         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U103S-49R                                                                             | Pediatric soft /1.1m   | C6         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U203S-49R                                                                             | Infant soft /1.1m      | C6         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U303-49R                                                                              | Neonate wrap /1.1m     | C6         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U603-49R                                                                              | Adult/Neo wrap/1.1m    | C4         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U603B-49R                                                                             | Pedia/infant wrap/1.1m | C4         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | 2258/<br>LNCS YI                                            | U803-49R                                                                              | Multi-site Y/1.1m      | C4         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | LNCS TC-I                                                   | U903-49R                                                                              | Adult ear clip /1.1m   | C6         |             |                                                                                     |                                                                                       |
| M Tech<br>*Refurbished<br>Only    |  | LNOP DCSC                                                   | U402-15R                                                                              | Adult finger /0.9m     | C8         | No need     |                                                                                     |                                                                                       |
|                                   |                                                                                       | LNOP DC-12                                                  | U410-15R                                                                              | Adult finger /3.0m     | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U110-15R                                                                              | Pediatric finger /3.0m | D1         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | LNOP DC-12                                                  | U410S-15R                                                                             | Adult soft /3.0m       | D1         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U110S-15R                                                                             | Pediatric soft /3.0m   | D1         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U210S-15R                                                                             | Infant soft /3.0m      | D1         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U310-15R                                                                              | Neonate wrap /3.0m     | D1         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U610-15R                                                                              | Adult/Neo wrap/3.0m    | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U610B-15R                                                                             | Pedia/infant wrap/3.0m | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U810-15R                                                                              | Multi-site Y/3.0m      | C8         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       |                                                             | U910-15R                                                                              | Adult ear clip /3.0m   | D1         |             |                                                                                     |                                                                                       |
|                                   |                                                                                       | Radical-7,Rad-87<br>Rad-57,Pronto-7,<br>Rad-8Masimo Rainbow |  | DCI-dc3                | U410-87    |             |                                                                                     |                                                                                       |
|                                   | U410S                                                                                 |                                                             |                                                                                       | Adult soft/3.0m        | D3         |             |                                                                                     |                                                                                       |
|                                   | U810-87                                                                               |                                                             |                                                                                       | Multi-site Y/3.0m      | D2         |             |                                                                                     |                                                                                       |


| OEM Manufacturer & Models                                 | Sensor Connector                                                                    | REF# OEM          | REF# UNIMED   | Sensor Descriptions | Price Code    | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|-----------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------|---------------|---------------------|---------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| MDE/Invivo                                                |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| Escort 1<br>Escort 300<br>(*Non Oximax)                   |    | DS100A            | U403-01       | Adult finger /1.1m  | C3            | U710-85     |    |    |
|                                                           |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| 20300,20401,20403,<br>Prism<br>(*Non Oximax)              |    | See the above     | See the above | See the above       | See the above | U708-01     |    |    |
| Millenia, Escort,(M)<br>M10TM,M12TM,<br>M8TM<br>(*M Tech) |    | 1863/<br>LNCS DCI | U403-49R      | Adult finger /1.1m  | C4            | U708M-15R   |    |    |
|                                                           |                                                                                     | 1863/<br>LNCS DCI | U403S-49R     | Adult soft /1.1m    | C6            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     | 2258/<br>LNCS YI  | U803-49R      | Multi-site Y/1.1m   | C4            |             |                                                                                       |                                                                                       |
|                                                           |    | LNOP DC-12        | U410-15R      | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                           |                                                                                     | LNOP DC-12        | U410S-15R     | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U810-15R      | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| 4500, ,Omni Track<br>Omeg 1445Millennia                   |    |                   | U410-14       | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U410S-14      | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U810-14       | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| 4500, ,Omni Track<br>Omeg 1445Millennia<br>M Tech         |    |                   | U410-14M      | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U410S-14M     | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U810-14M      | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| Meditronic - Physio Control                               |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| Lifepack 12<br>(*non Oximax)                              |   | DS100A            | U403-01       | Adult finger /1.1m  | C3            | U708-01     |   |   |
|                                                           |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| Lifepack 12, 20(M)<br>(*M Tech)                           |  | 1863/<br>LNCS DCI | U403-49R      | Adult finger /1.1m  | C4            | U708M-15R   |  |  |
|                                                           |                                                                                     | 1863/<br>LNCS DCI | U403S-49R     | Adult soft /1.1m    | C6            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     | 2258/<br>LNCS YI  | U803-49R      | Multi-site Y/1.1m   | C4            |             |                                                                                       |                                                                                       |
|                                                           |  | LNOP DC-12        | U410-15R      | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                           |                                                                                     | LNOP DC-12        | U410S-15R     | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U810-15R      | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |
| Medlab                                                    |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| NANOX 10 POX 10<br>Perarl 10<br>(*BCI)                    |  | 3044              | U403-06       | Adult finger /1.1m  | C3            | U708-76     |  |  |
|                                                           |                                                                                     |                   | U403S-06      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     | 3043              | U803-06       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| MEK                                                       |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
| MP 100,110,<br>400,1000                                   |  |                   | U410-26       | Adult finger /3.0m  | C6            | No need     |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U410S-26      | Adult soft /3.0m    | C7            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U810-26       | Multi-site Y/3.0m   | C6            |             |                                                                                       |                                                                                       |
| (*non Oximax)                                             |  | DS100A            | U403-01       | Adult finger /1.1m  | C3            | U708-56     |  |  |
|                                                           |                                                                                     | DS100A            | U403S-01      | Adult soft /1.1m    | C5            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     | D-YS              | U803-01       | Multi-site Y/1.1m   | C3            |             |                                                                                       |                                                                                       |
| MEK                                                       |                                                                                     |                   |               |                     |               |             |                                                                                       |                                                                                       |
|                                                           |  |                   | U410-56       | Adult finger /3.0m  | C8            | No need     |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U410S-56      | Adult soft /3.0m    | D1            |             |                                                                                       |                                                                                       |
|                                                           |                                                                                     |                   | U810-56       | Multi-site Y/3.0m   | C8            |             |                                                                                       |                                                                                       |

| OEM Manufacturer & Models                                                                                           | Sensor Connector                                                                      | REF# OEM                                                                              | REF# UNIMED   | Sensor Descriptions    | Price Code         | Ext - Cable                | Cable Plug                                                                            | Sensor Conn.                                                                          |  |
|---------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------|------------------------|--------------------|----------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--|
| Mennen                                                                                                              |                                                                                       |                                                                                       |               |                        |                    |                            |                                                                                       |                                                                                       |  |
| Envoy<br>(*non Oximax)                                                                                              |    | DS100A                                                                                | U403-01       | Adult finger /1.1m     | C3                 | U708-54                    |    |    |  |
|                                                                                                                     |                                                                                       | DS100A                                                                                | U403S-01      | Adult soft /1.1m       | C5                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       | D-YS                                                                                  | U803-01       | Multi-site Y/1.1m      | C3                 |                            |                                                                                       |                                                                                       |  |
| Horizon 1100,2000XL<br>Mercury,MR1300<br>(*non Oximax)                                                              |    | See the above                                                                         | See the above | See the above          | See the above      | U710-85                    |    |    |  |
| Horizon XL,Mercury<br>1200,MR1300<br>(*non Oximax)                                                                  |    | See the above                                                                         | See the above | See the above          | See the above      | U708-24                    |    |    |  |
| Horizon XL<br>(*Ohmeda)                                                                                             |    | OXY-F4-H                                                                              | U410-02       | Adult finger /3.0m     | C8                 | No need                    |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U410S-02      | Adult soft /3.0m       | D1                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       | OXY-W4-H                                                                              | U810-02       | Multi-site Y/3.0m      | C8                 |                            |                                                                                       |                                                                                       |  |
| Horizon,Envoy,<br>Vitalogik (M)<br>(*M Tech)                                                                        |    | 1863/<br>LNCS DCI                                                                     | U403-49       | Adult finger /1.1m     | C4                 | U708M-59                   |    |    |  |
|                                                                                                                     |                                                                                       |                                                                                       | U403S-49      | Adult soft /1.1m       | C6                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       | 2258/<br>LNCS YI                                                                      | U803-49       | Multi-site Y/1.1m      | C6                 |                            |                                                                                       |                                                                                       |  |
| Mindray (China)                                                                                                     |                                                                                       |                                                                                       |               |                        |                    |                            |                                                                                       |                                                                                       |  |
| 509B, MEC-1000,<br>MEC-1200<br>New PM6000<br>PM7000/8000/9000,<br>VS 800,PM600<br>PM9000 Express<br>(*Mindray Tech) |  |                                                                                       | U410-29       | Adult finger /3.0m     | C7                 | No need                    |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U110-29       | Pediatric finger /3.0m | C8                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U410S-29      | Adult soft /3.0m       | C8                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U110S-29      | Pediatric soft /3.0m   | C8                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U210S-29      | Infant soft /3.0m      | C8                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U310-29       | Neonate wrap /3.0m     | C8                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U610-29       | Adult/Neo wrap/3.0m    | C7                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U610B-29      | Pedia/infant wrap/3.0m | C7                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U810-29       | Multi-site Y/3.0m      | C7                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U910-29       | Adult ear clip/3.0m    | C8                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |  | 512D                                                                                  | U403-01       | Adult finger /1.1m     | C3                 | U708-29<br>(0010-20-42594) |  |  |  |
|                                                                                                                     |                                                                                       | DS100A                                                                                | U403S-01      | Adult soft /1.1m       | C5                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       | D-YS                                                                                  | U803-01       | Multi-site Y/1.1m      | C3                 |                            |                                                                                       |                                                                                       |  |
| PM-5000,PM 6000<br>PM-8000<br>(*Mindray Tech)                                                                       |  | See the above                                                                         | See the above | See the above          | See the above      | U708-22<br>(0515-30-11221) |  |  |  |
|                                                                                                                     |                                                                                       |  |               | U410-22                | Adult finger /3.0m | C8                         | No need                                                                               |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       |               | U410S-22               | Adult soft /3.0m   | D1                         |                                                                                       |                                                                                       |  |
|                                                                                                                     |  |                                                                                       | U810-22       | Multi-site Y/3.0m      | C8                 | No need                    |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U410-48MR     | Adult finger /3.0m     | D1                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       | U410S-48MR    | Adult soft /3.0m       | D2                 |                            |                                                                                       |                                                                                       |  |
| Beneview T5/T8<br>(Masimo Rainbow)                                                                                  |  |                                                                                       | U810-48MR     | Multi-site Y/3.0m      | D1                 |                            |                                                                                       |                                                                                       |  |
|                                                                                                                     |                                                                                       |                                                                                       |               |                        |                    |                            |                                                                                       |                                                                                       |  |
| Mindray (China)                                                                                                     |                                                                                       |                                                                                       |               |                        |                    |                            |                                                                                       |                                                                                       |  |

| Mindray (China)                                                                         |                                                                                     |                                       |                                       |                        |                    |                             |                                                                                       |                                                                                       |  |
|-----------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------|---------------------------------------|------------------------|--------------------|-----------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--|
| OEM Manufacturer & Models                                                               | Sensor Connector                                                                    | REF# OEM                              | REF# UNIMED                           | Sensor Descriptions    | Price Code         | Ext - Cable                 | Cable Plug                                                                            | Sensor Conn.                                                                          |  |
| Mindray (China)                                                                         |                                                                                     |                                       |                                       |                        |                    |                             |                                                                                       |                                                                                       |  |
| PM-6000<br>(*M tech)                                                                    |    | 1863/ LNCS DCI                        | *U403-49                              | Adult finger /1.1m     | C4                 | U708M-22<br>(0010-30-12452) |    |    |  |
|                                                                                         |                                                                                     | 1863/ LNCS DCI                        | *U403S-49                             | Adult soft /1.1m       | C6                 |                             |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     | 2258/ LNCS YI                         | *U803-49                              | Multi-site Y/1.1m      | C4                 |                             |                                                                                       |                                                                                       |  |
|                                                                                         |    |                                       | *U410-22M                             | Adult finger /3.0m     | C8                 | No need                     |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | *U410S-22M                            | Adult soft /3.0m       | D1                 |                             |                                                                                       |                                                                                       |  |
| New PM-6000,<br>PM-8000<br>(*M tech)                                                    |    |                                       | *U810-22M                             | Multi-site Y/3.0m      | C8                 | No need                     |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | *U410-29M                             | Adult finger /3.0m     | C7                 |                             |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | *U410S-29M                            | Adult soft /3.0m       | C8                 |                             |                                                                                       |                                                                                       |  |
| New PM-6000<br>PM7000,8000,9000<br>PM-9000 Express<br>VS 800<br>New PM6000<br>(*M tech) |    | See the above<br>M sensor             | See the above<br>M sensor             | See the above M sensor | See<br>the above   | U708M-30                    |    |    |  |
|                                                                                         |                                                                                     |                                       |                                       |                        |                    |                             |                                                                                       |                                                                                       |  |
|                                                                                         |    |                                       |                                       |                        |                    | No need                     |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       |                                       |                        |                    |                             |                                                                                       |                                                                                       |  |
| New PM-6000<br>PM7000,8000,9000<br>VS 800<br>PM-9000 Express<br>(*Oximax)               |   | DS100A                                | U403-01XP                             | Adult finger /1.1m     | C5                 | U710X-180                   |   |   |  |
|                                                                                         |                                                                                     | DS100A                                | U403S-01XP                            | Adult soft /1.1m       | C7                 |                             |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     | D-YS                                  | U803-01XP                             | Multi-site Y/1.1m      | C5                 |                             |                                                                                       |                                                                                       |  |
|                                                                                         |  |                                       |                                       | U410-180P              | Adult finger /3.0m | C8                          | No need                                                                               |                                                                                       |  |
|                                                                                         |                                                                                     |                                       |                                       | U410S-180P             | Adult soft /3.0m   | D1                          |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | U810-180P                             | Multi-site Y/3.0m      | C8                 |                             |                                                                                       |                                                                                       |  |
| Beneview T5/T8<br>(*Mindray Tech)                                                       |  | See the above<br>Non Oximax<br>Sensor | See the above<br>Non Oximax<br>Sensor | 1.1m                   | See<br>the above   | U708-48R                    |  |  |  |
|                                                                                         |                                                                                     |                                       |                                       |                        |                    |                             |                                                                                       |                                                                                       |  |
|                                                                                         |  |                                       |                                       | U410-48                | Adult finger /3.0m | C8                          | No need                                                                               |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | U410S-48                              | Adult soft /3.0m       | D1                 |                             |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | U810-48                               | Multi-site Y/3.0m      | C8                 |                             |                                                                                       |                                                                                       |  |
| Beneview T5/T8<br>(Oximax*)                                                             |  | See Nellcor<br>Oximax                 | See Nellcor<br>Oximax                 | 1.1m                   | See<br>Oximax      | U710X-48                    |  |  |  |
| Beneview T5/T8<br>(*M tech)                                                             |  | See M                                 | See M                                 | 1.1m                   | See M              | U708M-48                    |  |  |  |
| 503,507,508,<br>509,PM9000<br>(*Criticare tech)                                         |  |                                       | U410-05                               | Adult finger /3.0m     | C8                 | No need                     |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | U410S-05                              | Adult soft /3.0m       | D1                 |                             |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | U810-05                               | Multi-site Y/3.0m      | C8                 |                             |                                                                                       |                                                                                       |  |
| 500A,509B,<br>PM9000<br>(*BCI)                                                          |  |                                       | U410-31                               | Adult finger /3.0m     | C8                 | No need                     |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | U410S-31                              | Adult soft /3.0m       | D1                 |                             |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     |                                       | U810-31                               | Multi-site Y/3.0m      | C8                 |                             |                                                                                       |                                                                                       |  |
| PM-50,PM-60<br>(*Non-Oximax)                                                            |  | DS100A                                | U403-01                               | Adult finger /1.1m     | C3                 | No need                     |                                                                                       |                                                                                       |  |
|                                                                                         |                                                                                     | DS100A                                | U403S-01                              | Adult soft /1.1m       | C5                 |                             |                                                                                       |                                                                                       |  |


| OEM Manufacturer & Models                                                                                          | Sensor Connector                                                                      | REF# OEM           | REF# UNIMED        | Sensor Descriptions    | Price Code    | Ext - Cable        | Cable Plug                                                                            | Sensor Conn.                                                                          |
|--------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--------------------|--------------------|------------------------|---------------|--------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Minolta                                                                                                            |                                                                                       |                    |                    |                        |               |                    |                                                                                       |                                                                                       |
| PULSOX-300<br>PULSOX-300i                                                                                          |    | SR-5C              | U401-105           | Adult finger /0.3m     | C8            | No need            |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | SR-5C              | U401S-105          | Adult soft /0.3m       | D1            |                    |                                                                                       |                                                                                       |
| M&B(China)                                                                                                         |                                                                                       |                    |                    |                        |               |                    |                                                                                       |                                                                                       |
| MB526,MB526+,MB526T8<br>MB526T12,MB526+tw<br>MB800,CD2000,CD2000+                                                  |    |                    | U410-92            | Adult finger /3.0m     | C7            | No need            |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U410S-92           | Adult soft /3.0m       | C8            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U810-92            | Multi-site Y/3.0m      | C7            |                    |                                                                                       |                                                                                       |
| CD2000,CD2000+<br>(* Oximax tech)                                                                                  |    |                    | U410-92P           | Adult finger /3.0m     | C8            | No need            |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U410S-92P          | Adult soft /3.0m       | D1            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U810-92P           | Multi-site Y/3.0m      | C8            |                    |                                                                                       |                                                                                       |
| Nellcor-Tyco                                                                                                       |                                                                                       |                    |                    |                        |               |                    |                                                                                       |                                                                                       |
| NPB 40, NPB20,<br>NPB75<br>NPB 180,<br>NPB 185,<br>NPB-190,<br>NPB-195,<br>NPB 3900,<br>NPB 4000<br>(*Non-Oximax)+ |    | DS100A             | U403-01            | Adult finger /1.1m     | C3            | U704-01<br>U708-01 |    |    |
|                                                                                                                    |                                                                                       | OXI-P/I            | U103-01            | Pediatric finger /1.1m | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | DS100A             | U403S-01           | Adult soft /1.1m       | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | OXI-P/I            | U103S-01           | Pediatric soft /1.1m   | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | OXI-P/I            | U203S-01           | Infant soft /1.1m      | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | OXI-A/N            | U303-01            | Neonate wrap /1.1m     | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | D-YS               | U803-01            | Multi-site Y/1.1m      | C3            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | OXI-A/N            | U603-01            | Adult/Neowrap/1.1m     | C3            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | OXI-P/I            | U603B-01           | Pedia/infant wrap/1.1m | C3            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       | D-YS+YSE           | U903-01            | Adult ear clip /1.1m   | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |  |                    | U410-01            | Adult finger /3.0m     | C5            | No need            |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U110-01            | Pediatric finger /3.0m | C7            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U410S-01           | Adult soft /3.0m       | C7            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U110S-01           | Pediatric soft /3.0m   | C7            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U210S-01           | Infant soft /3.0m      | C7            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U310-01            | Neonate wrap /3.0m     | C7            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U610-01            | Adult/Neowrap/3.0m     | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U610B-01           | Pedia/infant wrap/3.0m | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U810-01            | Multi-site Y/3.0m      | C5            |                    |                                                                                       |                                                                                       |
|                                                                                                                    |                                                                                       |                    | U910-01            | Adult ear clip /3.0m   | C7            |                    |                                                                                       |                                                                                       |
| N-100,N-200                                                                                                        |  | See the above      | See the above      | 1.1m                   | See the above | U710-85            |  |  |
| N550,N560,N595,<br>N600,N751, N5500                                                                                |  | See Nellcor Oximax | See Nellcor Oximax | 1.1m                   | See Oximax    | U710-70            |  |  |
| NPB3000-3100,<br>Symphony,<br>NPB-290,<br>NPB-295, NPB-395                                                         |  | See the above      | See the above      | 1.1m                   | See the above | U708-69            |  |  |
| Nellcor-Tyco                                                                                                       |                                                                                       |                    |                    |                        |               |                    |                                                                                       |                                                                                       |



| ↳ Nellcor-Tyco                                                                                                              |                                                                                     |                                                                                     |                           |                           |            |             |            |                                                                                       |
|-----------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------|---------------------------|------------|-------------|------------|---------------------------------------------------------------------------------------|
| OEM Manufacturer & Models                                                                                                   | Sensor Connector                                                                    | REF# OEM                                                                            | REF# UNIMED               | Sensor Descriptions       | Price Code | Ext - Cable | Cable Plug | Sensor Conn.                                                                          |
| Nellcor-Tyco                                                                                                                |                                                                                     |                                                                                     |                           |                           |            |             |            |                                                                                       |
|                                                                                                                             |    |                                                                                     | U410-69                   | Adult finger /3.0m        | C8         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U110-69                   | Pediatric finger /3.0m    | D1         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U410S-69                  | Adult soft /3.0m          | D1         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U110S-69                  | Pediatric soft /3.0m      | D1         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U210S-69                  | Infant soft /3.0m         | D1         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U310-69                   | Neonate wrap /3.0m        | D1         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U610-69                   | Adult/Neowrap/3.0m        | C8         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U610B-69                  | Pedia/infant wrap/3.0m    | C8         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U810-69                   | Multi-site Y/3.0m         | C8         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U910-69                   | Adult ear clip /3.0m      | D1         |             |            |                                                                                       |
| N-65<br>(*Oximax)<br>Ending “X”<br>can be replaced<br>by “XP” or “P”<br>for fully<br>compatible with<br>all Oximax oximetry |    | DS100A                                                                              | U403-01XP                 | Adult finger /1.1m        | C5         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     | OXI-P/I                                                                             | U103-01XP                 | Pediatric finger /1.1m    | C7         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     | OXI-A/N                                                                             | U403S-01XP                | Adult soft /1.1m          | C7         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     | OXI-P/I                                                                             | U103S-01XP                | Pediatric soft /1.1m      | C7         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     | OXI-P/I                                                                             | U203S-01XP                | Infant soft /1.1m         | C7         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     | OXI-A/N                                                                             | U303-01XP                 | Neonate wrap /1.1m        | C7         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     | OXI-A/N                                                                             | U603-01XP                 | Adult/Neowrap/1.1m        | C5         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     | OXI-P/I                                                                             | U603B-01XP                | Pedia/infant wrap/1.1m    | C5         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     | D-YS                                                                                | U803-01XP                 | Multi-site Y/1.1m         | C5         |             |            |                                                                                       |
|                                                                                                                             | D-YS+YSE                                                                            | U903-01XP                                                                           | Adult ear clip /1.1m      | C7                        |            |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |  | See Nellcor<br>Non-Oximax | See Nellcor<br>Non-Oximax | 1.1m       | See Nellcor | U701-01X   |  |
| Neusoft (China)                                                                                                             |                                                                                     |                                                                                     |                           |                           |            |             |            |                                                                                       |
| NSC series<br>(*Oximax Tech)                                                                                                |  |                                                                                     | U410-78P                  | Adult finger /3.0m        | C8         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U410S-78P                 | Adult soft /3.0m          | D1         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U810-78P                  | Multi-site Y/3.0m         | C8         |             |            |                                                                                       |
|                                                                                                                             |  |                                                                                     | U410-107                  | Adult finger /3.0m        | C7         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U410S-107                 | Adult soft /3.0m          | C8         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U810-107                  | Multi-site Y/3.0m         | C7         |             |            |                                                                                       |
| New tech (China)                                                                                                            |                                                                                     |                                                                                     |                           |                           |            |             |            |                                                                                       |
| Pulse Oximer<br>NT1A,NT1B,NT1D<br>(Digital Sensor)                                                                          |  |                                                                                     | U403-42D                  | Adult finger /1.1m        | C5         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U403S-42D                 | Adult soft /1.1m          | C7         |             |            |                                                                                       |
| Patient monitor<br>Digital sensor                                                                                           |  |                                                                                     | U410-65D                  | Adult finger /3.0m        | C8         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U410S-65D                 | Adult soft /3.0m          | D1         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U810-65D                  | Multi-site Y/3.0m         | C8         |             |            |                                                                                       |
| NT2A,NT2C                                                                                                                   |  |                                                                                     | U410-165D                 | Adult finger /3.0m        | D1         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U410S-165D                | Adult soft /3.0m          | D2         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U810-165D                 | Multi-site Y/3.0m         | D1         |             |            |                                                                                       |
| NT3A,NT3B,<br>NT3E,NT3F,<br>NT3B60                                                                                          |  |                                                                                     | U410-167                  | Adult finger /3.0m        | C8         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U410S-167                 | Adult soft 成/3.0m         | D1         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U810-167                  | Multi-site Y/3.0m         | C8         |             |            |                                                                                       |
| NT3A,NT3B,<br>NT3E,NT3F,NT3B60<br>(*Digital Tech)                                                                           |  |                                                                                     | U410-173D                 | Adult finger /3.0m        | D1         | No need     |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U410S-173D                | Adult soft /3.0m          | D2         |             |            |                                                                                       |
|                                                                                                                             |                                                                                     |                                                                                     | U810-173D                 | Multi-site Y/3.0m         | D1         |             |            |                                                                                       |






| OEM Manufacturer & Models                                                                                                                                                      | Sensor Connector                                                                      | REF# OEM             | REF# UNIMED          | Sensor Descriptions    | Price Code           | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------------------|----------------------|------------------------|----------------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Nihon Kohden                                                                                                                                                                   |                                                                                       |                      |                      |                        |                      |             |                                                                                       |                                                                                       |
| BSM-1101/2101<br>/23101<br>7100,7200,<br>OLV-1100<br>ALife scope                                                                                                               |    |                      | U410-16              | Adult finger /3.0m     | C8                   | No need     |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U110-16              | Pediatric finger /3.0m | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U410S-16             | Adult soft /3.0m       | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U110S-16             | Pediatric soft /3.0m   | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U210S-16             | Infant soft /3.0m      | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U310-16              | Neonate wrap /3.0m     | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U610-16              | Adult/Neowrap/3.0m     | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U610B-16             | Pedia/infant wrap/3.0m | C8                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U810-16              | Multi-site Y/3.0m      | C8                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U910-16              | Adult ear clip /3.0m   | D1                   |             |                                                                                       |                                                                                       |
| Lifescope S,<br>BSS-9800<br>Lifescope M,<br>BSM 9510<br>Lifescope P<br>BSM 4100<br>Life scope A,<br>BSM 5100<br>Lifescope L,<br>BSM 2350<br>Lifescope I<br>BSM 2301K,<br>2303K |    | TL-201T              | U405-16              | Adult finger /1.5m     | C7                   | U710-36     |    |    |
|                                                                                                                                                                                |                                                                                       |                      | U105-16              | Pediatric finger /1.5m | C8                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U405S-16             | Adult soft /1.5m       | C8                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U105S-16             | Pediatric soft /1.5m   | C8                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U205S-16             | Infant soft /1.5m      | C8                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U305-16              | Neonate wrap /1.5m     | C8                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U605-16              | Adult/Neowrap/1.5m     | C7                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U605B-15             | Pedia/infant wrap/1.5m | C7                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       | TL-220T              | U805-16              | Multi-site Y/1.5m      | C7                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U905-16              | Adult ear clip /1.5m   | C8                   |             |                                                                                       |                                                                                       |
| Lifescope S,<br>BSS-9800<br>Lifescope M,<br>BSM 9510<br>Lifescope P<br>BSM 4100<br>Life scope A,<br>BSM 5100<br>Lifescope L,<br>BSM 2350<br>Lifescope I<br>BSM 2301K,<br>2303K |  |                      | U410-36              | Adult finger /3.0m     | D1                   | No need     |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U110-36              | Pediatric finger /3.0m | D2                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U410S-36             | Adult soft /3.0m       | D2                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U110S-36             | Pediatric soft /3.0m   | D2                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U210S-36             | Infant soft /3.0m      | D2                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U310-36              | Neonate wrap /3.0m     | D2                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U610-36              | Adult/Neowrap/3.0m     | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U610B-36             | Pedia/infant wrap/3.0m | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U810-36              | Multi-site Y/3.0m      | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U910-36              | Adult ear clip /3.0m   | D2                   |             |                                                                                       |                                                                                       |
| OPV-1500                                                                                                                                                                       |  | The same<br>as above | The same<br>as above | 1.1m                   | The same<br>as above | U710-106    |  |  |
| ☆ SVM-7501                                                                                                                                                                     |  |                      | U410-224             | Adult finger /3.0m     | C8                   | No need     |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       | TL-651S              | U410S-224            | Adult soft /3.0m       | D1                   |             |                                                                                       |                                                                                       |
|                                                                                                                                                                                |                                                                                       |                      | U810-224             | Multi-site Y/3.0m      | C8                   |             |                                                                                       |                                                                                       |

| OEM Manufacturer & Models                                                                               | Sensor Connector                                                                                                                                                        | REF# OEM          | REF# UNIMED       | Sensor Descriptions     | Price Code           | Ext - Cable | Cable Plug                                                                          | Sensor Conn.                                                                        |
|---------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|-------------------|-------------------------|----------------------|-------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Nonin                                                                                                   |                                                                                                                                                                         |                   |                   |                         |                      |             |                                                                                     |                                                                                     |
| 8500<br>2500<br>9840<br>Avant 4000,9600, 9700<br>Avant 2120<br>PulseSense<br>Model 7500,7500FO (Nonin*) |                                                                                        | 8000AA            | U403-08           | Adult finger /1.1m      | C3                   | No need     |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         | 8000AP            | U103-08           | Pediatric finger /1.1m  | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         | 8000SL            | U403S-08          | Adult soft /1.1m        | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         | 8000SM            | U103S-08          | Pediatric soft /1.1m    | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U203S-08          | Infant soft /1.1m       | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U303-08           | Neonate wrap /1.1m      | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U603-08           | Adult/Neowrap /1.1m     | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U603B-08          | Pedia/infant wrap /1.1m | C3                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U803-08           | Multi-site Y /1.1m      | C3                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | 8000Q2            | U903-08                 | Adult ear clip /1.1m |             |                                                                                     |                                                                                     |
| Avant 4000 Wearable<br>Digital Pulse Oximetry                                                           |                                                                                        |                   | U410-08L          | Adult finger /3.0m      | C5                   | No need     |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U401-146          | Adult finger /0.3m      | D1                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U101-146          | Pediatric finger /0.3m  | D1                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U401S-146         | Adult soft /0.3m        | D1                   |             |                                                                                     |                                                                                     |
| 8600                                                                                                    |                                                                                        | The same as above | The same as above | 1.1m                    | The same as above    | NA          |                                                                                     |  |
|                                                                                                         |                                                                                                                                                                         |                   |                   |                         |                      |             |                                                                                     |                                                                                     |
| 8604                                                                                                    | <br> | The same as above | The same as above | 1.1m                    | The same as above    | U708-08     |  |  |
|                                                                                                         |                                                                                                                                                                         |                   | U410-08           | Adult finger /3.0m      | C8                   | No need     |                                                                                     |                                                                                     |
| Novametrix                                                                                              |                                                                                                                                                                         |                   |                   |                         |                      |             |                                                                                     |                                                                                     |
| 512/513/710/715, MARS                                                                                   |                                                                                      |                   | U403-73           | Adult finger /1.1m      | C5                   | No need     |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U103-73           | Pediatric finger /1.1m  | C7                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U403S-73          | Adult soft /1.1m        | C7                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U103S-73          | Pediatric soft /1.1m    | C7                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U203S-73          | Infant soft /1.1m       | C7                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U303-73           | Neonate wrap /1.1m      | C7                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U603-73           | Adult/Neowrap /1.1m     | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U603B-73          | Pedia/infant wrap /1.1m | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U803-73           | Multi-site Y /1.1m      | C5                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U903-73           | Adult ear clip /1.1m    | C7                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                      |                   | U403-103          | Adult finger /1.1m      | C8                   | NA          |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U103-103          | Pediatric finger /1.1m  | D1                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U403S-103         | Adult soft /1.1m        | D1                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U103S-103         | Pediatric soft /1.1m    | D1                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U203S-103         | Infant soft /1.1m       | D1                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U303-103          | Neonate wrap /1.1m      | D1                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U603-103          | Adult/Neowrap /1.1m     | C8                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U603B-103         | Pedia/infant wrap /1.1m | C8                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U803-103          | Multi-site Y /1.1m      | C8                   |             |                                                                                     |                                                                                     |
|                                                                                                         |                                                                                                                                                                         |                   | U903-103          | Adult ear clip /1.1m    | D1                   |             |                                                                                     |                                                                                     |
| Novametrix                                                                                              |                                                                                                                                                                         |                   |                   |                         |                      |             |                                                                                     |                                                                                     |


Novametrix


| OEM Manufacturer & Models                                    | Sensor Connector                                                                      | REF# OEM               | REF# UNIMED            | Sensor Descriptions     | Price Code  | Ext - Cable | Cable Plug                                                                           | Sensor Conn.                                                                         |
|--------------------------------------------------------------|---------------------------------------------------------------------------------------|------------------------|------------------------|-------------------------|-------------|-------------|--------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| Novametrix                                                   |                                                                                       |                        |                        |                         |             |             |                                                                                      |                                                                                      |
| 505/510/511 /515/B/C/520<br>7000/7100,Oxypleth (Novametrix*) |    |                        | U410-03                | Adult finger /3.0m      | C8          | No need     |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U110-03                | Pediatric finger /3.0m  | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U410S-03               | Adult soft /3.0m        | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U110S-03               | Pediatric soft /3.0m    | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U210S-03               | Infant soft /3.0m       | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U310-03                | Neonate wrap /3.0m      | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U610-03                | Adult/Neowrap /3.0m     | C8          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U610B-03               | Pedia/infant wrap /3.0m | C8          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U810-03                | Multi-site Y/3.0m       | C8          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U910-03                | Adult ear clip /3.0m    | D1          |             |                                                                                      |                                                                                      |
| PaceTech/Generra                                             |                                                                                       |                        |                        |                         |             |             |                                                                                      |                                                                                      |
|                                                              |    | See Nellcor Non-Oximax | See Nellcor Non-Oximax | 1.1m                    | See Nellcor | U708-38     |   |   |
|                                                              |    |                        | U410-38                | Adult finger /3.0m      | C7          | No need     |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U410S-38               | Adult soft /3.0m        | C8          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U810-38                | Multi-site Y/3.0m       | C7          |             |                                                                                      |                                                                                      |
|                                                              |   | See Nellcor Non-Oximax | See Nellcor Non-Oximax | 1.1m                    | See Nellcor | U708-39     |  |  |
|                                                              |                                                                                       |                        | U410-39                | Adult finger /3.0m      | C8          | No need     |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U410S-39               | Adult soft /3.0m        | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U810-39                | Multi-site Y/3.0m       | C8          |             |                                                                                      |                                                                                      |
| Mediaid/Palco                                                |                                                                                       |                        |                        |                         |             |             |                                                                                      |                                                                                      |
| M300/M305/340, M5305/M5340                                   |  |                        | U410-18                | Adult finger /3.0m      | C8          | No need     |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U110-18                | Pediatric finger /3.0m  | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U410S-18               | Adult soft /3.0m        | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U110S-18               | Pediatric soft /3.0m    | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U210S-18               | Infant soft /3.0m       | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U310-18                | Neonate wrap /3.0m      | D1          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U610-18                | Adult/Neowrap /3.0m     | C8          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U610B-18               | Pedia/infant wrap /3.0m | C8          |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U810-18                | Multi-site Y/3.0m       | D1          |             |                                                                                      |                                                                                      |
| M30/40,M100,M130, M900,M960P, M6 Aster,M4 Nova, M1000        |  |                        | U410-148               | Adult finger /3.0m      | D1.3        | No need     |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U110-148               | Pediatric finger /3.0m  | D1.6        |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U410S-148              | Adult soft /3.0m        | D1.6        |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U110S-148              | Pediatric soft /3.0m    | D1.6        |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U210S-148              | Infant soft /3.0m       | D1.6        |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U310-148               | Neonate wrap /3.0m      | D1.6        |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U610-148               | Adult/Neowrap /3.0m     | D1.3        |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U610B-148              | Pedia/infant wrap /3.0m | D1.3        |             |                                                                                      |                                                                                      |
|                                                              |                                                                                       |                        | U810-148               | Multi-site Y/3.0m       | D1.6        |             |                                                                                      |                                                                                      |

| OEM Manufacturer & Models                                                                   | Sensor Connector                                                                    | REF# OEM               | REF# UNIMED            | Sensor Descriptions    | Price Code          | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|---------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|------------------------|------------------------|------------------------|---------------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Philips/HP                                                                                  |                                                                                     |                        |                        |                        |                     |             |                                                                                       |                                                                                       |
| 78352C,78354C, 78834C,M1020A, M1025B,M1025A, M1350B,M1772A/B, M1732A/B,M2475B (*Non-Oximax) |    | M1190A                 | <b>U410-20</b>         | Adult finger /3.0m     | C8                  | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U110-20</b>         | Pediatric finger /3.0m | D1                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     | M1190A                 | <b>U410S-20</b>        | Adult soft /3.0m       | D1                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U110S-20</b>        | Pediatric soft/3.0m    | D1                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U210S-20</b>        | Infant soft /3.0m      | D1                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U310-20</b>         | Neonate wrap /3.0m     | D1                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U610-20</b>         | Adult/Neowrap/3.0m     | C8                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U610B-20</b>        | Pedia/infant wrap3.0m  | C8                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U810-20</b>         | Multi-site Y/3.0m      | C8                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U910-20</b>         | Adult ear clip /3.0m   | D1                  |             |                                                                                       |                                                                                       |
| C1,C3,M3,M4, VirDia, MP20,MP30, MP40/50/60/70/90. VM4,VM6,VM8 (Philips*)                    |    | See Philips' sensor    | See Philips' sensor    | See Philips' sensor    | See Philips' sensor | U708-41     |    |    |
|                                                                                             |   | See Nellcor Non-Oximax | See Nellcor Non-Oximax | 1.1m                   | See Nellcor         | U708-43     |   |   |
| A1, A3 (*Non-Oximax)                                                                        |  | See Nellcor Non-Oximax | See Nellcor Non-Oximax | 1.1m/3.0m              | See Nellcor         | U708-01     |  |  |
| MP20/30/40/50 PM60/70/90 VN4,VM6,VM8 (*Oximax)                                              |  | DS100A                 | <b>U403-01P</b>        | Adult finger /1.1m     | C5                  | U710X-43    |  |  |
|                                                                                             |                                                                                     | DS100A                 | <b>U403S-01P</b>       | Adult soft /1.1m       | C7                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     | D-YS                   | <b>U803-01P</b>        | Multi-site Y/1.1m      | C5                  |             |                                                                                       |                                                                                       |
| MP20/30/40/ PM50/60/70/90 VN4,VM6,VM8 (*M tech)                                             |  | See M                  | See M                  | 1.1m                   | See M               | U708M-43R   |  |  |
|                                                                                             |                                                                                     |                        | <b>U410-91MR</b>       | Adult finger /3.0m     | D1                  | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U410S-91MR</b>      | Adult soft /3.0m       | D2                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U810-91MR</b>       | Multi-site Y/3.0m      | D1                  |             |                                                                                       |                                                                                       |
| MP20/30/40/50/60 70/90/VN4,VM6,VM8 (*Oximax)                                                |  |                        | <b>U410-91P</b>        | Adult finger /3.0m     | D1                  | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U410S-91P</b>       | Adult soft/3.0m        | D2                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U810-91P</b>        | Multi-site Y/3.0m      | D1                  |             |                                                                                       |                                                                                       |
| (Non Oximax)                                                                                |  |                        | <b>U410-01</b>         | Adult finger /3.0m     | C5                  | U710X-43    |  |  |
|                                                                                             |                                                                                     |                        | <b>U410S-01</b>        | Adult soft/3.0m        | C7                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U810-01</b>         | Multi-site Y/3.0m      | C5                  |             |                                                                                       |                                                                                       |
| (*Oximax)                                                                                   |  |                        | <b>U410-01XP</b>       | Adult finger /3.0m     | C7                  | U708-41     |  |  |
|                                                                                             |                                                                                     |                        | <b>U410S-01XP</b>      | Adult soft/3.0m        | C8                  |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                     |                        | <b>U810-01XP</b>       | Multi-site Y/3.0m      | C7                  |             |                                                                                       |                                                                                       |




| OEM Manufacturer & Models                                                                   | Sensor Connector                                                                      | REF# OEM | REF# UNIMED          | Sensor Descriptions     | Price Code | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------|----------------------|-------------------------|------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Philips/HP                                                                                  |                                                                                       |          |                      |                         |            |             |                                                                                       |                                                                                       |
| 78352C,78354C, 78834C,M1020A, M1025B,M1025A, M1350B,M1772A/B, M1732A/B,M2475B (*Non-Oximax) |    | DS100A   | U403-01              | Adult finger /1.1m      | C3         | U708-20     |    |    |
|                                                                                             |                                                                                       | DS100A   | U403S-01             | Adult soft /1.1m        | C5         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       | D-YS     | U803-01              | Multi-site Y/1.1m       | C3         |             |                                                                                       |                                                                                       |
|                                                                                             |    | M1196A   | U403-91<br>U410-91   | Adult finger /1.1m 3.0m | C8         | U708-40     |    |    |
|                                                                                             |                                                                                       | M1192A   | U103-91              | Pediatric finger /1.1m  | D1         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       | M1191A   | U403S-91<br>U410S-91 | Adult soft /1.1m 3.0m   | D1         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       | M1192A   | U103S-91             | Pediatric soft /1.1m    | D1         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       | M1195A   | U203S-91             | Infant soft /1.1m       | D1         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       | M1193A   | U303-91              | Neonate wrap /1.1m      | D1         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U603-91              | Adult/Neowrap/1.1m      | C8         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U603B-91             | Pedia/infant wrap/1.1m  | C8         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U803-91              | Multi-site Y/1.1m       | D1         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       | M1194A   | U903-91              | Adult ear clip /1.1m    | D1         |             |                                                                                       |                                                                                       |
| Schiller ODM/Bruker                                                                         |                                                                                       |          |                      |                         |            |             |                                                                                       |                                                                                       |
| Defigard4000, Defigard5000<br>ARGUS PRO LifeCare<br>Masimo SET                              |   |          | U410-63M             | Adult finger /3.0m      | D1         | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U410S-63M            | Adult soft /3.0m        | D2         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U810-63M             | Multi-site Y/3.0m       | D1         |             |                                                                                       |                                                                                       |
| ARGUS PRO LifeCare 2<br>DEFIGARD Touch 7<br>Masimo Rainbow                                  |  |          | U410-63MR            | Adult finger /3.0m      | D2         | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U410S-63MR           | Adult soft /3.0m        | D3         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U810-63MR            | Multi-site Y/3.0m       | D2         |             |                                                                                       |                                                                                       |
| Spacelabs                                                                                   |                                                                                       |          |                      |                         |            |             |                                                                                       |                                                                                       |
| Ultraview<br>(*Masimo Rainbow)                                                              |  |          | U410-74MR            | Adult finger /3.0m      | D2         | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U410S-74MR           | Adult soft /3.0m        | D3         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U810-74MR            | Multi-site Y/3.0m       | D2         |             |                                                                                       |                                                                                       |
| IBIO R&DT                                                                                   |                                                                                       |          |                      |                         |            |             |                                                                                       |                                                                                       |
| RD9012                                                                                      |  |          | U410-100             | Adult finger /3.0m      | C7         | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U410S-100            | Adult soft /3.0m        | C8         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U810-100             | Multi-site Y/3.0m       | C7         |             |                                                                                       |                                                                                       |
| ARGUS PRO LifeCare 2<br>DEFIGARD Touch 7<br>Masimo Rainbow                                  |  | DS100A   | U410-01              | Adult finger /3.0m      | C5         | No need     |  |  |
|                                                                                             |                                                                                       | FlexMAX  | U410S-01             | Adult soft /3.0m        | C7         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       | D-YS     | U810-01              | Multi-site Y/3.0m       | C5         |             |                                                                                       |                                                                                       |
| ☆Triton                                                                                     |                                                                                       |          |                      |                         |            |             |                                                                                       |                                                                                       |
| Defigard4000, Defigard5000<br>ARGUS PRO LifeCare<br>Masimo SET                              |  |          | U410-135             | Adult finger /3.0m      | C7         | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U410S-135            | Adult soft /3.0m        | C8         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U810-135             | Multi-site Y/3.0m       | C7         |             |                                                                                       |                                                                                       |
| ARGUS PRO LifeCare 2<br>DEFIGARD Touch 7<br>Masimo Rainbow                                  |  |          | U410-236             | Adult finger /3.0m      | C8         | No need     |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U410S-236            | Adult soft /3.0m        | D1         |             |                                                                                       |                                                                                       |
|                                                                                             |                                                                                       |          | U810-236             | Multi-site Y/3.0m       | C8         |             |                                                                                       |                                                                                       |

| OEM Manufacturer & Models | Sensor Connector                                                                    | REF# OEM | REF# UNIMED | Sensor Descriptions  | Price Code | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|---------------------------|-------------------------------------------------------------------------------------|----------|-------------|----------------------|------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| ☆UTAS                     |                                                                                     |          |             |                      |            |             |                                                                                       |                                                                                       |
| ☆UM-300                   |    |          | U410-129    | Adult finger /3.0m   | C8         |             |                                                                                       |                                                                                       |
|                           |                                                                                     |          | U410S-129   | Adult soft /3.0m     | D1         |             |                                                                                       |                                                                                       |
|                           |                                                                                     |          | U810-129    | Adult ear clip /3.0m | C8         |             |                                                                                       |                                                                                       |
| Unicare(China)            |                                                                                     |          |             |                      |            |             |                                                                                       |                                                                                       |
| ☆F8/S                     |    |          | U410-120    | Adult finger /3.0m   | C7         |             |                                                                                       |                                                                                       |
|                           |                                                                                     |          | U410S-120   | Adult soft /3.0m     | C8         |             |                                                                                       |                                                                                       |
|                           |                                                                                     |          | U810-120    | Adult ear clip /3.0m | C7         |             |                                                                                       |                                                                                       |
| Risingmed                 |                                                                                     |          |             |                      |            |             |                                                                                       |                                                                                       |
| RPM-9000A                 |    |          | U410-217    | Adult finger /3.0m   | C7         | No need     |                                                                                       |                                                                                       |
|                           |                                                                                     |          | U410S-217   | Adult soft/3.0m      | C8         |             |                                                                                       |                                                                                       |
|                           |                                                                                     |          | U810-217    | Multi-site Y/3.0m    | C7         |             |                                                                                       |                                                                                       |
|                           |    | DS100A   | U403-01     | Adult finger /1.1m   | C3         | U708-217    |    |    |
|                           |                                                                                     | DS100A   | U403S-01    | Adult soft/1.1m      | C5         |             |                                                                                       |                                                                                       |
|                           |                                                                                     | D-YS     | U803-01     | Multi-site Y/1.1m    | C3         |             |                                                                                       |                                                                                       |
| RGB medical               |                                                                                     |          |             |                      |            |             |                                                                                       |                                                                                       |
| (*BCI)                    |  | 3044     | U403-06     | Adult finger /1.1m   | C3         | U708-76     |  |  |
|                           |                                                                                     |          | U403S-06    | Adult soft /1.1m     | C3         |             |                                                                                       |                                                                                       |
|                           |                                                                                     | 3078     | U803-06     | Adult ear clip /1.1m | C3         |             |                                                                                       |                                                                                       |



| OEM Manufacturer & Models                       | Sensor Connector                                                                      | REF# OEM               | REF# UNIMED            | Sensor Descriptions | Price Code      | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|-------------------------------------------------|---------------------------------------------------------------------------------------|------------------------|------------------------|---------------------|-----------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Schiller ODM/Bruker                             |                                                                                       |                        |                        |                     |                 |             |                                                                                       |                                                                                       |
| Propaq CS, Encore Series, Atlas (* non-Oximax)  |    | DS100A                 | <b>U403-01</b>         | Adult finger /1.1m  | C3              | U708-01     |    |    |
|                                                 |                                                                                       | DS100A                 | <b>U403S-01</b>        | Adult soft /1.1m    | C5              |             |                                                                                       |                                                                                       |
|                                                 |                                                                                       | D-YS                   | <b>U803-01</b>         | Multi-site Y/1.1m   | C3              |             |                                                                                       |                                                                                       |
| TM 910-Defigard 500(M) (* non-Oximax)           |    | See Nellcor Non-Oximax | See Nellcor Non-Oximax | 1.1m/3.0m           | See Nellcor     | U708-63     |    |    |
| Argus VCM (*Oximax)                             |    | DS100A                 | <b>U403-01XP</b>       | Adult finger /1.1m  | C5              | U710-70     |    |    |
|                                                 |                                                                                       | DS100A                 | <b>U403S-01XP</b>      | Adult soft /1.1m    | C7              |             |                                                                                       |                                                                                       |
|                                                 |                                                                                       | D-YS                   | <b>U803-01XP</b>       | Multi-site Y/1.1m   | C5              |             |                                                                                       |                                                                                       |
| CM8 SM783,784,785, TM910                        |    | See RGB medical        | See RGB medical        | See RGB medical     | See RGB medical | U708-17     |    |    |
|                                                 |                                                                                       |                        | <b>U410-17</b>         | Adult finger /3.0m  | C8              | No need     |                                                                                       |                                                                                       |
|                                                 |                                                                                       |                        | <b>U410S-17</b>        | Adult soft /3.0m    | D1              |             |                                                                                       |                                                                                       |
|                                                 |    |                        | <b>U810-17</b>         | Multi-site Y/3.0m   | C8              |             |                                                                                       |                                                                                       |
| TM7,Argus LCM, (*M tech)                        |  | 1863/ LNCS DCI         | <b>U403-49</b>         | Adult finger /1.1m  | C4              | U708M-17    |  |  |
|                                                 |                                                                                       | 1863/ LNCS DCI         | <b>U403S-49</b>        | Adult soft /1.1m    | C6              |             |                                                                                       |                                                                                       |
|                                                 |                                                                                       | 2258/ LNCS YI          | <b>U803-49</b>         | Multi-site Y/1.1m   | C4              |             |                                                                                       |                                                                                       |
|                                                 |  |                        | <b>U410-17M</b>        | Adult finger /3.0m  | C8              | No need     |                                                                                       |                                                                                       |
|                                                 |                                                                                       |                        | <b>U410S-17M</b>       | Adult soft /3.0m    | D1              |             |                                                                                       |                                                                                       |
|                                                 |                                                                                       |                        | <b>U810-17M</b>        | Multi-site Y/3.0m   | C8              |             |                                                                                       |                                                                                       |
| Propaq CS, Encore, Micropaq (New) (M) (*M tech) |  | See M                  | See M                  | 1.1m                | See M           | U708M-15R   |  |  |
|                                                 |                                                                                       | LNOP DC-12             | <b>U410-15R</b>        | Adult finger /3.0m  | C8              | No need     |                                                                                       |                                                                                       |
|                                                 |  | LNOP DC-12             | <b>U410S-15R</b>       | Adult soft /3.0m    | D1              |             |                                                                                       |                                                                                       |
|                                                 |                                                                                       |                        | <b>U810-15R</b>        | Multi-site Y/3.0m   | C8              |             |                                                                                       |                                                                                       |

| OEM Manufacturer & Models                                | Sensor Connector                                                                    | REF# OEM               | REF# UNIMED            | Sensor Descriptions    | Price Code  | Ext - Cable | Cable Plug                                                                            | Sensor Conn.                                                                          |
|----------------------------------------------------------|-------------------------------------------------------------------------------------|------------------------|------------------------|------------------------|-------------|-------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Spacelabs                                                |                                                                                     |                        |                        |                        |             |             |                                                                                       |                                                                                       |
| 90419/01,90651A/07, 90501, Ultraview 90496 (*Novametrix) |    |                        | U410-03                | Adult finger /3.0m     | C8          | No need     |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U110-03                | Pediatric finger /3.0m | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U410S-03               | Adult soft /3.0m       | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U110S-03               | Pediatric soft /3.0m   | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U210S-03               | Infant soft /3.0m      | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U310-03                | Neonate wrap /3.0m     | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U610-03                | Adult/Neo wrap/3.0m    | C8          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U610B-03               | Pedia/infant wrap/3.0m | C8          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U810-03                | Multi-site Y/3.0m      | C8          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U910-03                | Adult ear clip /3.0m   | D1          |             |                                                                                       |                                                                                       |
| 90351-0/6,90465, 90466, 90456,PC exoress (*Non-Oximax)   |    | DS100A                 | U403-01                | Adult finger /1.1m     | C3          | U708-04     |    |    |
|                                                          |                                                                                     | DS100A                 | U403S-01               | Adult soft /1.1m       | C5          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     | D-YS                   | U803-01                | Multi-site Y/1.1m      | C3          |             |                                                                                       |                                                                                       |
|                                                          |  |                        | U410-04                | Adult finger /3.0m     | C8          | No need     |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U110-04                | Pediatric finger /3.0m | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U410S-04               | Adult soft /3.0m       | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U110S-04               | Pediatric soft /3.0m   | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U210S-04               | Infant soft /3.0m      | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U310-04                | Neonate wrap /3.0m     | D1          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U610-04                | Adult/Neo wrap/3.0m    | C8          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U610B-04               | Pedia/infant wrap/3.0m | C8          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U810-04                | Multi-site Y/3.0m      | C8          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U910-04                | Adult ear clip /3.0m   | D1          |             |                                                                                       |                                                                                       |
| Ultraview (*Non-Oximax)                                  |  | See Nellcor Non-Oximax | See Nellcor Non-Oximax | 1.1m                   | See Nellcor | U708-74     |  |  |
|                                                          |                                                                                     |                        | U410-74                | Adult finger /3.0m     | D1          |             |                                                                                       |                                                                                       |
|                                                          |  |                        | U410S-74               | Adult soft /3.0m       | D2          | No need     |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U810-74                | Multi-site Y/3.0m      | D1          |             |                                                                                       |                                                                                       |
| Ultraview (M Tech)                                       |  |                        | U410-74MR              | Adult finger /3.0m     | D1          | No need     |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U410S-74MR             | Adult soft /3.0m       | D2          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U810-74MR              | Multi-site Y/3.0m      | D1          |             |                                                                                       |                                                                                       |
| Ultraview (Oximax Tech)                                  |  |                        | U410-74P               | Adult finger /3.0m     | D1          | No need     |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U410S-74P              | Adult soft /3.0m       | D2          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U810-74P               | Multi-site Y/3.0m      | D1          |             |                                                                                       |                                                                                       |
| PRO2 (*BCI)                                              |  | 3044                   | U403-06                | Adult finger /1.1m     | C3          | No need     |                                                                                       |                                                                                       |
|                                                          |                                                                                     |                        | U403S-06               | Adult soft /1.1m       | C5          |             |                                                                                       |                                                                                       |
|                                                          |                                                                                     | 3043                   | U803-06                | Multi-site Y/1.1m      | C3          |             |                                                                                       |                                                                                       |

| OEM Manufacturer & Models                                        | Sensor Connector                                                                      | REF# OEM          | REF# UNIMED | Sensor Descriptions | Price Code | Ext - Cable | Cable Plug                                                                          | Sensor Conn.                                                                        |
|------------------------------------------------------------------|---------------------------------------------------------------------------------------|-------------------|-------------|---------------------|------------|-------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Takaoka                                                          |                                                                                       |                   |             |                     |            |             |                                                                                     |                                                                                     |
| Takaoka 9504,                                                    |    |                   | U410-99     | Adult finger /3.0m  | C7         | No need     |                                                                                     |                                                                                     |
| 9504ECG, Oxipani ,                                               |                                                                                       |                   | U410S-99    | Adult soft /3.0m    | C8         |             |                                                                                     |                                                                                     |
| Oxifast                                                          |                                                                                       |                   | U810-99     | Multi-site Y/3.0m   | C7         |             |                                                                                     |                                                                                     |
| Uray                                                             |                                                                                       |                   |             |                     |            |             |                                                                                     |                                                                                     |
| UR-8000                                                          |    |                   | U410-143    | Adult finger /3.0m  | C7         | No need     |                                                                                     |                                                                                     |
| UR-9000                                                          |                                                                                       |                   | U410S-143   | Adult soft /3.0m    | C8         |             |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U810-143    | Multi-site Y/3.0m   | C7         |             |                                                                                     |                                                                                     |
| Zoll                                                             |                                                                                       |                   |             |                     |            |             |                                                                                     |                                                                                     |
| M series AED<br>(*M tech)                                        |    | 1863/<br>LNCS DCI | U403-49R    | Adult finger /1.1m  | C4         | U708M-15R   |  |  |
|                                                                  |                                                                                       | 1863/<br>LNCS DCI | U403S-49R   | Adult soft /1.1m    | C6         |             |                                                                                     |                                                                                     |
|                                                                  |                                                                                       | 2258/<br>LNCS YI  | U803-49R    | Multi-site Y/1.1m   | C4         |             |                                                                                     |                                                                                     |
|                                                                  |    | LNOP DC-12        | U410-15R    | Adult finger /3.0m  | C8         | No need     |                                                                                     |                                                                                     |
|                                                                  |                                                                                       | LNOP DC-12        | U410S-15R   | Adult soft /3.0m    | D1         |             |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U810-15R    | Multi-site Y/3.0m   | C8         |             |                                                                                     |                                                                                     |
| Zondan (China)                                                   |                                                                                       |                   |             |                     |            |             |                                                                                     |                                                                                     |
| Apollo N2, Apollo N4<br>Apollo N5<br>(M tech optional)<br>DS120C |   |                   | U410-45     | Adult finger /3.0m  | C7         | No need     |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U410S-45    | Adult soft /3.0m    | C8         |             |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U810-45     | Multi-site Y/3.0m   | C7         |             |                                                                                     |                                                                                     |
| ZD120B<br>ZD120D<br>ZD120E<br>(*Digital Tech)                    |  |                   | U410-66D    | Adult finger /3.0m  | C8         | No need     |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U410S-66D   | Adult soft /3.0m    | D1         |             |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U810-66D    | Multi-site Y/3.0m   | C8         |             |                                                                                     |                                                                                     |
| ZD120D/ZD120E                                                    |  |                   | U410-45     | Adult finger /3.0m  | C7         | No need     |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U410S-45    | Adult soft /3.0m    | C8         |             |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U810-945    | Multi-site Y/3.0m   | C7         |             |                                                                                     |                                                                                     |
| Apollo<br>N1/N2/N3/N4/N5                                         |  |                   | U410-210M   | Adult finger /3.0m  | D1         | No need     |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U410S-210M  | Adult soft /3.0m    | D2         |             |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U810-210M   | Multi-site Y/3.0m   | D1         |             |                                                                                     |                                                                                     |
| Unicare(China)                                                   |                                                                                       |                   |             |                     |            |             |                                                                                     |                                                                                     |
| (*Digital tech)                                                  |  |                   | U410-66D    | Adult finger /3.0m  | C8         | No need     |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U410S-66D   | Adult soft /3.0m    | D1         |             |                                                                                     |                                                                                     |
|                                                                  |                                                                                       |                   | U810-66D    | Multi-site Y/3.0m   | C8         |             |                                                                                     |                                                                                     |

Please be Reminded : The other types of sensors like U100, U100S, U200S, U300, U600, U600B and U900 of the listed respected brands are available upon demands.










### ③ SpO2 Interface Cable





#### Specification

|                |                                |
|----------------|--------------------------------|
| Cable material | TPU Jacket, Light gray, Blue   |
| Connector      | Nylon plug &Golden plated pins |
| Cable length   | 1.1m/2.2m/3.0m                 |
| Cable diameter | 4.0mm/5.0mm                    |
| Package        | 1 pcs per bag                  |





| REF#      | OEM REF#                                    | Sensor Tech.           | Cable length | Monitor end Conn. | Sensor end Conn. | Price code |
|-----------|---------------------------------------------|------------------------|--------------|-------------------|------------------|------------|
| Artema    |                                             |                        |              |                   |                  |            |
| U708-80   | 59831-0                                     | Nellcor,non-Oximax     | 2.2m         |                   |                  | C3         |
| BCI       |                                             |                        |              |                   |                  |            |
| U708-76   |                                             | BCI                    | 2.2m         |                   |                  | B8         |
| U708-06   |                                             | BCI                    | 2.2m         |                   |                  | B8         |
| Biolight  |                                             |                        |              |                   |                  |            |
| U710X-51  |                                             | Nellcor Oximax         | 3.0m         |                   |                  | C3         |
| U708-68   |                                             | Digital sensor         | 2.2m         |                   |                  | C3         |
| Bionet    |                                             |                        |              |                   |                  |            |
| U708-58   |                                             | Nellcor,non-Oximax     | 2.2m         |                   |                  | B8         |
| U708-108  |                                             | Bionet                 | 2.2m         |                   |                  | C3         |
| Criticare |                                             |                        |              |                   |                  |            |
| U708-05   | 518LD                                       | CSI                    | 2.2m         |                   |                  | C3         |
| U708-75   | 518DD                                       | CSI                    | 2.2m         |                   |                  | B8         |
| Datascope |                                             |                        |              |                   |                  |            |
| U708-07   | 0012-00-0516-02                             | DATASCOPE              | 2.2m         |                   |                  | C3         |
| Dolphin   |                                             |                        |              |                   |                  |            |
| U708-50   |                                             | Dolphin digital sensor | 2.2m         |                   |                  | C3         |
| Drager    |                                             |                        |              |                   |                  |            |
| U708-28   | 5720084                                     | Nellcor,non-Oximax     | 2.2m         |                   |                  | D4         |
| EDAN      |                                             |                        |              |                   |                  |            |
| U708X-191 | MS1-30131<br>(M9,M9A,M8,M8A,M8B)            | Nellcor-Oximax         | 2.2m         |                   |                  | C3         |
| U708-29   | M13-30691<br>(M9,M9A,M8,M8A,M8B,<br>M3,M3B) | EDAN Oximetry          | 2.2m         |                   |                  | B8         |








| REF#                 | OEM REF#                                   | Sensor Tech.                                                  | Cable length | Monitor end Conn.                                                                    | Sensor end Conn.                                                                      | Price code |
|----------------------|--------------------------------------------|---------------------------------------------------------------|--------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------------|
| GE healthcare        |                                            |                                                               |              |                                                                                      |                                                                                       |            |
| U708-21              | 2006644-001<br>E9004GE                     | Nellcor,non-Oximax                                            | 2.2M         |    |    | C3         |
| U710X-21             | 2021406-001                                | Nellcor,Oximax                                                | 3.0M         |    |    | C7         |
| U708-21M             | M technology                               | Nellcor Non-Oximax                                            | 2.2M         |    |    | C3         |
| U708M-21R            | LNC-10-GE<br>2027263-002<br>(REF 2016)     | M tech.(Masimo LNCS)                                          | 2.2M         |    |    | C3         |
| U708M-21P            | 20169041-001                               | M tech.(Masimo LNOP)                                          | 2.2M         |    |    | C7         |
| U708-09              | OXY-C3                                     | Datex                                                         | 2.2m         |    |    | C3         |
| U708-81              | Datex>>Oxytip<br>Used with<br>Oxy-C3 cable | Oxytip                                                        | 2.2m         |    |    | C3         |
| U708-55              | 30344358                                   | Hellige(Nellcor)                                              | 2.2m         |    |    | C3         |
| U708-25              | 4033CAX                                    | Corometrics(Nellcor)                                          | 2.2m         |    |    | C3         |
| U708-72*             | All Ohmeda<br>Oximeters                    | Ohmeda Used with<br>Unimed N500 or U500<br>Disposable sensors | 2.2m         |   |   | C3         |
| U708-160             | OXY-MC3                                    | GE OxyTip+ sensor                                             | 2.2m         |  |  | C7         |
| U708-102<br>U704-102 | OXY-OL3                                    | GE OxyTip+ sensor                                             | 2.2m<br>1.1m |  |  | C5         |
| U708-209             | OXY-SL3                                    | GE OxyTip+ sensor                                             | 2.2m         |  |  | C5         |
| U708-121<br>U710-121 | OXY-ES3                                    | GE OxyTip+ sensor                                             | 2.2m<br>3.0m |  |  | C7         |
| U708-81<br>U701-81   | OXY-SLA<br>OXY-SLC                         | GE OxyTip+ sensor                                             | 2.2m<br>0.3m |  |  | C5<br>C3   |
| U708-111             | TS-G3                                      | GE OxyTip+ sensor                                             | 2.2m         |  |  | C5         |
| U708-112             | TS-H3                                      | GE OxyTip+ sensor                                             | 2.2m         |  |  | C5         |
| U708-116             | TS-M3                                      | GE OxyTip+ sensor                                             | 2.2m         |  |  | C7         |
| U708-119             | TS-N3                                      | GE OxyTip+ sensor                                             | 2.2m         |  |  | C5         |
| Goldway              |                                            |                                                               |              |                                                                                      |                                                                                       |            |
| U708-11              |                                            | Nellcor Non-Oximax                                            | 2.2m         |  |  | C3         |
| U708-37              |                                            | Nellcor Non-Oximax                                            | 2.2m         |  |  | B8         |

| REF#         | OEM REF#                                                                      | Sensor Tech.                     | Cable length | Monitor end Conn.                                                                     | Sensor end Conn.                                                                      | Price code |
|--------------|-------------------------------------------------------------------------------|----------------------------------|--------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------------|
| Goldway      |                                                                               |                                  |              |                                                                                       |                                                                                       |            |
| U708-38      |                                                                               | Nellcor Non-Oximax               | 2.2m         |    |    | B8         |
| U710X-38     |                                                                               | Nellcor Oximax                   | 3.0m         |    |    | C3         |
| U708-11      |                                                                               | Nellcor Non-Oximax               | 2.2m         |    |    | C3         |
| U708-37      |                                                                               | Nellcor Non-Oximax               | 2.2m         |    |    | B8         |
| Invivo       |                                                                               |                                  |              |                                                                                       |                                                                                       |            |
| U708-14      |                                                                               | Nellcor,non-Oximax               | 2.2m         |    |    | C3         |
| U708-14M     |                                                                               | (Masimo)<br>Nellcor, non-Oximax  | 2.2m         |    |    | C3         |
| U708M-14     |                                                                               | Masimo LNCS sensor               | 2.2m         |    |    | C3         |
| Kontron      |                                                                               |                                  |              |                                                                                       |                                                                                       |            |
| U708-13      |                                                                               | Nellcor,non-Oximax               | 2.2m         |    |    | C3         |
| U708-53      |                                                                               | Nellcor,non-Oximax               | 2.2m         |    |    | B8         |
| M Technology |                                                                               |                                  |              |                                                                                       |                                                                                       |            |
| U708-15      | M technology                                                                  | (Masimo)<br>Nellcor, non-Oximax  | 2.2m         |  |  | C3         |
| U708M-15R    | LNCS Patient Cable<br>LNC-4(REF 2017)<br>LNC-10(REF 1814)<br>LNC-14(REF 2013) | Masimo LNCS sensor               | 2.2m         |  |  | C3         |
| U708M-15PR   | PC 04 / PC 08 / PC 12                                                         | Masimo LNOP sensor               | 2.2m         |  |  | C3         |
| U704M-49     | LNC-4-Ext<br>(REF 2021)<br>For LNCS Cable                                     | Masimo LNCS sensor               | 1.1m         |  |  | B8         |
| U701M-125    | LNCS to LNOP<br>LNCS to PC series<br>(REF 1816)                               | Masimo LNCS sensor               | 0.3m         |  |  | C3         |
| MEK          |                                                                               |                                  |              |                                                                                       |                                                                                       |            |
| U708-56      |                                                                               | MEK<br>(Nellcor non-oximax)      | 2.2m         |  |  | C3         |
| Mennen       |                                                                               |                                  |              |                                                                                       |                                                                                       |            |
| U708-54      |                                                                               | Nellcor,non-Oximax               | 2.2m         |  |  | C3         |
| U708-24      |                                                                               | Nellcor,non-Oximax               | 2.2m         |  |  | D4         |
| U708-59M     | M technology                                                                  | M tech.<br>(Nellcor, non-Oximax) | 2.2m         |  |  | C3         |
| U708M-59     | M technology                                                                  | Masimo LNCS sensor               | 2.2m         |  |  | C3         |
| U710X-59     |                                                                               | Nellcor Oximax                   | 3.0m         |  |  | C3         |

| REF#            | OEM REF#                                                                                                                                       | Sensor Tech.           | Cable length | Monitor end Conn.                                                                    | Sensor end Conn.                                                                      | Price code |
|-----------------|------------------------------------------------------------------------------------------------------------------------------------------------|------------------------|--------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------------|
| Mindray         |                                                                                                                                                |                        |              |                                                                                      |                                                                                       |            |
| U708-29         | 0010-20-42594<br>Accutorr V, DPM 3, DPM 4,<br>DPM 5, Passport V                                                                                | Nellcor,non-Oximax     | 2.2M         |    |    | B8         |
| U708-29M        | M technology<br>For old PM 6000 only                                                                                                           | Nellcor,non-Oximax     | 2.2M         |    |    | B8         |
| U708M-29        | M technology<br>For old PM 6000 only                                                                                                           | Masimo LNCS            | 2.2M         |    |    |            |
| U708-30M        | M technology<br>(0010-30-42625)<br>DPM 3                                                                                                       | Nellcor,non-Oximax     | 2.2M         |    |    | B8         |
| U708M-30        | DPM 4<br>DPM 5                                                                                                                                 | Masimo LNCS            | 2.2M         |    |    |            |
| U708-22         | 0515-30-11221                                                                                                                                  | Nellcor,non-Oximax     | 2.2M         |    |    | C3         |
| U708-22M        | M technology<br>0010-30-12452                                                                                                                  | non-Oximax Nellcor     | 2.2M         |    |    | C3         |
| U708M-22        | M technology<br>0010-30-12452                                                                                                                  | Masimo LNCS            | 2.2M         |    |    | C3         |
| U708-31         | CSI Technology                                                                                                                                 | CSI                    | 2.2M         |   |   | C3         |
| U708M-48R       | 0010-30-42738<br>T5/T8<br>DPM 6, DPM 7                                                                                                         | Masimo LNCS            | 2.2m         |  |  | C3         |
| U708-48         | 0010-20-42710 (C)<br>T5/T8<br>DPM 6, DPM 7<br>Nellcor Non-Oximax                                                                               | Nellcor,non-Oximax     | 2.2m         |  |  | C3         |
| U710X-48        | 0010-20-42712<br>T5/T8, DPM 6, DPM 7<br>Nellcor Oximax                                                                                         | Nellcor Oximax         | 3.0m         |  |  | C4         |
| U708M-15R       | (0012-00-1652),1.2m<br>(0012-00-1599),3.1m<br>(0012-00-1653),4.3m<br>Accutorr V, Duo, Trio,<br>Passport 2, PassportV,<br>Spectrum/OR, V12, V21 | Masimo LNCS            |              |  |  | C3         |
| U708-69         | 0012-00-1254<br>Trio, Passport 2,<br>Spectrum                                                                                                  | Nellcor,non-Oximax     | 2.2m         |  |  | C3         |
| U710-70         | 0012-00-1464<br>Accutorr V, Duo, Trio,<br>Passport 2, Passport V,<br>Spectrum/OR, V12, V21                                                     | Nellcor Oximax         | 3.0m         |  |  | C6         |
| U710X-180       | 0010-20-42595<br>DPM4                                                                                                                          | Nellcor Oximax         | 3.0m         |  |  | C3         |
| Nellcor/Coviden |                                                                                                                                                |                        |              |                                                                                      |                                                                                       |            |
| U708-69         | SCP-10/MC-10                                                                                                                                   | Nellcor,<br>Non-Oximax | 2.2m         |  |  | C3         |
| U710-85         | M-200-13                                                                                                                                       |                        | 3.0m         |  |  | D7         |
| U704-71         | DEC-4                                                                                                                                          | Nellcor-Oximax         | 1.1m         |  |  | B8         |
| U708-71         | DEC-8                                                                                                                                          |                        | 2.2m         |  |  | B8         |
| U710-70         | DOC-10                                                                                                                                         | Nellcor-Oximax         | 3.0m         |  |  | C6         |
| U710-70X        | DOC-10 +<br>Oximax adapter                                                                                                                     | Nellcor- Non Oximax    | 3.0m         |  |  | C8         |

| REF#                 | OEM REF#                               | Sensor Tech.        | Cable length | Monitor end Conn.                                                                     | Sensor end Conn.                                                                      | Price code |
|----------------------|----------------------------------------|---------------------|--------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------------|
| Nellcor/Coviden      |                                        |                     |              |                                                                                       |                                                                                       |            |
| U701-01X<br>U708-01X | Oximax adapter                         | Nellcor, Non-Oximax | 0.3m<br>2.2m |    |    | C3         |
| Nihon Kohden         |                                        |                     |              |                                                                                       |                                                                                       |            |
| U710-36              | JL-900P<br>BSM/PVM                     | NIHON KOHDEN        | 3.0m         |    |    | C7         |
| U710X-36             | JL-650P                                | Oximax              | 3.0m         |    |    | C7         |
| U710M-36R            | JL-631P                                | M Tech              | 3.0m         |    |    | C7         |
| U710-106             | JL-302T<br>OPV-1500,OLV                | NIHON KOHDEN        | 3.0m         |    |    | C7         |
| Nonin                |                                        |                     |              |                                                                                       |                                                                                       |            |
| U708-08              |                                        | Nonin               | 2.2m         |    |    | C3         |
| Novametrix           |                                        |                     |              |                                                                                       |                                                                                       |            |
| U708-103             | 8853                                   | Novametrix          | 2.2m         |    |    | C3         |
| Pace tech/Generra    |                                        |                     |              |                                                                                       |                                                                                       |            |
| U708-37              |                                        | Nellcor non-Oximax  | 2.2m         |   |   | B8         |
| U708-38              |                                        | Nellcor non-Oximax  | 2.2m         |  |  | B8         |
| U708-39              |                                        | Nellcor non-Oximax  | 2.2m         |  |  | C3         |
| Philips/HP           |                                        |                     |              |                                                                                       |                                                                                       |            |
| U708-40              | M1940A/B                               | Philips             | 2.2m         |  |  | C3         |
| U708-41              | M1941A/B                               |                     | 2.2m         |  |  | C3         |
| U708-20              | M1900A/B                               | Nellcor non-Oximax  | 2.2m         |  |  | C3         |
| U708M-20P            | M1020-61102<br>CMS M1020A module.      | Masimo LNOP         | 2.2m         |  |  | C3         |
| U708-43<br>U710-43   | M1943A<br>M1943AL                      | Nellcor non-Oximax  | 2.2m<br>3.0m |  |  | C3         |
| U710X-43             | M1943NL                                | Nellcor,Oximax      | 3.0m         |  |  | C7         |
| U708-43MR            | M Technology                           | Nellcor non-Oximax  | 2.2m         |  |  | C3         |
| U708M-43R            | 989803148221<br>LNC MP10<br>(REF 2281) | Masimo LNCS sensor  | 2.2m         |  |  | C3         |
| U708M-43PR           | M1020-61100                            | Masimo LNOP sensor  | 2.2m         |  |  | C7         |

4 Accessories

| REF#                   | OEM REF#                                                            | Sensor Tech.                   | Cable length | Monitor end Conn.                                                                    | Sensor end Conn.                                                                      | Price code |
|------------------------|---------------------------------------------------------------------|--------------------------------|--------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|------------|
| Schiller               |                                                                     |                                |              |                                                                                      |                                                                                       |            |
| U708-17                |                                                                     | BCI                            | 2.2m         |    |    | C3         |
| U708-17M               | M tech.                                                             | Nellcor non-Oximax             | 2.2m         |    |    | C3         |
| U708M-17               | M tech.                                                             | Masimo LNCS sensor             | 2.2m         |    |    | C3         |
| U708-63<br>ODAM/Bruker |                                                                     | BCI                            | 2.2m         |    |    | C3         |
| Siemens/Dräger         |                                                                     |                                |              |                                                                                      |                                                                                       |            |
| U704-23                | 3368433                                                             | Nellcor non-Oximax             | 1.1m         |    |    | C6         |
| U708-23                | 3375834                                                             | Nellcor non-Oximax             | 2.2m         |    |    | C6         |
| U710X-23               |                                                                     | Nellcor Oximax                 | 3.0m         |    |    | C7         |
| U704M-23P              | MS18680                                                             | Masimo <sup>®</sup> LNOP       | 1.1m         |    |    | C6         |
| U708M-23PR             | MS17041                                                             |                                | 2.2m         |   |   | C6         |
| U704M-23R              | MS24303                                                             | Masimo <sup>®</sup> LNCS       | 1.1m         |  |  | C6         |
| U708M-23R              | MS17522                                                             |                                | 2.2m         |  |  | C6         |
| U708-33                | 4533840                                                             | Nellcor non-Oximax             | 2.2m         |  |  | C3         |
| Spacelabs              |                                                                     |                                |              |                                                                                      |                                                                                       |            |
| U708-04                | 700-0002-00<br>90351-D/G,<br>90465, 90466,<br>90467, 90489,M77      | Nellcor non-Oximax             | 2.2m         |  |  | C3         |
| U708-74                | 700-0030-00<br>90496,91496-U                                        | Nellcor non-Oximax             | 2.2m         |  |  | C6         |
| U710X-74               | 700-0792-00<br>91496-N,91220                                        | Nellcor-Oximax                 | 3.0m         |  |  | C7         |
| *U708M-74              | 700-0906-00<br>91496-M                                              | Masimo LNCS                    | 2.2m         |  |  | C6         |
| *U708M-74P             | 700-0789-00<br>91496-M                                              | Masimo LNOP                    | 2.2m         |  |  | C6         |
| U701-04NV              | 175-0646-00<br>90351-D/G,90465,<br>90466,90467,<br>90489, 90651-08) | Novametrix-<br>Hypertronic Con | 0.3m         |  |  | C3         |
| U701-74NV              | 700-0029-00<br>(90496,91496-U)                                      | Novametrix-<br>Hypertronic Con | 0.3m         |  |  | C6         |
| Welchallyn             |                                                                     |                                |              |                                                                                      |                                                                                       |            |
| U704X-104              | 001966-E<br>(PIC 40TM Defibrillator<br>PIC 50TM Defibrillator)      | Nellcor Oximax                 | 2.2m         |  |  | C7         |
| U704-104               | 001955-U<br>(PIC40 Defibrillator<br>PIC50 Defibrillator)            | Nellcor                        | 1.1m         |  |  | C7         |

\*Refurbish only

|  |                                                              |                                                                               |                   |          |         |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------|-------------------------------------------------------------------------------|-------------------|----------|---------|
| REF#                                                                                | Product Name                                                 | Purpose of Use                                                                | Package           | Price    | Figures |
| Y04-E041                                                                            | Ear sensor hanger                                            | Used with Unimed ear clip sensor U900 series or Other brands ear sensor(Gray) | 1 pcs             | A5       | A       |
| UEC-01                                                                              | Adult ear clip Or as animal small size lingua clip           | Used with D-YS or Unimed U800 sensor(Black)                                   | 1 pcs             | A5       | B       |
| UVC-01                                                                              | Animal lingua clip,Big size                                  | Used with D-YS or Unimed U800 sensor(Gray)                                    | 1 pcs             | A8       | C       |
| UVC-02                                                                              | Animal lingua clip,Middle size                               | Used with D-YS or Unimed U800 sensor(Gray)                                    | 1 pcs             | A7       | D       |
| Y04-E019                                                                            | Reusable Silicon Y sensor wrap tape                          | Used with Nellcor D-YS or Unimed U800 series sensors                          | 1 pcs             | A6       | E       |
| E023-12                                                                             | Reusable 3 holes foam tape<br>For All patient sizes          | Used with Nellcor D-YS or Unimed U800 /U600 series sensors                    | 12 pcs<br>Per bag | B5/bag   | F       |
| E024-25                                                                             | Disposable 3M adhesive textile tape<br>For Adult and neonate | Used with Nellcor D-YS or Unimed U800 series sensors                          | 25 pcs<br>Per bag | B7/bag   | G       |
| E025-25                                                                             | Disposable 3M adhesive textile tape for Pediatric and infant | Used with Nellcor D-YS or Unimed U800 series sensors                          | 25 pcs<br>Per bag | B7/bag   | H       |
| T002-12                                                                             | Velcro foam wrap tape                                        | To fix the sensor cable onto the patient wrist                                | 12 pcs<br>per bag | C1.4/bag | I       |