




1.2 Temp. Adapter Cable for Reusable Probe

Reference Number Rule Description

TSM-AD/2	
T	Temperature probe
SM	Factory/model (See other on previous page)
2	Dual channel(Default is single channel)
AD	Adapter, used with reusable YSI probes




OEM & Models	Unimed REF#	Application	Cable length	Monitor end plug	Probe end Connector	Price code
Artema-S&W						
	TSW-AD	Used with reusable YSI400 compatible probes	0.3m			C3
Biolight						
V6 Vital Sign AnyView A2/A3/A5/A6/A8	TBLT-AD	Used with reusable YSI400 compatible probes	0.3m			C3
GE Marquette						
Eagle, Tram, Solar, Mac-lab Cath lab	TMQ-AD	Used with reusable YSI400 compatible probes	0.3m			C3
	TMQ400-AD/2 402015-004	Used with reusable YSI400 compatible probes	0.3m			C4
	TMQ700-AD/2 402015-004	Used with reusable YSI700 compatible probes	0.3m			C4
Philips Medical						
All	THP-AD 989803163281	Used with reusable YSI400 compatible probes	0.3m			C3
Mindray						
T5/T8 DPM 6, DPM 7	TMR-AD	Used with reusable YSI400 compatible probes	0.3m			C3
Datascope Passport V, V12, V21	TDT-AD	Used with reusable YSI400 compatible probes	0.3m			C3
Siemens/Drager						
	TSM-AD 5198333	Used with reusable YSI400 compatible probes	0.3m			C3
Spacelabs						
	TSL-AD	Used with reusable YSI400 compatible probes	0.3m			C5
90496, 91496	TSL-AD/2 700-0031-00	Used with reusable YSI400 or YSI700 compatible probes	0.3m			D5
90496/91496 306655-002/306644-001	CO-SL-AD 700-0027-00	CO adapter cable	0.3m			D3



2.5 Temp. Adapter Cable For Disposable probe

Reference Number Rule Description

THP-30-AD	
T	Temperature Adapter Cable
HP	Factory/model (See others on the previous page)
30	3.0m cable length
AD	used with YSI 400 series disposable probe(PH: Philips disposable probe)

OEM & Models	Probe end Conn	UNIMED REF#	Probe end Conn	UNIMED REF#	Cable length	Monitor end plug	Price code
	For 400 Disposable		For Philips Disposable				
Artema-S&W							
		TSW-30-AD		TSW-30-PH	3.0m		C3
Biolight							
		TBLT-30-AD		TBLT-30-PH	3.0m		C3
GE Marquette							
Eagle, Tram, Solar, Mac-lab Cath lab		TMQ-30-AD 2021701-001		TMQ-30-PH	3.0m		C3
Philips Medical							
All		THP-30-AD		THP-30-PH (21082A),3.0m (21082B),1.5m	3.0m		C3
Mindray							
T5/T8 DPM 6, DPM 7		TMR-30-AD		TMR-30-PH	3.0m		C3
Passport V, V12, V21		TDT-30-AD 040-000091-00		TDT-30-PH	3.0m		C3
Siemens							
		TSM-30-AD		TSM-30-PH	3.0m		C3
Spacelabs							
		TSL-30-AD		TSL-30-PH	3.0m		C5
YSI 400							
		TYSI-30-AD (165640),1.3m (165641),2.8m (9445-801),3.6m		TYSI-30-PH (989803162591)3.0m (989803162601)1.5m	3.0m		B8