

Fetal Transducers

Philips 50 Series



Philips Ultrasound

Compatible Monitors:
1350A, Viridia Series 50XM

OEM P/N: M1356A
UNIMED P/N: UFU200-20

Connector
RED



Philips Tocodynamic

Compatible Monitors:
1350A, Viridia Series 50XM

OEM P/N: M1355A
UNIMED P/N: UFT200-10

Connector
TAN



Philips Avalon FM Series



Philips Avalon Ultrasound

Compatible Monitors:
FM20, FM30

OEM P/N: M2734A
Original

Connector
RED



Philips Avalon Tocodynamic

Compatible Monitors:
FM20, FM30

OEM P/N: M2736A
Original

Connector
RED



Fetal Transducers

Corometrics 170 Series

GE Corometrics Ultrasound

Compatible Monitors:
115, 118, 120, 170

OEM P/N: 5700HAX (button style)
UNIMED P/N: UFU300-20

Connector
GREY



GE Corometrics Tocodynamic

Compatible Monitors:
115, 118, 120, 170

OEM P/N: 2264HAX (button style)
UNIMED P/N: UFT300-20

Connector
WHITE



Fetal Transducers

Accessories

P/N: E044
Fetal Transducer Belt



Features

- Reusable
- Stretchable with holes & locking button
- Size: 1300mm (L) x 60mm (W)
- Colour: Blue
- Latex Free