

Resuscitators

Reusable - Silicone



BS5401R

Adult

1400ml

Size 5 mask



BS3601

Paediatric

550ml

Size 3 mask



BS1601

Infant

330ml

Size 1 mask

Features

- Textured surface for easy grip
- Pressure Release Valve for Patient Safety
- CE, ISO & FDA approved
- 100% Latex Free
- Silicone - Reusable
- Autoclavable at 134°C

Resuscitators

Disposable - PVC



BV5401R

Adult

1400ml

Size 5 mask




BV3601

Paed

550ml

Size 3 mask



BV1601

Infant

330ml

Size 1 mask

Features

- Textured surface for easy grip
- Pressure Release Valve for Patient Safety
- CE, ISO & FDA approved
- 100% Latex Free
- PVC - Disposable

Resuscitators

Optional Accessories - Reusable



Reservoir Bag

BP3100	BP3500
Adlt/ Paed	Infant
2.9L	810ml
PVC - Not Autoclavable	



PEEP Valve Adaptor

BP620S
Autoclavable



PEEP Valve

BP6320	BP6310
Adult	Paed/Infant
5-20 H2O cm	0-10 H2O cm